

Analysts Report Amendment

Study No.: 275.30

Amendment No.: 1


Effective Date: 3 Feb 2011

Study Title: Provide data on various arsenic species present in broilers treated with roxarsone: Comparison with untreated birds


Changes:

The following minor errors were noted by the Quality Assurance Unit in review of the Analysts Report after it had been signed:


1. On page 10, line 203, the Jackson and Bertsch reference number should be (8), not (7).
2. On page 20, line 359, should read "Tables 3a though 3d" not 3c.
3. On page 23, line 424, should read "(Table 3b)" not 3a.




Analyst Mary C. Carson (CVM)



Date




Analyst Sean D. Conklin (CFSAN)




Date




Philip J. Kijak, Director, Division of Residue Chemistry (CVM)



Date



Stephen G. Capar, Chief, Chemical Contaminants Branch (CFSAN)



Date

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30

Office of Research Final Report:

Analysts' Report for Liver, Feed and Premix Analyses

*Provide data on various arsenic species present in broilers treated with roxarsone:
Comparison with untreated birds.*

OR Study 275.30

Study Director: Joseph Kawalek, Ph.D.
OR Division: HFV-520, Division of Animal Research

Analysts: Mary C. Carson, Ph.D. HFV-510
Sean D. Conklin, Ph.D. HFS-716
Karyn Howard HFV-520

Analytical Phase Initiation Date: 13 Oct 2009

Facility Addresses: Division of Residue Chemistry
Office of Research
Center for Veterinary Medicine
Food and Drug Administration
8401 Muirkirk Road
Laurel, Maryland 20708

and

Chemical Contaminants Branch
Division of Bioanalytical Chemistry
Office of Regulatory Science
Center for Food Safety and Applied Nutrition
Food and Drug Administration
5100 Paint Branch Parkway
College Park, Maryland 20740

31 INDEX

32 I. Abstract 4

33 II. Narrative..... 5

34 Introduction..... 5

35 Personnel..... 8

36 Procedures..... 8

37 Total As determination..... 8

38 Speciation in Liver 9

39 Analysis Dates..... 16

40 Statistical Methods and Calculations 16

41 Summary of Results 17

42 Total Arsenic Determination..... 17

43 Arsenic Speciation in Liver..... 18

44 Detection and Quantification Limits of the Speciation Method 18

45 Speciation Method Evaluation 20

46 Speciation Analysis of Roxarsone-Treated and Control Livers..... 22

47	Extraction Efficiency and Mass Balance of Arsenic.....	23
48	Speciation in Control Feed, Medicated Feed, and Type A Medicated Articles.....	24
49	Conclusions.....	26
50	Signatures	27
51	III. Tables	28
52	Table 1. Dates for total arsenic analyses.....	28
53	Table 2. Total arsenic in liver results.....	29
54	Tables 3a-d. Validation results for arsenic speciation method	34
55	Tables 4a-b. Arsenic species in liver.....	39
56	Tables 5a-b. Arsenic species in feed and type A medicated articles	46
57	IV. List of Excel Spreadsheets Contributing to Tables 2, 3, 4, and 5	48
58	V. References.....	51
59	VI. Attachment 1: SOP for ICP-MS analysis for total arsenic	
60	VII. Attachment 2: SOP for Speciation of arsenic compounds related to roxarsone use in	
61	chickens	
62	VIII. Electronic files	

63 **I. Abstract**

64 We were tasked with developing and applying methods to characterize the arsenic (As) profiles
65 in tissue and excreta from control chickens and those treated with roxarsone. A literature search
66 showed that ion chromatography-inductively coupled plasma-mass spectrometry (IC-ICP-MS) is
67 the most sensitive and versatile instrumentation for this purpose. Validated methods to
68 accomplish As speciation in roxarsone treated poultry were not available. The animal phase of
69 Study 275.30 was completed before we could develop such methods. We determined total As in
70 the breast and liver samples from the study birds using a standard method. Total As
71 concentrations were much higher (ca 40x) in liver than in breast, so method development
72 focused on liver. A preliminary speciation method developed using standards and fortified
73 control tissue proved inadequate when used with the study samples due to the presence of
74 unknown As compounds (presumed roxarsone metabolites) which interfered with the
75 chromatography. The IC-ICP-MS speciation method was refined with the primary purpose of
76 identifying and quantifying trace levels of arsenite (AsIII) or arsenate (AsV) in liver in the
77 presence of large amounts of roxarsone and other unknown organic As species. This method
78 was used to characterize the study liver samples. Livers from untreated birds did not have any
79 As species above the method's lower limit of quantification. Livers from treated birds all had
80 roxarsone present, as well as several other As-containing compounds. A modified version of the
81 tissue speciation method was also applied to feed and type A medicated articles to further
82 characterize them.

83 **II. Narrative**

84 *Introduction*

85 The fact that arsenic comes in a variety of forms (with greatly differing toxicological effects) has
86 been known for over 100 years. Arsenic speciation was limited to determination of arsenate and
87 arsenite, until the development of hydride-generation (HG) techniques coupled with atomic
88 absorption detection in the 1970s, which enabled measurement of AsIII, AsV, dimethylarsinic
89 acid (DMA) and monomethylarsenic acid (MMA). Hydride generation is limited, though,
90 because many arsenicals are not amenable to formation of hydrides. Liquid chromatography
91 (LC), and in particular ion chromatography (IC), was more suitable as a method of separating a
92 larger number of species, but success was limited for lack of a suitable detector. This changed
93 with the commercialization of inductively coupled plasma-mass spectrometers (ICP-MS) in the
94 1980s. ICP-MS represented a new detector with superb detection limits which was easy to
95 interface with LC, a powerful tool for separating a large number of compounds in the same
96 solution. The element-specific ICP-MS is especially well suited for detection of arsenic species,
97 as the As heteroatom (m/z 75) allows discrimination from all the other organic material in
98 biological samples. The combination of the two—LC-ICP-MS—enabled a proliferation of
99 arsenic speciation research from the 1990s to the present. A number of other techniques are still
100 used, in part because of the expense associated with ICP-MS. The main weakness of ICP-MS is
101 that it offers no molecular information, meaning that species identification is primarily based on
102 retention time matching with known standards. (*History based on Francesconi & Kuehnelt's*
103 *review on the state of arsenic speciation*)(1)

104 Previous arsenic speciation in chicken tissue

105 Falnoga (2) used LC-HG-AAS to analyzed extracts of freeze-dried liver tissue from hens given
106 feed enriched with arsenic trioxide (30 µg/g), and found AsIII (29 ng/g) and DMA (189 ng/g).
107 Species were separated using a Hamilton PRP-X100 column with 15 mM potassium phosphate
108 pH 6 mobile phase.

109 In 2004, Polatajko and Szpunar (3) reported on the speciation of arsenic in chicken meat, in the
110 context of characterizing the species identification and stability in a freeze-dried chicken meat
111 candidate reference material (source not specified). The authors used methanol-water extraction
112 and LC-ICP-MS detection to find 106 ng/g DMA, 37 ng/g arsenobetaine (AsB, which could
113 come from the fish meal which was a feed ingredient), and ~15 ng/g unknown As (60-65%
114 extraction efficiency reported, 157 ng/g total As in tissue). Species were separated using an AS7
115 column with 0.5 and 50 mM nitric acid mobile phase with 0.5% methanol added.

116 In her Ph.D. dissertation (4), Dr. Tyre Grant developed an LC-ICP-MS method for extraction and
117 analysis of arsenic species in chicken tissue. Several different chromatographic conditions were
118 compared in her work, including PRP-X100 column with 10mM ammonium phosphate + 10mM
119 ammonium nitrate pH 6.3, PRP-X100 column with 20mM ammonium carbonate, and AS7
120 column with 10mM nitric acid. While the PRP-X100 column was good at separating AsIII,
121 DMA, MMA and AsV, she found that roxarsone did not elute from the PRP. Roxarsone was
122 successfully detected using the AS7 column for separation. In tetramethylammonium hydroxide
123 (TMAH) extracts of freeze-dried store-bought chicken liver she found AsIII (0.043 mg/kg), AsV
124 (0.116 mg/kg), MMA (0.106 mg/kg) and roxarsone (0.644 mg/kg), claiming confirmation of
125 roxarsone in edible tissues “for the first time.” Several unidentified peaks presumed to be

126 roxarsone metabolites were also reported. Total arsenic values ranged from 0.02 to 0.057 mg/kg
127 in muscle, and 0.606 to 1.999 mg/kg in liver. The average mass balance (sum of total species
128 found chromatographically divided by total As) was 51%.

129 Another group's report focused on extraction methods for liberating arsenic from chicken meat
130 (5). Using LC-UV-HG-AFS analysis of methanol-water (+ heat + ultrasonics) extracts of freeze-
131 dried store-bought chicken breast, the authors found AsB (48 µg/kg) and Nitarsonone (227 µg/kg)
132 (reported extraction efficiency 80-100%, 270 ng/g total As). Species were separated using a
133 PRP-X100 column with 25mM potassium phosphate pH 5.8 mobile phase.

134 Pizzaro, et al., focused primarily on extraction of arsenic species from various matrices (6).
135 They determined arsenic species in chicken tissue extracts, again as part of a characterization of
136 candidate reference materials. Using methanol/water, they achieved 70-75% extraction
137 efficiency of 0.168 mg/g total arsenic. Of the arsenic extracted, ~15% was AsB, 50% DMA and
138 34% was an unidentified arsenical. Anionic species were separated using a PRP-X100 column
139 with 10mM phosphate pH 6 mobile phase.

140 Using protease digestion and an ultrasonic probe extraction, Sanz et al (7) found AsB (4.7
141 µg/kg), AsIII (2.3 µg/kg) and DMA (133 µg/kg) in freeze-dried chicken muscle tissue. These
142 species represented an 83% mass balance compared to the 169 µg/kg total. Of the matrices
143 examined (rice, fish, chicken muscle, and soil), only chicken required enzymatic digestion to
144 optimize extraction. Speciation was performed using the PRP-X100 anion exchange column
145 with 10 mM ammonium phosphate mobile phase pH 8.5 and 2% methanol and ICP-MS
146 detection. The chicken muscle analyzed in that study was from 70 day old cockerels given an
147 AsIII-enriched diet.

148 For the analytical portion of Study 275.30, we evaluated existing literature methods, combined or
149 modified them as necessary to meet study needs, validated the resulting method and applied it to
150 the analysis of study tissues. Total As levels were much lower in muscle than liver, and muscle
151 was more problematic to extract, so initial efforts focused on speciation of As in liver. This
152 report is limited to the analysis of liver tissues from the control and treated chickens.

153 *Personnel*

154 **Mary Carson** has over 20 years experience developing and validating methods for veterinary
155 drug residue analysis. She is also trained in biochemistry, and is familiar with mass
156 spectrometry and multiple modes of liquid chromatography. However, at the start of this study
157 she had no experience with As speciation. She contacted a colleague at CFSAN, Stephen Capar
158 an expert in trace metal analysis, for advice. He referred her to a member of his branch. **Sean**
159 **Conklin** has several years experience in IC-ICP-MS, with emphasis on As speciation in fruit
160 juice and seafoods. Throughout this project, he was primarily responsible for instrumentation
161 maintenance, calibration and tuning, and took the lead on chromatographic optimization. Mary
162 Carson took the lead in sample preparation and data processing, and as she learned ICP-MS,
163 routine instrument operation. **Karyn Howard** assisted with sample preparation.

164 *Procedures*

165 *Total As Determination*

166 Total As concentration was determined in sliced portions of liver and muscle as described in
167 Attachment 1, "ICP-MS analysis for total arsenic." This method is based on an FDA Elemental
168 Analysis Manual standard analytical procedure for total metals determination. The nitric acid

169 digestion procedure results in a very stable solution which can be stored at room temperature for
170 extended periods of time with no noticeable effect. Instrumental calibrants varied from one set
171 to the next, depending on the expected As concentrations of the samples. Performance of the
172 method was assured by analysis of a NIST standard reference material (SRM 1577c - Bovine
173 Liver, 19.6 ± 1.4 ng As/g certified value, found values in acceptable data sets 15.9 to 22.8 ng/g),
174 as well as by analysis of negative and positive (fortified) controls, with each set. The amount of
175 nitric acid leftover after digestion can affect signal response, so nitric acid concentration in the
176 calibrants needs to match that of the sample digests. Some sets of digests required reanalysis
177 with calibrants prepared in a different concentration of nitric acid.

178 *As Speciation in Liver*

179 Speciation was conducted as described in Attachment 2, “Speciation of arsenic compounds
180 related to roxarsone use in chickens.” The method was developed for this project. Once
181 developed and optimized, method performance was validated concomitantly with sample
182 analysis by inclusion of control and fortified samples with each set. We had standards available
183 for AsIII, AsV, DMA, MMA, Rox, and possible metabolites 3-amino-4-hydroxyphenylarsonic
184 acid (3-Amino) and N-acetyl-4-hydroxyphenylarsonic acid (N-Acetyl). Mixed standards and
185 mixed fortification solutions had equal amounts of As for all seven of these species.

186 Several extraction schemes and chromatographic separations were evaluated during the course of
187 method development. Extraction with water and methanol, similar to that previously reported
188 (3-5) was promising with fortified muscle samples, but could not distinguish a treated muscle
189 sample (total As = 93 ppb) from a control muscle sample (total As = 3 ppb). Subsequent flow
190 injection analysis of the treated muscle extract suggested that extraction efficiency was very low.

191 Extraction with an aqueous alkaline solution, TMAH, (as in Dr. Tyre Grant's dissertation)
192 appeared more efficient. TMAH is a strong base and causes AsIII to be oxidized to AsV. It can
193 also strongly affect ion chromatography (IC) results.


194 IC separation is based on ionic interactions between the analyte, the solid phase, and the mobile
195 phase. These interactions are strongly affected by pH and ionic strength of the mobile phase;
196 ideally samples are in a solvent that is identical to the mobile phase. Many As compounds have
197 a pKa value below 7, making them negatively charged under neutral or basic pH conditions and
198 therefore amenable to separation via anion exchange. A number of ion chromatography methods
199 have been applied to arsenic speciation in a wide variety of samples, but as indicated above, a
200 survey of papers describing arsenic speciation in chicken tissue extracts reveals only 2 different
201 columns have been used for this purpose—the Hamilton PRP-X100 and the Dionex AS7. Other
202 methods and columns have been applied to speciation of As in poultry litter and manure—
203 Jackson and Bertsch (7) used AS14 and AS16 as well as AS7. We confirmed early on the
204 observation by Grant (4) that roxarsone is not easily eluted from the PRP-X100 and therefore
205 that column is not the best choice for this project. However, since the PRP-X100 is often used
206 for determination of AsV (among other species), this column was useful for confirming the
207 presence of AsV in liver extracts.

208 We evaluated various dilution, neutralization, and cleanup steps following extraction of 0.5 g
209 samples with 3 mL 0.625% aqueous TMAH. This concentration of TMAH solubilized most of
210 the liver tissue, with little solid material left. The resultant extract was quite crude.

211 Trichloroacetic acid both neutralized TMAH and precipitated protein, enabling separation on an
212 AS7 column with a nitric acid mobile phase, similar to Grant's method. This method looked

213 promising with fortified samples, though inter-conversion between AsIII and AsV was highly
214 variable, and MMA tended to split between two peaks, MMA1 and Unknown, as shown in the
215 chromatogram of a 10 ppb (mix of 7 compounds) fortified liver extract below:


216 **Figure 1 Chromatography on AS7 with acid mobile phase, 10 ppb stds**



217

218 Unfortunately, the first analysis with an incurred roxarsone liver sample resulted in a very large
219 peak at the front, which was suspected to be unretained roxarsone (or a metabolite closely related
220 to roxarsone) and not AsIII:

221 **Figure 2 AS7 chromatogram of a treated liver extract.**




222

223 Once the AS7 separation was found to be unsuitable, we evaluated a Dionex AS16 column,
224 based on the work of Jackson and Bertsch (8) which showed good separation of several As
225 species in under 10 min. On the AS16, good separation was achieved for most of the species in
226 our standard mix, but two species (N-Acetyl and 3-Amino) co-eluted and the peak shape of AsV
227 was poor.

228 The Dionex AS18 column was the next choice, and it showed a good ability to separate all 7
229 species in the standard mix within 20 min. Fortuitously, the separation was achieved using 45
230 mM tetramethylammonium hydroxide with 1% methanol mobile phase, so it was completely
231 compatible with the TMAH tissue extracts without need for neutralization. In the figure below,
232 the green trace is a 1 ng/g mixed standard in water (equivalent to 20 ppb in tissue), and the black
233 trace is a 20 ppb fortified liver extract.


234 **Figure 3 Standards and fortified liver extract on an AS18 column with 45 mM TMAH**
235 **mobile phase. Retention times: DMA (2.85), AsIII (3.05), MMA (3.3), AsV (7.1), 3-Amino**
236 **(8.9), N-Acetyl (9.9), Rox (18).**




237

238 Injections of incurred tissues showed As-containing compounds that were strongly retained on
239 the column and eluted in subsequent runs. A comparison of the two treated liver extracts below
240 shows the difference between the first treated sample in a set and a subsequent (several runs
241 later) treated sample with interfering late eluters.

242 **Figure 4 AS18 chromatogram of a treated liver extract at the beginning of a set.**



244 **Figure 5 AS18 chromatogram of a treated liver extract later in the set, showing late-eluting**
245 **peaks from a previous injection.**




247 In order to get all peaks, including late eluters, out in the same run, the mobile phase was ramped
248 from 45 mM TMAH to 70 mM for 25 min after roxarsone eluted. It was also necessary to

249 introduce an ultrafiltration cleanup of the extracts, as the unfiltered extracts quickly caused
250 column back pressure to increase, which eventually resulted in column leakage and system
251 shutdown.


252 The chromatograms below illustrate chromatographic results obtained with the final method.
253 The known standards elute in the order DMA, AsIII, MMA (a tight grouping in the first 4 min),
254 AsV, 3-Amino, N-Acetyl, and Roxarson. The first chromatogram is a standard equivalent to 20
255 ppb, the second is a commercial control liver sample, the third is a 20 ppb fortified liver sample,
256 and the bottom one is from a liver extract from treated bird #18. Several of the peaks in this
257 chromatogram do not correspond to known As compounds, including the large peaks at 6 min
258 and after 40 min. We did inject a number of species that were not included in the multi-species
259 standard, including arsenobetaine (RT = 2.8 min), trimethylarsine oxide (RT = 2.8min), 4-
260 arsanilic acid (RT = 3.6min), tetramethylarsonium ion (RT = 3.9 min), arsenocholine (RT = 3.92
261 min), and nitarson (4-nitrobenzenearsonic acid, RT = 7.6 min). None of these compounds is a
262 retention time match with any of the unknowns mentioned above, although arsenobetaine and
263 tetramethylarsine oxide both partially co-elute with DMA. Most of these species are of little
264 concern, because species eluting in less than 4 min from treated bird tissue extracts were present
265 only at very low levels, mostly less than the lower limit of quantification (LLOQ).

266 **Figure 6 Gradient AS18 chromatograph of a standard equivalent to 20 ppb in tissue.**




267

268 **Figure 7 Chromatogram of a liver extract from (commercial) control livers. The AsV peak**
269 **shown is below the LLOQ.**




270

271 **Figure 8 Chromatogram of an extract from 20 ppb fortified liver.**



272

273 **Figure 9 Chromatogram of an extract from a treated liver.**



275 *Analysis Dates*

276 Total As analyses were conducted on dates listed in Table 1.

277 Attempts at speciation analysis began in Nov 2009. However, the speciation method was not
278 finalized until Mar 2010. The reported results come from analyses conducted from 22 Mar 2010
279 through 23 Aug 2010.

280 *Statistical Methods and Calculations*

281 All results are reported as mass elemental As per mass wet tissue.

282 Total As calculations are described in Attachment 1. Values for each sample are generally from
283 replicate analyses with mean and standard deviation reported, though most muscle results are
284 from a single analysis.

285 Species concentrations of individual extracts were determined by comparison to an external
286 standard curve for the known analytes, and by reference to known analytes for unknown
287 concentration (using the slopes of the adjacent known analytes as a response factor for the

288 unknowns). The dilution factor [(9.5+ sample wt)/sample wt] was included in the sequence
289 information, so the numbers reported by the chromatographic data system were in ng As/g tissue.
290 Each extraction set included a Method Blank (water substituted for the tissue, so essentially just
291 the TMAH and water), Method Blank + fortification standards, and quality control (QC) samples
292 consisting of control liver (purchased at a local store and found to have very low levels of As)
293 and fortified control liver samples. If an analyte value in the Method Blank significantly
294 exceeded the LLOQ (or >20% above the lowest calibrant included), the set was repeated as the
295 critical reagent criterion had not been met (see section J.5.a of the liver speciation method SOP).
296 Otherwise, all positive Method Blank values were subtracted from the tissue values for that set.
297 Only AsV ever had a positive Method Blank value, and only in some of the sets. Since As is a
298 naturally occurring element and common environmental contaminant, method blank subtraction
299 is routinely used (9). Reported results are the average of two speciation analyses for each tissue.

300 **Summary of Results**

301 *Total As Determination*

302 Table 2 summarizes the findings for total As determination using microwave assisted nitric acid
303 digestion and ICP-MS analysis. The feed had to be ground in an IKA mill prior to sampling as
304 the original analyses were too variable, which was likely due to heterogeneity of the feed. The
305 entire laboratory portion of both the control and medicated feeds were ground to a fine powder
306 before reanalysis. This laboratory sample had been subsampled from the Control and Test
307 Materials in September 2009. The control feed had 156 ppb of As and the medicated feed had
308 11.3 ppm, which is equivalent to 39.6 ppm of roxarsone.

309 Day 3 muscle samples were not analyzed. All of the livers and all but 3 of the muscles from the
310 untreated birds had total As concentrations below 5 ppb at each time point. Many of the muscle
311 concentrations were below the lowest point of the calibration curves. Total As values for the
312 livers from the treated birds were much higher and varied considerably from bird to bird. They
313 ranged from a high of over 5 ppm for a 0 day withdrawal bird to a low of 275 ppb for one of the
314 5 day withdrawal birds. Values for total As in muscle from treated birds were much lower than
315 liver values. There was approximately 40 times as much total As in liver as in muscle for both 0
316 and 5 day withdrawal treated birds. All muscle values at both 0 and 5 days were below the
317 tolerance value of 0.5 ppm. Some birds at each time point exceeded the liver tolerance of 2 ppm.

318 *Arsenic Speciation in Liver*

319 *Detection and Quantification Limits of the Speciation Method*

320 CVM policy strongly discourages extrapolation beyond upper or lower limits of a calibration
321 curve(10), therefore the lowest level below which AsV was not reliably quantifiable (the lower
322 limit of quantification, LLOQ) was by default defined as the lowest level calibration standard
323 used in the analysis. The lowest calibration mixed standard used was either 0.1 (earlier analyses)
324 or 0.03 ng As/g solvent for each analyte (most analyses used in this report). These standard
325 solutions were equivalent to 2 or 0.6 ppb in tissue.

326 The Limits of Detection (LOD) and LLOQ were statistically estimated by a variety of methods.
327 All methods were based on estimating noise, multiplying it by factors for LOD and LLOQ, and
328 converting the values obtained to equivalent tissue ppb. The initial estimate was derived while
329 still using 0.1 ng/mL (2 ppb tissue equivalent) as the lowest calibrant. A 0.2 ng/mL mixed
330 standard was prepared and analyzed ten (10) times. The standard deviation of the found

331 concentrations was multiplied by $2 \times 1.833 \times \sqrt{1.1}$ for LOD and by 10 for LLOQ. The results are
332 shown in Table 3a.

333 Shortly after this experiment, we discovered that soaking the autosampler vials in 2% Trace
334 Metal Grade Nitric Acid significantly reduced the method blank values for AsV, enabling lower
335 detection and quantification limits. We added a lower calibrant (0.03 ng/mL, equivalent to 0.6
336 ppb in tissue) to the method. With assistance from Dr. Idowu of the CVM Office of New
337 Animal Drug Evaluation, LOD and LLOQ for AsV were estimated for both the older calibration
338 curves (range 0.1 to 100 ng/mL) and the new curves (0.03 to 100 ng/mL) by determining the
339 standard errors of the intercepts (SE Intercept) of the calibration curves used for analysis. A total
340 of 21 curves for AsV were evaluated—10 with the 0.03 ng/mL calibrant and 11 early curves only
341 going as low as 0.1 ng/mL. The evaluation also compared no weighting, 1/x, and 1/x² weighting
342 and verified that 1/x² weighting consistently resulted in the best fit. The SE Intercepts for each
343 1/x² weighted curve were multiplied by 3.28 (LOD) or 10 (LLOQ) and converted to tissue ppb.
344 The average results are shown in Table 3a. With the 0.03 calibrant included, the estimated
345 LLOQ for AsV was less than 0.5 ppb.

346 The data shown in Table 3b provide a point of reference to the accuracy of analyses near these
347 low levels, with 5 replicate analyses of 1 ppb AsV only fortified liver giving an average reading
348 of 1.02 ppb with standard deviation of 0.13 ppb (13% RSD). For AsV, the LOD and LLOQ
349 based on 3× or 10× this standard deviation (0.4 or 1.3 ppb, respectively) are comparable to the
350 LLOQ values of 0.6 or 2 ppb based on lowest point of the calibration curve. For the other
351 analytes, 2 ppb was used as the LLOQ even if the lowest calibrant was 0.03 ng As/g solvent.
352 Roxarsone, 3-Amino and N-Acetyl were not always detectable in the 0.03 ng/g standard, and

353 DMA and MMA showed a pronounced rise in analytical variance when fortification
354 concentrations decreased from 2 ppb to 1 ppb.

355 Speciation Method Evaluation

356 Due to time constraints, the speciation method performance evaluation was conducted using the
357 fortified QC samples prepared and analyzed with each set, along with some additional analyses
358 performed after completion of the study liver samples. The results are summarized in Tables 3a
359 through 3c. AsIII and AsV readily interconvert via oxidation and reduction, and the extraction
360 procedure oxidizes most AsIII to AsV (single analyte fortified sample analyses results, data not
361 shown). This is evident in Table 3b, where accuracy for AsIII is quite low, while accuracy of
362 AsV approaches 200% in all the mixed standard fortified samples. If one sums AsIII and AsV
363 nominal values in the mixed fortified QC samples, the accuracy for AsV determination ranges
364 from 78% to 102% across all concentrations tested, with relative standard deviations (RSDs) of
365 13% to 23%. These results meet CVM criteria (60-110% accuracy and RSD < 20%) for a
366 determinative procedure for residue concentrations less than 100 ppb.

367 The method is clearly incapable of distinguishing AsIII from AsV. The high RSDs (~30%) also
368 mean the method is not suitable for quantifying Rox or its two metabolites 3-Amino and N-
369 Acetyl. The accuracy of the Rox determination was also too low (~70%) for a residue whose
370 concentration is ≥ 100 ppb (CVM criteria 80-110% accuracy and <10% RSD).

371 The performance characteristics for DMA and MMA were generally acceptable for
372 quantification between 2 and 20 ppb based on CVM criteria for a determinative procedure for
373 residue concentrations less than 100 ppb. The analytical variance for MMA is slightly high.


374 Because the determination of iAs was the most crucial measurement, and the method was

375 optimized for this measurement, the less than ideal performance characteristics for DMA, MMA,
376 3-Amino, N-acetyl, and Rox were deemed acceptable. The speciation method should only be
377 used to estimate concentrations of 3-Amino, N-Acetyl, and Rox. MMA and DMA should only
378 be determined above 2 ppb.

379 A liver from a treated bird was selected for replicate analysis to evaluate precision in a “real”
380 sample. Five replicate portions of liver tissue from the same bird were extracted and analyzed.
381 The results are in Table 3c and are consistent with precision results in the fortified samples.

382 AsV is clearly one of the smallest peaks in a liver extract chromatogram from a treated bird. To
383 demonstrate that this peak is truly AsV, we conducted two additional experiments. First, we
384 took portions of an extract from a putative AsV-containing liver, added amounts of AsV standard
385 equivalent to roughly 4x, 2x, and 1x the estimated AsV concentration, and analyzed them. The
386 resulting AsV peak was still a single, symmetrical peak, with area increased as expected, as
387 shown in the chromatograms below:

388 **Figure 10 AsV standard addition to a treated liver extract.**



389

390 AsV is the peak at 7.4 min. The peak at 9.3 min is 3-Amino and shows there is a slight shift
391 forward of retention times over the course of this experiment. The amount of AsV in the original
392 extract was quantified as 4.4 ppb by external calibration and as 5.0 ppb by standard addition.

393 The second experiment was to chromatograph the extracts from a treated liver on an alternate
394 system. In this case, we set up a second IC-ICP-MS instrument with a PRP-X100 column and
395 isocratic mobile phase of pH 9 ammonium carbonate. AsIII elutes before DMA in this system,
396 and AsV elutes at 30 min rather than 7.4 min as on the AS18 system. Some of the extracts from
397 the replicate incurred analyses were divided and analyzed concomitantly on both systems to
398 avoid any concerns of extract stability. The number of injections of Rox-containing extracts had
399 to be limited on the PRP to avoid interference. The results are shown in Table 3d. The amount
400 of AsV found, 4.4 ppb, agrees very well with the 4.2 ppb found on the AS18 analyses in Table
401 3c. Together these two experiments provide strong evidence that the peak we identified as AsV
402 truly is AsV.

403 *Speciation Analysis of Roxarsone-Treated and Control Livers*

404 The speciation results for the study livers are tabulated in Tables 4a and 4b. Values that were
405 below the LLOQ are indicated by greyed-out cells. The only analyte above the LLOQ in any
406 untreated bird at any time point was Rox in Day 3 bird #67, which was due to residual
407 contamination of the system with Rox following speciation analysis of the 3-Nitro 20% Premix
408 (see below). The AsV concentrations are higher in the roxarsone-treated birds than in the control
409 birds.

410 In addition to Rox, 3-Amino was readily detected in all treated samples and a probable peak for
411 N-Acetyl was seen in most. However, a striking feature of the speciation analyses was the

412 presence of several unknown As-containing compounds. Until they can be identified, they are
413 simply labeled by their approximate retention time. Some of these unknowns comprised >10%
414 of the total detectable As in the speciation chromatogram.

415 *Extraction Efficiency and Mass Balance of Arsenic*

416 We compared the sum of all As species found chromatographically with the total As determined
417 in each sample. The sums in general do not add up to the total. The average mass balance
418 (sum/total) was 39%, but the values ranged from 8 to 98%. We investigated this discrepancy by
419 determining total As at each step of the extraction for two incurred liver samples (Day 0 #11 and
420 Day 5 #27). Of the total As concentrations of 1760 and 2940 µg/kg, 93% and 108% were
421 extracted by the TMAH, respectively, indicating nearly quantitative extraction of As compounds
422 from the tissue. The ultrafiltrates contained 67% and 69% of the total tissue As. This loss
423 through ultrafiltration was consistent with the loss seen for roxarsone recovery in fortified
424 experiments (Table 3a). However, the species for these two samples, including estimated
425 concentrations of the unknown As-containing metabolites, added to only 53% and 32%,
426 suggesting additional loss of As analytes during chromatography. To verify this possibility, we
427 removed resin from a used column and guard cartridge, determined the total As concentration by
428 microwave-assisted nitric acid digestion, and compared those values with one obtained from new
429 guard cartridge resin. Used resin had approximately 2 mg/kg As at the inlet end of the column,
430 and 0.5 mg/kg at the outlet end of the same column, while new guard cartridge resin had 0.02
431 mg/kg As. Clearly, some As compound(s) bound very tightly to the AS18 column under the
432 conditions we used.

433 The low and variable mass balances suggest that we do not have a full accounting of all the As
434 species that may exist in livers from roxarsone-treated birds, and that the numbers we are
435 reporting should be treated as minimum numbers. Actual concentrations may be higher.

436 *Speciation in Control Feed, Medicated Feed, and Type A Medicated Articles*

437 When inorganic arsenic in livers from chickens treated with roxarsone medicated feed was
438 found, we wanted to be sure it was not because the medicated feed itself contained inorganic
439 arsenic. Control feed, the medicated feed, and the premix used to make the medicated feed were
440 tested for the two inorganic arsenic species and several organic arsenic species.

441 An *ad hoc* method based on the instrumental procedure used for the chicken tissue samples was
442 used. Analysis consisted of mixing a small portion of each feed and premix with water or
443 TMAH to dissolve the drug and any water soluble contaminants, and then analyze the extract for
444 arsenic species by IC-ICP-MS using the same chromatographic conditions employed for the liver
445 speciation. The water extraction has the potential to distinguish AsIII from AsV, while the
446 TMAH extraction converts most AsIII to AsV. The *ad hoc* method is not the same as the official
447 feed analysis procedure. The official feed method will not distinguish inorganic arsenic from
448 roxarsone.

449 The results are shown in Table 5a. In the first analysis of the original 20% roxarsone premix, a
450 random sample was taken from the top of the 50lb premix bag and analyzed. The premix
451 contained 1 part per million AsIII, 12 parts per million AsV, and 38,000 parts per million
452 roxarsone (equivalent to 14% roxarsone) using the water extraction. Based on the total arsenic in
453 the premix, 0.03% was inorganic arsenic and 99.97% was roxarsone or other organic arsenic.

454 The control feed had 0.024 parts per million AsIII and 0.032 parts per million AsV. The premix

455 is diluted 4000-fold with control feed to make medicated feed. The amount of inorganic As
456 contributed by the Rox premix to the medicated feed is in the vicinity of 3-4 ppb, or $\leq 10\%$
457 above control feed levels. We therefore find it unlikely that the AsV found in the treated bird
458 livers is a result of direct inorganic As contamination of the medicated feed.

459 It was decided to further investigate the inorganic arsenic concentrations in organic arsenical
460 premixes. A penicillin premix (a presumed negative control), a premix containing 50%
461 nitarsonsone, and a second lot of a 20% roxarsone premix were purchased. Greater care was taken
462 to ensure that representative samples of the premixes were obtained using a “core sampling
463 probe” that is approximately one (1) inch in diameter and about three (3) feet long. The probe
464 was pushed into the bags of medicated premix until it reached the bottom of the bag and was
465 rotated several times to enable the test article to fill the probe. The probe’s outer tube was
466 rotated to close off the probe and the collected materials were poured into a sample collection
467 bag. This process was repeated about 12-15 times by sampling in a clockwise rotation within the
468 bag, thereby collecting material from all parts of the bags of each of the test articles. During the
469 sampling process, it was noted that despite the fact that these four (4) preparations were
470 “supposed” to have about the similar proportions of inert ingredients, there was a marked
471 difference in the texture of the various products.

472 The results are shown in Table 5b. The amount of inorganic arsenic in the penicillin premix was
473 very low, in the parts per billion range. The results for the original roxarsone premix were
474 consistent with the initial analysis. However, the second lot of roxarsone premix had a much
475 higher percentage of inorganic arsenic--0.4%. The percent of inorganic arsenic in the nitarsonsone
476 premix was also high, but within the manufacturing specifications for this compound.

477 **Conclusions**

478 The CVM study wherein chickens were raised on feed medicated with roxarsone according to
479 label directions along with chickens raised on unmedicated feed was specifically designed to
480 determine which arsenic species are present in chicken tissue as a direct result of roxarsone
481 treatment. The liver tissues collected from these birds were found to contain higher levels of
482 total arsenic than had been reported in store-bought tissues in previous work (4), and more
483 arsenic species were found in our samples than in previous work. This was at least partially due
484 because we included more species in our multi-species standard than in previous work, *i.e.*, we
485 found more peaks because we were looking for them. While some inorganic arsenic had been
486 reported in chicken tissues (4) (7), the source of the arsenic contamination cannot be ascertained
487 since the tissues were purchased from stores and the birds' treatment prior to the point of
488 purchase was unknown to the researchers. By comparing treated and untreated birds raised
489 under carefully controlled conditions, it should be possible to determine whether roxarsone
490 treatment represents a food safety risk based on increased inorganic arsenic concentration in the
491 treated birds.

492 **Signatures**

493

494 Mary C. Carson

495 Analyst Mary C. Carson (CVM)

1/24/2011

Date

496

497 Sean Conklin

498 Analyst Sean D. Conklin (CFSAN)

1/24/2011

Date

499

500 Philip J. Kijak

501 Philip J. Kijak, Director, Division of Residue Chemistry (CVM)

01/24/11

Date

502

503 Stephen G. Capar

504 Stephen G. Capar, Chief, Chemical Contaminants Branch (CFSAN)

1/24/2011

Date

505

III. Tables

506 Table 1. Dates for digestion and total As determination in study samples

Digestion Date	Analysis Date	Sample Set	Comments
14 Oct 2009	15 Oct 2009	Feed and standards	Sequence SDC1015b
12 Nov 2009	12 Nov 2009	Day 0 Breast	Prewighed breast, instrumental problems
12 Nov 2009	13 Nov 2009	Day 0 Breast	Prewighed breast, reanalysis
12 Nov 2009	13 Nov 2009	Day 0 Liver Ground Feed	Prewighed liver
30 Nov 2009	2 Dec 2009	Day 5 Liver	Prewighed liver
2 Dec 2009	3 Dec 2009	Day 5 Breast	Prewighed breast
12 Nov 2009	4 Dec 2009	Day 0 Breast	Reanalysis with optimized tune file
12 Nov 2009	4 Dec 2009	Day 0 Liver Ground Feed	Reanalysis with optimized tune file
2 Dec 2009	4 Dec 2009	Day 5 Breast	Reanalysis with optimized tune file
9 Apr 2010	9 Apr 2010	Selected Day 0 and Day 5 liver repeats	Freshly weighed samples
12 Apr 2010	12 Apr 2010	Selected Day 5 liver re-repeats	Freshly weighed samples
1 Jun 2010 and 2 Jun 2010	2 Jun 2010	Day 3 Liver	Freshly weighed, duplicate or triplicate
4 Oct through 12 Oct 2010	4 Oct through 1 Nov 2010	Day 0 and Day 5 Liver	Freshly weighed, triplicates, reanalysis with calibrants in different % nitric

507 Table 2. Total Arsenic Concentrations in Feed, Liver, and Muscle Samples

Bird #	Treatment	Withdrawal time (days)	Total As (ng/g wet wt)	n	Total As (ng/g wet wt)	n
Feed	Control	N/A	156	3		
Feed	Roxarsone	N/A	11300	3		
			LIVER		MUSCLE	
66	Control	0	1.3±0.2	3	5.0	1
69	Control	0	1.8±0.2	3	3.9	1
75	Control	0	3.1±0.5	4	3.6	1
77	Control	0	4.6±1.3	3	3.5	1
82	Control	0	3.2±0.6	3	2.7±0.1*	3
84	Control	0	4.4±0.5	4	4.7	1
97	Control	0	3.5±0.2	4	10.7	1
99	Control	0	4.9±0.8	4	5.4	1
11	Roxarsone	0	1760±60	4	72.9±17.1	3
12	Roxarsone	0	2920±260	4	55.4	1
21	Roxarsone	0	5430±220	5	68.9	1
26	Roxarsone	0	4200±580	4	72.7	1
28	Roxarsone	0	4140±240	4	126.	1
36	Roxarsone	0	1950±330	4	64.4	1
40	Roxarsone	0	899±106	4	54.2	1
46	Roxarsone	0	4300±120	4	78.7	1
48	Roxarsone	0	2920±410	4	80.7	1

Bird #	Treatment	Withdrawal time (days)	Total As (ng/g wet wt)	n	Total As (ng/g wet wt)	n
52	Roxarsone	0	2460±360	4	81.6	1
53	Roxarsone	0	3560±650	4	91.8	1
54	Roxarsone	0	4410±150	4	71.2	1
56	Roxarsone	0	652±35	5	54.1	1
58	Roxarsone	0	2180±420	5	61.1	1
59	Roxarsone	0	1360±170	4	57.7	1
60	Roxarsone	0	2240±230	4	38.6	1
67	Control	3	3.7±0.5	6		
68	Control	3	4.7±1.1	3		
74	Control	3	4.3±0.3	3		
78	Control	3	2.3±0.2	3		
79	Control	3	3.5±0.1	3		
83	Control	3	3.4±0.2	3		
91	Control	3	4.5±0.2	3		
100	Control	3	2.9±0.1	3		
17	Roxarsone	2	3170±110‡	3		
34	Roxarsone	2	2250±210	3		
1	Roxarsone	3	4590±410	3		
2	Roxarsone	3	1600±760	3		

Bird #	Treatment	Withdrawal time (days)	Total As (ng/g wet wt)	n	Total As (ng/g wet wt)	n
5	Roxarsone	3	988±42	2		
9	Roxarsone	3	2510±290	2		
13	Roxarsone	3	3090±850	3		
18	Roxarsone	3	843±24	3		
23	Roxarsone	3	3270±150	2		
30	Roxarsone	3	2520±940	3		
31	Roxarsone	3	941±135	3		
33	Roxarsone	3	1720±350	3		
38	Roxarsone	3	1070±170	2		
47	Roxarsone	3	1110±580	3		
50	Roxarsone	3	469±38	3		
61	Roxarsone	3	2010±1120	3		
62	Roxarsone	3	2640±15	2		
65	Roxarsone	3	1930±290	2		
70	Control	5	4.3±0.6	4	3.2	1
71	Control	5	4.6±1.2	4	3.4	1
72	Control	5	2.8±0.3	3	2.5	1
73	Control	5	2.8±0.4	4	2.5	1
81	Control	5	3.6±0.7	4	3.1	1
86	Control	5	3.2±2.3	3	3.2	1
87	Control	5	3.0±0.3	4	3.3	1

Bird #	Treatment	Withdrawal time (days)	Total As (ng/g wet wt)	n	Total As (ng/g wet wt)	n
89	Control	5	2.5±1.2	4	1.8	1
92	Control	5	3.0±0.8	4	4.9	1
94	Control	5	2.7±1.0	4	2.9	1
98	Control	5	3.3±0.5	4	3.6	1
3	Roxarsone	5	1060±310	4	30.6	1
4	Roxarsone	5	1150±170	4	31.5	1
6	Roxarsone	5	655±190	4	24.7	1
7	Roxarsone	5	1910±540	4	44.7	1
8	Roxarsone	5	1590±230	4	23.1	1
14	Roxarsone	5	691±416	6	19.1	1
15	Roxarsone	5	427±44	5	23.2	1
16	Roxarsone	5	480±45	6	17.4	1
19	Roxarsone	5	301±12	4	13.9	1
24	Roxarsone	5	1640±180	6	39.9	1
27	Roxarsone	5	2940±360	6	46.6	1
29	Roxarsone	5	581±61	8	21.5	1
32	Roxarsone	5	1220±210	7	19.8	1
37	Roxarsone	5	2170±220	5	48.4	1
39	Roxarsone	5	275±22	6	21.4	1
44	Roxarsone	5	633±136	4	14.0	1
49	Roxarsone	5	388±53	4	39.0	1

Bird #	Treatment	Withdrawal time (days)	Total As (ng/g wet wt)	n	Total As (ng/g wet wt)	n
51	Roxarsone	5	907±160	4	25.0	1
55	Roxarsone	5	499±23	6	23.9	1
63	Roxarsone	5	1030±160	4	21.6	1
64	Roxarsone	5	1745±52	4	30.5	1

508 *Two replicates below that set's lowest calibrant—this value was censored in the statistical
509 analysis.

510 ‡One replicate above that set's highest calibrant

511 Table 3a. Estimates of LOD and LLOQ by various methods. Values are in tissue ppb.*

Analyte	DMA	AsIII	MMA	AsV	3-Amino	N-Acetyl	Rox
Method	10 analyses of 0.2 ng/mL standard; lowest calibrant 0.1 ng/mL						
LOD	0.8	1.6	0.6	1.1	1.5	0.4	1.4
LLOQ	2.0	4.0	1.6	3.0	3.9	1.0	3.5
Method	SE Intercept of 11 calibration curves, range 0.1 to 100 ng/mL						
LOD				0.7			
LLOQ				2.1			
Method	SE Intercept of 10 calibration curves, range 0.03 to 100 ng/mL						
LOD				0.15			
LLOQ				0.45			

512 *The dilution factor from equivalent tissue ppb ($\mu\text{g As/kg tissue}$) to standard ng/mL is 20.

513 Table 3b. Validation Results for Arsenic Speciation Method—Fortified Liver Accuracy and
514 Precision

Fortification Level		DMA	AsIII*	MMA	AsV*	3-Amino	N-Acetyl	Rox
2000 ppb Rox	Average Found ppb	0.3			0.9	2.7		1490
n = 16	Accuracy							74%
	RSD	97%			103%	173%		28%
20 ppb mix	Average Found ppb	21.3	2.1	16.4	38.7	16.6	14.9	14.5
n = 12	Accuracy	106%	10%	82%	194%	83%	75%	72%
	RSD	10%	275%	23%	15%	33%	33%	33%
4 ppb mix	Average Found ppb	4.6	0.0	3.77	8.0	4.1	3.8	3.4
n = 7	Accuracy	116%	0%	94%	201%	103%	95%	86%
	RSD	9%		17%	15%	32%	30%	28%
2 ppb mix	Average Found ppb	2.3	0.04	1.7	3.6	1.7	1.6	1.4
n = 9	Accuracy	117%	2%	86%	180%	86%	82%	70%
	RSD	7%	300%	21%	18%	32%	29%	33%
1 ppb mix	Average Found ppb	1.2	0.00	0.8	1.6	1.1	0.8	0.6
n = 10	Accuracy	119%	0%	85%	157%	113%	78%	62%
	RSD	23%		32%	23%	25%	42%	86%

Fortification Level		DMA	AsIII*	MMA	AsV*	3-Amino	N-Acetyl	Rox
1 ppb AsV	Average Found ppb	0.3			1.0			
n = 5	Accuracy				102%			
	RSD	16%			13%			

515 *The extraction procedure oxidizes most AsIII to AsV.

516 Table 3c. Validation Results for Arsenic Speciation Method—Precision with an Incurred Roxarsone Liver Sample. Cells that are
 517 grey indicate a censored value below the LLOQ. Concentration units are ppb or µg As/kg wet weight liver.

Set	Sample	DMA	AsIII	MMA	Unk 3.5	Unk 4.5	Unk 5.6	AsV	3-Amino	N-Acetyl	Unk 13	Rox	Unk 21	Unk 32	Unk 36
10H23 #10	Treated #11	1.2	1.1	1.8	5.0	1.7	62	4.5	65	9.7	2.1	543	8.0	8.4	121
10H23 #11	Treated #11	1.3	1.0	1.4	5.4	1.2	63	3.2	44	7.9	0.0	377	4.9	7.0	81
10H23 #12	Treated #11	1.3	1.3	1.8	6.0	1.9	60	4.7	71	10.1	2.2	605	7.6	10.9	122
10H23 #13	Treated #11	1.1	1.2	1.4	3.7	1.8	42	4.6	41	6.9	2.0	328	6.1	8.4	96
10H23 #14	Treated #11	1.2	0.9	1.6	4.7	1.8	51	4.2	49	7.4	2.4	375	5.4	9.3	98
	Average	1.2	1.0	1.6	4.9	1.7	56	4.2	54	8.4	1.7	446	6.4	8.8	104
	Std Dev	0.1	0.1	0.2	0.9	0.3	8.8	0.6	13	1.4	1.0	121	1.3	1.4	17.8
	RSD	8%	12%	13%	18%	17%	16%	14%	24%	17%	57%	27%	21%	16%	17%

518 Table 3d. Validation Results for Arsenic Speciation Method—Verification of AsV in Incurred
 519 Sample by Alternate Chromatography. Concentration units are ppb or µg As/kg wet weight
 520 liver.

521 Column: PRP-X100 anion exchange column

522 Mobile phase: 20 mM Ammonium Carbonate pH 9 (adjusted with ammonium hydroxide)

523 Conditions: 1 mL/min flow rate, 100 µL injections, 40 min run time

524 Instrument: Agilent 7500 ICP-MS, operated in helium collision mode

525 Analysis within 24 hrs of AS18 column analysis.

526

	AsIII	DMA	MMA	AsV
19Aug #12 Bird #11 liver extract	0	0.95	1.90	3.22
19Aug #13 Bird #11 liver extract	0	1.17	2.14	5.25
23Aug #10 Bird #11 liver extract	0	0.93	2.04	3.72
23Aug #10 Bird #11 liver extract	0	0.93	2.04	4.27
23Aug #11 Bird #11 liver extract	0	0.83	2.28	3.93
23Aug #11 Bird #11 liver extract	0	1.04	2.07	4.76
23Aug #12 Bird #11 liver extract	0	0.96	2.11	4.81
23Aug #14 Bird #11 liver extract	0	0.60	2.19	5.17
Average	0	0.92	2.10	4.39
RSD		18%	5%	17%

527 Table 4a. As Species in Liver Samples from Birds Treated with Roxarsone Medicated Feed.

528 Concentrations are in ppb ($\mu\text{g As per kg wet weight liver}$). Results are the average of duplicate analyses.

529 See the text for a description of LLOQ estimation. Cells that are grey indicate a censored value below the LLOQ.

530 The extraction solvent (tetramethylammonium hydroxide in water) oxidizes most AsIII to AsV, so AsIII is rarely found.

531 Compounds are listed in chromatographic elution order.

532 Unknown As compound concentrations were estimated by using the slope of the adjoining known compounds as a response factor.

With- drawal Day	Bird	DMA	AsIII	MMA	Unk 3.5	Unk 4.5	Unk 5.5	AsV	3- Amino	N- Acetyl	Unk 13	Rox	Unk 21	Unk 32	Unk 36
0	#11	0.4	0.8	1.7	7	2	54	8.2	97.7	12.8	1	588	19	n.d.†	137
0	#12	1.0	0.3	2.0	10	1	71	21.5	245	8.4	3	1640	37	n.d.	216
0	#21	0.0	0.5	0.2	3	1	12	18.1	73.8	7.1	2	333	24	n.d.	90
0	#26	1.2	0.4	2.3	12	2	93	19.3	215	15.4	2	1920	29	n.d.	234
0	#28	0.8	0.4	1.7	10	2	155	8.4	66.0	15.9	2	2810	13	n.d.	107
0	#36	0.8	0.4	1.6	5	2	84	3.8	39.6	9.5	1	1430	7	n.d.	66
0	#40	0.9	0.1	0.8	2	1	70	0.1	3.7	0.5	0	284	0	n.d.	58
0	#46	1.9	0.0	1.3	5	1	123	5.2	48.4	12.5	2	3560	8	n.d.	222
0	#48	1.6	0.0	1.3	4	1	76	4.2	33.3	5.4	0	2040	6	n.d.	164
0	#52	1.8	1.2	1.5	3	2	66	5.4	31.9	8.0	1	865	9	n.d.	204

With- drawal Day	Bird	DMA	AsIII	MMA	Unk 3.5	Unk 4.5	Unk 5.5	AsV	3- Amino	N- Acetyl	Unk 13	Rox	Unk 21	Unk 32	Unk 36
0	#53	1.7	0.0	1.6	12	2	175	11.4	118	15.5	3	2560	13	n.d.	240
0	#54	2.1	0.3	1.0	3	1	70	4.0	23.4	4.8	1	1290	2	n.d.	125
0	#56	0.4	0.4	0.5	2	1	15	0.3	8.1	1.4	0	43.1	0	n.d.	19
0	#58	1.7	0.5	1.1	3	2	44	2.2	19.8	3.1	0	367	4	n.d.	57
0	#59	1.3	0.0	1.0	2	1	51	2.3	16.3	3.0	0	322	2	n.d.	41
0	#60	1.3	1.6	1.0	6	3	32	6.8	52.9	9.2	0	306	13	n.d.	88
2	#17	1.2	0.5	0.9	4	1	125	2.5	18.4	10.5	3	1520	5	0	210
2	#34	1.6	0.5	0.9	1	1	51	1.5	8.4	3.9	0	838	0	0	83
3	#1	1.3	0.0	0.1	11	2	105	3.2	103	34.6	1	2150	13	6.5	246
3	#2	0.6	0.0	0.0	2	1	31	0.0	15.9	4.5	0	270	3	0	67
3	#5	1.0	0.0	0.0	1	1	77	0.0	1.3	0.5	0	112	1	0	45
3	#9	1.2	0.0	0.0	2	1	70	0.1	7.2	2.2	0	585	2	0	77
3	#13	1.0	0.0	0.2	5	1	79	2.2	49.3	5.2	0	1490	5	6.5	215
3	#18	1.5	0.5	0.4	1	1	50	0.3	1.6	0.6	0	99.2	0	0	35

With- drawal Day	Bird	DMA	AsIII	MMA	Unk 3.5	Unk 4.5	Unk 5.5	AsV	3- Amino	N- Acetyl	Unk 13	Rox	Unk 21	Unk 32	Unk 36
3	#23	1.5	0.3	0.8	5	2	145	3.2	34.8	8.3	2	1670	6	7.7	197
3	#30	0.8	0.0	0.3	5	1	72	1.2	34.8	7.8	2	771	4	7.2	92
3	#31	1.2	0.2	0.1	1	1	34	0.1	1.9	0.5	0	90.9	0	0	20
3	#33	1.0	0.2	0.2	0	1	47	0.3	1.7	0.5	0	299	0.8	0	42
3	#38	1.1	0.3	0.3	1	1	54	0.3	2.3	0.8	0	172	0	0	36
3	#47	1.2	0.5	0.4	1	1	49	0.3	2.0	0.6	0	139	0	0	33
3	#50	1.0	0.4	0.3	1	1	53	0.3	1.4	0.3	0	96.1	0	0	29
3	#61	1.1	0.5	0.8	1	1	53	0.5	1.1	0.2	0	66.9	0	0	35
3	#62	0.7	0.1	1.2	10	0	141	4.5	97.8	11.8	4	1350	10	9.0	217
3	#65	1.2	0.5	0.4	1	1	71	0.5	4.6	1.8	0	310	0	0	60
5	#3	1.1	0.6	0.2	1	1	28	0	1.1	0.9	0	72.0	0	1	25
5	#4	0.9	0.4	0.7	3	1	54	1.7	41.2	8.4	1	630	2	16	105
5	#6	0.0	0.0	0.6	9	1	36	0.5	4.9	12.3	1	32.2	2	17	29
5	#7	0.2	0.3	1.5	21	1	63	5.2	117	31.7	5	443	17	55	187
5	#8	0.7	0.5	0.3	2	1	25	0	7.2	2.4	0	345	0	19	41

With- drawal Day	Bird	DMA	AsIII	MMA	Unk 3.5	Unk 4.5	Unk 5.5	AsV	3- Amino	N- Acetyl	Unk 13	Rox	Unk 21	Unk 32	Unk 36
5	#14	0.9	0.1	0.2	0	1	30	0	0.6	0.2	0	28.6	0	0	19
5	#15	1.2	0.2	0.3	1	1	37	0.1	0.6	0.7	0	23.0	1	0	20
5	#16	0.8	0.3	0.6	10	1	57	0	3.6	4.0	0	38.8	0	9	20
5	#19	0.8	0.1	0.5	0	1	20	0.3	0.6	0.6	0	20.8	1	1	20
5	#24	1.2	0.6	0.4	1	2	17	0.8	10.4	2.6	0	316	1	5	75
5	#27	0.0	0.0	0.2	9	1	50	9.1	93.0	15.8	3	592	14	13	134
5	#29	1.0	0.2	0.7	3	1	102	1.3	7.3	4.2	1	175	2	8	74
5	#32	1.2	0.2	0.4	1	1	39	0.3	2.4	1.5	0	287	1	0	47
5	#37 (extra)	1.5	0.3	0.5	1	1	52	0.9	13.3	6.6	0	1010	2	5	103
5	#39 (extra)	0.8	0.1	0.3	0	1	32	0	0.0	0.0	0	10.3	0	4	20
5	#44	1.0	0.5	0.5	0	1	27	0.3	2.2	0.7	0	100	0	0	22
5	#49 (extra)	3.6	1.1	0.7	1	2	49	0	1.1	0.5	0	41.9	1	0	20
5	#51	0.7	0.4	0.5	1	1	31	0.9	8.5	1.8	1	162	1	5	46
5	#55 (extra)	1.4	0.6	0.5	1	1	56	0	1.1	0.6	0	28.8	1	4	27

With- drawal Day	Bird	DMA	AsIII	MMA	Unk 3.5	Unk 4.5	Unk 5.5	AsV	3- Amino	N- Acetyl	Unk 13	Rox	Unk 21	Unk 32	Unk 36
5	#63	0.8	0.5	0.5	1	1	27	0.6	9.6	2.0	0	196	1	1	46
5	#64	1.1	0.8	0.6	1	1	35	1.0	14.1	2.4	0	313	1	3	60

533 †Not done. This unknown compound was not added to the list of unknowns in the processing method until later.

534 Table 4b. As Species in Liver Samples from Untreated Birds. Conditions the same as in Table 4a.

With drawal Day	Bird	DM A	AsIII	MM A	Unk 3.5	Unk 4.5	Unk 5.5	AsV *	3- Amino	N- Acetyl	Unk 13	Rox	Unk 21	Unk 32	Unk 36
0	#66	0	0	0.0	0	0	0	0	0	0	0	0	0	n.d.	0
0	#69	0.1	0	0.2	0	0	0	0.0	0	0	0	0	0	n.d.	0
0	#75	0.4	0	0.4	0	1	0	1.1	0	0	0	0	0	n.d.	0
0	#77	0.3	0	0.4	0	0	0	0	0	0	0	0	0	n.d.	0
0	#82	0.2	0	0.3	0	0	0	0	0	0	0	0	0	n.d.	0
0	#84	0.5	0	0.6	0	1	0	0.5	0	0	0	0	0	n.d.	0
0	#97	0.3	0	0.5	1	1	0	0	0	0	0	0	0	n.d.	0
0	#99	0.6	0	0.6	0	1	0	0.2	0	0	0	0	0	n.d.	0
3	#67	0.9	0	0	0	0	0	0	0	0	0	2.6	0	0	0
3	#68	0.8	0	0	0	0	0	0	0	0	0	1.0	0	0	0
3	#74	0.8	0	0	0	1	0	0	0	0	0	1.8	0	0	0
3	#78	0.8	0.1	0.6	0	0	0	0.1	0	0	0	0	0	0	0
3	#79	0.6	0	0.5	0	0	0	0.1	0	0	0	0	0	0	0

With drawal Day	Bird	DM A	AsIII	MM A	Unk 3.5	Unk 4.5	Unk 5.5	AsV *	3- Amino	N- Acetyl	Unk 13	Rox	Unk 21	Unk 32	Unk 36
3	#83	0.8	0.0	0.5	0	1	0	0	0	0	0	0	0	0	0
3	#91	0.9	0.3	1.0	0	1	0	0	0	0	0	0	0	0	0
3	#100	0.7	0.1	0.6	0	0	0	0	0	0	0	0	0	0	0
5	#70	0.6	0.1	0.5	0	1	0	0	0	0	0	0	0	0	0
5	#71	0.4	0.3	0.6	0	1	0	0	0	0	0	0	0	0	0
5	#72	0.3	0.1	0.1	0	0	0	0	0	0	0	0	0	0	0
5	#73	0.4	0.0	0.4	0	0	0	0	0	0	0	0	0	0	0
5	#81	0.9	0.1	0.7	0	1	0	0	0	0	0	0	0	0	0
5	#86	0.7	0.0	0.4	0	1	0	0	0	0	0	0	0	0	0
5	#87	0.4	0.0	0.1	0	0	0	0	0	0	0	0	0	0	0
5	#89	0.5	0.0	1.3	0	1	0	0	0	0	0	0	0	0	0
5	#92	1.0	0.5	0.7	0	0	0	0	0	0	0	1.3	0	0	0
5	#94	0.6	0.3	0.5	0	0	0	0	0	0	0	0	0	0	0
5	#98	0.8	0.3	0.8	0	0	0	0	0	0	0	0	0	0	0

536 Table 5a. Initial Speciation Analysis Of Control Feed, Medicated Feed, And Premix. **Concentration is in ppm (mg As per kg).**

Sample	DMA	AsIII	MMA	Unk 3.5	Unk 4.5	Unk 5.6	AsV	3-Amino	N-Acetyl	Unk 13	Rox*	Unk 21	Unk 32	Unk 36
Control feed (water extract)	<LLOQ	0.024	<LOD	<LOD	<LOD	<LOD	0.0316	<LOD	<LOD	<LOD	0.00364	<LOD	<LOD	<LOD
Control feed (TMAH extract)	<LLOQ	<LOD	<LLOQ	<LOD	<LOD	<LOD	0.0381	<LOD	<LOD	<LOD	0.0154	<LOD	<LOD	<LOD
Medicated feed (water extract)	<LLOQ	0.0228	<LOD	<LOD	<LOD	<LOD	0.0362	<LLOQ	<LOD	<LLOQ	11.3	0.003	<LOD	0.15
Medicated feed (TMAH extract)	<LLOQ	<LLOQ	<LLOQ	<LOD	<LLOQ	<LOD	0.0470	0.0098	<LLOQ	<LLOQ	8.00	0.003	0.046	0.051
Premix (diluted) (water extract)	<LOD	1.02	<LOD	<LOD	0.44	<LLOQ	12.2	2.03	1.64	0.110	38100	22	810	53
Premix (diluted) (TMAH extract)	<LOD	0.731	<LOD	<LLOQ	0.45	<LLOQ	12.4	2.07	1.63	0.182	38200	12	460	58
Premix (direct) (water extract)	<LOD	1.00	<LOD	0.094	0.11	0.18	11.2	1.72	1.45	<LOD	Off scale	11	670	40
Premix (direct) (TMAH extract)	<LOD	0.770	<LOD	0.074	0.13	0.36	12.5	1.82	1.22	<LOD	Off scale	15	420	52

537 *The conversion from Rox concentration as ppm As to percent roxarsone is: $\text{ppm}/10^6 \times 263/75 \times 100$

538 Table 5b. Analysis of Type A Medicated Articles for iAs using Water Extraction

Product	API	Number of assays	API Conc.		AsIII Conc.		AsV Conc.	
			NADA spec.	Found	Found	NADA spec.**	Found	NADA spec.**
3-Nitro 20 (Bag 1/test 1)	Rox	3	18 – 22 %	13.4%	1 ppm	≤50 ppm	12 ppm	≤100 ppm
3-Nitro 20 (Bag 1/test 2)	Rox	5	18 – 22 %	21.1%	1 ppm	≤50 ppm	22 ppm	≤100 ppm
3-Nitro 20 (Bag 2)	Rox	5	18 – 22 %	21.3%	31 ppm	≤50 ppm	865 ppm	≤100 ppm
Histostat 50	Nitarstone	5	45 – 50 %	47.5%	686 ppm	≤5000 ppm	207 ppm	≤5000 ppm
Penicillin 100*	Penicillin G	5	NA	NA	-	NA	30 ppb	NA

539 *contains similar inactive carriers/extenders as arsenical products

540 **Roxarsone API impurity limit is 0.025% for AsIII and 0.05% for AsV. Nitarstone has single
 541 specification for limit (1% of API) of inorganic arsenic which would represent the sum of As
 542 (III) and As (V).

543	IV. List of Excel Files Contributing to Tables 2, 3, 4, and 5	
544	Table 2	0- 3- and 5-day liver totals- adjusted for NR audit.xls (muscle and feed)
545		All liver totals results.xls (liver)
546		Contributing:
547		Day 3 liver digest worksheet.xls
548		09I02 5-day birdsTOTALS As worksheet.xls
549		09I09 0-Day birds TOTAL As worksheet.xls
550		10D09 and 10D12 treated liver MWdigest TOTAL As checks.xls
551		10J04 digestion form 434-142.xls
552		10J04 digestion form.xls
553		10J05 digestion form 434-143.xls
554		10J05 digestion form.xls
555		10J06 digestion form 434-146 HNO3.xls
556		10J06 digestion form.xls
557		10J06pm digestion form 434-147.xls
558		10J06pm digestion form.xls
559		10J07 digestion form 434-148.xls
560		10J12am digestion form 434-150.xls
561		10J12pm digestion form 434-151.xls
562		10J13 digestion form 434-152.xls
563		434-143 -146 -151 -152 Totals with GOOD SRM.xls
564	Table 3a	10C11 replicates QUANTTAB.xls
565		arsenic_calibration_data_Table_Remi-MCC comments.xls

566 Contributing: This file summarizes all valid calibration data acquired from 22
567 Mar to 20 Sep 2010, so pulls data from all QUANTTAB.xls files listed below for
568 Tables 3b, 3c, 4, and 5, as well as these additional sets:

569 10C30QUANTTAB.xls (tissue results not used because reanalyzed later with
570 lower calibrant)

571 10D20QUANTTAB.xls (tissue results not used because reanalyzed later with
572 lower calibrant)

573 10D22QUANTTAB.xls (tissue results not used because reanalyzed later with
574 lower calibrant)

575 10H19QUANTTAB.xls (Method Blank high for this set—but calibration ok)

576 10I20QUANTTAB.xls

577 Table 3b MethodEval + ASDL and LOD + orthogonal-MBcorrected.xls, worksheet
578 Fortification recoveries

579 Contributing: All files contributing to Table 4, plus

580 10E20QUANTTAB.xls

581 10F15QUANTTAB.xls

582 10H23QUANTTAB.xls

583 Table 3c MethodEval + ASDL and LOD + orthogonal-MBcorrected.xls, worksheet
584 Incurred replicates

585 Contributing:

586 10H23QUANTTAB.xls

587 Table 3d MethodEval + ASDL and LOD + orthogonal-MBcorrected.xls, worksheet
588 “Orthogonal” results

589 Contributing:

590 10H19 QUANTTAB PRP-X100.xls

591 10H24 QUANTTAB treated extracts on PRP-X100.xls

592 Table 4 Chicken liver speciation summary-rev 20110114.xls

593 Contributing:

594 10C22QUANTTAB.xls

595		10C23QUANTTAB.xls
596		10C25QUANTTAB.xls
597		10D01QUANTTAB.xls
598		10E06QUANTTAB.xls
599		10E17QUANTTAB.xls
600		10E19QUANTTAB.xls
601		10F24QUANTTAB.xls
602		10F29QUANTTAB.xls
603		10F30QUANTTAB.xls
604	Data for extraction efficiency:	
605		10I16 liver extract filt + unfilt TOTALS worksheet vs avg totals.xls
606	Table 5a	Feed speciation LOQ.xls
607		Contributing:
608		10E03QUANTTAB.xls
609	Table 5b	10I23QUANTTAB.xls
610	Generally contributing to all extraction data:	extractionform.xls

611

REFERENCES

612

613 (1) Francesconi, K.A. & Kuehnelt, D. (2004) *Analyst*. **129**, 373-395

614 (2) Falnoga, I., Stibilj, E., Tusek-Znidaric, M., Slejkovec, Z., Mazej, D., Jacimovic, R., &

615 Scancar, J. (2000) *Biol.Trace Elem.Res.* **78**, 241-254

616 (3) Polatajko, A. & Szpunar, J. (2004) *JAOAC Int* **87**, 233-237

617 (4) Grant, Tyre. Dissertation (5-26-2004) *Assessing The Environmental And Biological*

618 *Implications Of Various Elements Through Elemental Speciation Using Inductively*

619 *Coupled Plasma Mass Spectrometry: Chapter 4: Characterization of arsenic species in*

620 *poultry tissue: Identification of 3-nitro-4-hydroxyphenylarsonic acid.* University of

621 Cincinnati

622 (5) Sanchez-Rodas, D., Gomez-Ariza, J.L., & Oliveira, V. (2006) *Analytical and Bioanalytical*

623 *Chemistry* **385**, 1172-1177

624 (6) Pizarro, I., Gomez, M., Camara, C., & Palacios, M.A. (2003) *Analytica Chimica Acta* **495**,

625 85-98

626 (7) Sanz, E., Munoz-Olivas, R., & Camara, C. (2005) *Journal of Chromatography A* **1097**, 1-8

627 (8) Jackson, B.P. & Bertsch, P.M. (2001) *Environmental Science & Technology* **35**, 4868-4873

628 (9) **SOP Number T012(001): Determination of Trace Elements in Aqueous Solutions by**

629 **Simultaneous Inductively Coupled Plasma Atomic Emission Spectrometry, D**

630 Heitkemper, B Barnes, J Urban, and B Zimmer, last updated Feb. 2000,

631 <http://inside.fda.gov:9003/ORA/CentralRegion/ForensicChemistryCenter/ucm040360.htm>

632 accessed on

633 (10) Guidance for Industry: **Bioanalytical Method Validation**, FDA CDER and FDA CVM,

634 last updated May 2001,

635 <http://www.fda.gov/downloads/Drugs/GuidanceComplianceRegulatoryInformation/Guidan>

636 [ces/UCM070107.pdf](http://www.fda.gov/downloads/Drugs/GuidanceComplianceRegulatoryInformation/Guidances/UCM070107.pdf) accessed on 12/6/2010

637

638

Microwave digestion

Accurately weigh ~0.5-1g of sample into clean, dry tared microwave digestion vessel (record actual mass of sample). Include at least two method blanks (0.5g DIW) for each batch, and at least one SRM in each batch. Fortified blank and SRM samples should also be included (fortify by adding a small amount of 10ppm As standard to the vessel along with the sample).

Add 10g (7mL) concentrated (70%) nitric acid (Tracemetal grade or better) to each vessel.

Seal each vessel using the torque device and place it in the numbered carousel slot.

Place carousel containing vessels in microwave. Heat according to the following program: 20 min ramp to 200°C, 20 min hold at 200°C.

Allow vessels to cool to ~50°C before opening. **Use caution when opening vessels, as contents may be under pressure, and hazardous spray may occur when pressure is released.**

Transfer contents of each vessel to a **tared** 50mL polypropylene tube. Rinse each vessel with several portions (~10mL) of deionized (DI) water, and add each rinse to the 50mL tube. Add DI water up to 50g total weight (record final mass of diluted digest solution).

Prepare calibration standards, blanks and check solutions. All standards should be matrix-matched to the digested samples. A final dilution weight of 50g for the 7mL (10g) of 100% nitric acid in the microwave vessel results in 20%(wt/wt) nitric acid. Assuming that about half of the acid is used up in the digestion process, standards should be made up in 10% (wt/wt) HNO₃. The blank should be 10% HNO₃. A reasonable set of calibration standards might be 0.05, 0.2, 1, 4, 20, and 50 ng/g As, with a 2ng/g check solution.

ICP-MS analysis


(Note that this SOP is based on the Agilent 7500ce ICP-MS system. Other systems may be used for the analysis, with details of this protocol adapted for that instrument.)

Open main valve on Argon supply.


Turn on chiller.

Tighten peristaltic pump tubing for drain, internal standard and sample introduction lines, and tighten clamps on ICP peristaltic pump.

Open “Configuration” software, and setup system for the appropriate autosampler configuration.



Click "save" and exit. Open ICP-MS software, open the Instrument Control panel and turn Plasma ON.




Wait for the plasma to ignite. If it does not come on, make a note of the error message at the bottom of instrument control window and troubleshoot.

Allow 30mins for the instrument to warm up and the plasma to stabilize. This is a good time to edit the sequence table with your sample list and place samples in the racks.

The sample introduction line should be connected to a mixing tee, along with an internal standard line, with the output of the tee going to the nebulizer. The internal standard should be ~100ng/g germanium in 7%(vol/vol) HNO₃, with a flow rate about 1/20 of the flow rate from the autosampler. This can be accomplished by using peristaltic pump tubing with three black+white tabs for the sample flow, and tubing with orange/red/orange tabs for the internal standard line.

In the tune window, load the most recent no-gas mode tune file. Check tuning by uptaking a 1ppb tune solution (1ng/g each Li, Y, Tl, Ce and Co in 2% HNO₃ + 0.5% HCl). Under the "Tune" dropdown, select "Autotune" and tune the torch horizontal and vertical position (select Hot Plasma only). Ensure that signals meet the requirements specified by the FCC ICP tuning SOP (attached- Hydrogen mode steps may be disregarded, as it is not used), adjusting tune parameters as needed. Generate a tune report, and save the tune file (for example, filename 091204ng.u for the no-gas mode tune file saved on Dec 4, 2009). Load the most recent helium-mode tune file. Check that tuning SOP requirements are met (print a tune window with about 50 data points plotted). Save the tune file. If tuning requirements are not met after adjusting tune parameters, turn off plasma and clean the cones (follow steps in Agilent maintenance video on PC). Make sure to run the analysis with the He mode tune file open.



Close the tune window, and load the method Arsenic.m (this method is set up to acquire data for m/z 72, 75 and 77, 0.1 s dwell time and 3 repetitions per mass).

Edit the Sample Log Table. Under “sequence”, click “edit sample log table”, and fill in the table (example below). Make sure the proper method is selected, make sure calibration standards are entered in the right place (calibration stds in the CALIB sequence, samples in the SMPL sequence), and check that the Vial number matches the sample’s placement in the rack. After the calibration standards, and after every ten samples thereafter, run a blank, a check solution, and another blank.

The screenshot shows the 'Edit Sample Log Table' window in the ICP-MS software. The window title is 'Edit Sample Log Table - 091202.S'. The main table has the following columns: Method, Type, Vial, Data File, Sample, Comment, Dil Lvl, ISTD Conc, and Action on Failure. The table contains 7 rows of calibration standards and 7 rows of blank samples. The 'Sample' column contains values like 'Blank', '434-46F 0.025ppb As', '434-46E 0.1ppb As', '434-46D 0.5ppb As', '434-46C 2.5ppb As', '434-46B 10ppb As', and '434-46A 50ppb As'. The 'Action on Failure' column contains 'NextSmpl'.

	Method	Type	Vial	Data File	Sample	Comment	Dil Lvl	ISTD Conc	Action on Failure	Sk
1	JCPCHEM\1\METHODS\Arsenic.M	CalStd	1101 01		Blank	10% HNO3	Level 1	Level 1	NextSmpl	
2	JCPCHEM\1\METHODS\Arsenic.M	CalStd	4101 02		434-46F 0.025ppb As	in 10% HNO3	Level 2	Level 1	NextSmpl	
3	JCPCHEM\1\METHODS\Arsenic.M	CalStd	4102 03		434-46E 0.1ppb As	in 10% HNO3	Level 3	Level 1	NextSmpl	
4	JCPCHEM\1\METHODS\Arsenic.M	CalStd	4103 04		434-46D 0.5ppb As	in 10% HNO3	Level 4	Level 1	NextSmpl	
5	JCPCHEM\1\METHODS\Arsenic.M	CalStd	4104 05		434-46C 2.5ppb As	in 10% HNO3	Level 5	Level 1	NextSmpl	
6	JCPCHEM\1\METHODS\Arsenic.M	CalStd	4105 06		434-46B 10ppb As	in 10% HNO3	Level 6	Level 1	NextSmpl	
7	JCPCHEM\1\METHODS\Arsenic.M	CalStd	4106 07		434-46A 50ppb As	in 10% HNO3	Level 7	Level 1	NextSmpl	
8										
9										
10										
11										
12										
13										
14										
15										
16										
17										
18										
19										
20										
21										
22										
23										
24										
25										
26										
27										
28										
29										
30										
31										
32										

At the bottom of the window, there are buttons for 'List Method-Cal', 'Print', 'OK', 'Cancel', and 'Help'. The 'Sample' field at the bottom left contains the text 'Blank'.

Switch from CALIB to SMPL to Sequence Flow here

The screenshot shows the 'Edit Sample Log Table - 091202.S' dialog box. The 'Method' is set to 'Arsenic.M' and the 'Sequence' is '091202.S'. The 'SMPL' dropdown menu is open, showing a list of methods including 'TERM', 'CVY', 'ICS', 'LRS', 'USER', 'CALIB', and 'BLANK'. The main table lists 33 sample entries with columns for 'method', 'Type', 'Vial', 'Data File', 'Sample', 'Comment', 'Dil/Lvl', and 'ISTD Conc'. Buttons for 'List Method-Cal', 'Print', 'OK', 'Cancel', and 'Help' are visible at the bottom.

Set up the Sequence Flow to run the CALIB stds first, followed by the SMPL sequence.

The screenshot shows the 'Edit Sample Log Table - 091202.S' dialog box with the 'Sequence Flow' section expanded. The 'Sequence Flow' dropdown is set to 'Sequence Flow'. Below it, a table lists the sequence blocks:

Block	Method
1	CALIB
2	SMPL
3	
4	
5	
6	
7	
8	
9	
10	


Buttons for 'List Method-Cal', 'Print', 'OK', 'Cancel', and 'Help' are visible at the bottom. The 'Block' dropdown at the bottom left is set to 'CALIB'.

Once all the samples have been entered, click the “sequence” dropdown and select “run”. While the sequence is running, samples may be added to or removed from the table.

Quantitation

After the sequence has finished, open the offline data analysis program. First, load the Arsenic.m method. Open one of the files from the data set to be processed. Click on the calibration button. Enter the exact concentration of each calibration standard as indicated below.

Calibration



Edit Calibration(Offline) - C:\ICPCHEM\1\CALIB\091204.c
 Calibration Type: External Calibration Method
 Calibration Title: Sample calibration curve
 Weighting Method: 1/(SD*SD)
 Mass Interpolation Fit for VIS: Point to Point

Blank Conc. Subtraction

Table | Standard Files | Graph Summary | Graph Detail

<Analyte>

	Step	Mass	Element	Curve Fit	ISTD	Weight	Min Conc	Units	01.d	02.d	03.d	04.d	06.d
									Level 1	Level 2	Level 3	Level 4	Level 5
1	1	75	As	Y=aX+blank	72	OFF	0.000		0.000	4.300E-02	2.200E-01	1.030	19.15

$Y=aX+blank$

Insert CAL STD concentrations here

<ISTD>

Step	Mass	Element	VIS	Units	All Levels
1	72	Ge			1.000

Use this slider to input level 6, 7, etc

Add Element
 Step: 1, Mass: 2, Element: , Analyte: Analyte ISTD

Configure Analyte/ISTD... Add

New... Print... OK Cancel Help

Click the “Standard Files” tab, and select the appropriate data files for each calibration level.

Edit Calibration(Offline) - C:\ICPCHEM\1\CALIB\091204.c
 Calibration Type: External Calibration Method
 Calibration Title: Sample calibration curve
 Weighting Method: 1/(SD*SD)
 Mass Interpolation Fit for VIS: Point to Point

Blank Conc. Subtraction

Table | Standard Files | Graph Summary | Graph Detail

Directory: C:\...CHEM\1\DATA\CVMSTUDY41505\09L04\00.B\

Standard Files Reload all masses to calib table


Data Correction: Bkg File: ..., Rejected Masses: ..., Interference Correction: OFF

	Data File	Sample Name	Date Acquired
1	c:\icpchem\1\data\cvmstudy41505\09l04\00.b\01.d	Blank	Dec 4 2009 10:54 am
2	c:\icpchem\1\data\cvmstudy41505\09l04\00.b\02.d	434-50F 0.043ppb As	Dec 4 2009 10:56 am
3	c:\icpchem\1\data\cvmstudy41505\09l04\00.b\03.d	434-50E 0.22pb As	Dec 4 2009 10:59 am
4	c:\icpchem\1\data\cvmstudy41505\09l04\00.b\04.d	434-50D 1.03ppb As	Dec 4 2009 11:01 am
5	c:\icpchem\1\data\cvmstudy41505\09l04\00.b\06.d	434-50B 19.15ppb As	Dec 4 2009 11:06 am
6	c:\icpchem\1\data\cvmstudy41505\09l04\00.b\07.d	434-50A 100.16ppb As	Dec 4 2009 11:08 am
7
8
9
10
11
12
13
14
15
16
17
18
19
20


Select-> <-Clear

New... Print... OK Cancel Help

Click the “Graph Detail” tab, select As, and check the calibration plot for linearity. Minimum correlation coefficient for standard curve is 0.995.




Under the “Calibrate” dropdown, save the calibration file. Under the “Tools” dropdown, select “do list”, select option for “fullQ-summary- screen”. Click OK. In the left-most section of the window, open the folder containing the data to be quantitated. In the middle section, select the data files to be processed, and “Add” them over to the right-most window. Click “Process”.



Quantitation results will open for each data file. Once the processing has finished, these can all be closed.

Open FileView32 software. Select the files to be processed (procedure similar to the previous “do list” file selection). Click “Process”.



The resulting data table will contain counts per second data for each mass in each sample. Click the “CMC” button to convert to concentration.

The screenshot shows a software window titled "Quantitation Results" with a menu bar (File, Count Info, Quant Info, Tools, Configure, Help) and a toolbar. Below the toolbar is a table with the following columns: File, Sample, Ge / 72, As / 75, and / 77. The table contains 59 rows of data, including blank samples, various methods (DIW, NIST, Fort), and different sample types (La, Lc).

File	Sample	Ge / 72	As / 75	/ 77
08.D	Blank	1	.004787	---
09.D	434-46G	1	4.007	---
10.D	Blank	1	.001516	---
11.D	DIW Method	1	.0235	---
12.D	NIST 1722c	1	.09651	---
13.D	70-5d-La	1	.08695	---
14.D	71-5d-La	1	.04479	---
15.D	73-5d-La	1	.06885	---
16.D	81-5d-La	1	.04537	---
17.D	87-5d-La	1	.04818	---
18.D	89-5d-La	1	.04858	---
19.D	92-5d-La	1	.07411	---
20.D	94-5d-La	1	.034	---
21.D	Blank	1	.001516	---
22.D	434-46G	1	3.882	---
23.D	Blank	1	.004042	---
24.D	98-5d-La	1	.04788	---
25.D	Fort method	1	8.72	---
26.D	03-5d-La	1	11.41	---
27.D	04-5d-La	1	23.58	---
28.D	06-5d-La	1	10.03	---
29.D	07-5d-La	1	20.26	---
30.D	08-5d-La	1	34.14	---
31.D	14-5d-La	1	7.304	---
32.D	15-5d-La	1	22.26	---
33.D	16-5d-La	1	11.85	---
34.D	Blank	1	.004042	---
35.D	434-46G	1	3.908	---
36.D	Blank	1	.002906	---
37.D	70-5d-Lc	1	10.38	---
38.D	19-5d-La	1	5.618	---
39.D	24-5d-La	1	43.51	---
40.D	27-5d-La	1	65.63	---
41.D	29-5d-La	1	20.04	---
42.D	32-5d-La	1	24.52	---
43.D	37-5d-La	1	71.63	---
44.D	39-5d-La	1	4.099	---
45.D	44-5d-La	1	8.277	---
46.D	49-5d-La	1	8.663	---
47.D	Blank	1	.2794	---
48.D	434-46G	1	3.954	---
49.D	Blank	1	.1981	---
50.D	51-5d-La	1	17.92	---
51.D	55-5d-La	1	2.039	---
52.D	Fort Method	1	28.86	---
53.D	Fort NIST	1	38.93	---
54.D	63-5d-La	1	9.188	---
55.D	64-5d-La	1	10.52	---
56.D	Fort Method	1	21.42	---
57.D	NIST 1577c	1	19.82	---
58.D	Blank	1	.0643	---
59.D	434-46G	1	3.958	---

Under “Tools” dropdown, this data can be exported to a CSV file, which can then be opened in Excel. Finally, calculate dilution-corrected concentrations, and check recovery of QC samples: Check Solution Recovery must be $100 \pm 10\%$. Control limits for the reference material (True Value Recovery) is $100 \pm 20\%$ from certified value. Control limits for the fortified analytical portion (FAP recovery) is $100 \pm 30\%$.





STANDARD OPERATING PROCEDURE

U.S. FOOD AND DRUG ADMINISTRATION
FORENSIC CHEMISTRY CENTER

Agilent 7500ce #2 and #3 Inductively Coupled Plasma Mass Spectrometers

Author(s): Barbara S. Barnes
John R. Urban
Nohora V. Shockey

Version Date: 06/28/2005

Approvals:		Date
Inorganic Branch Director	 Karen A. Wolnik	6/29/05
Organic Branch Director	 R. Duane Satzger	6/29/2005
Center Director	 Fred L. Fricke	6/30/05
QA Director	 Robert D. Riley	7/5/2005 Issued

This document, when viewed electronically from the network file is the current issue. If you are using a printed copy, please confirm that it is the latest issue before use.

PURPOSE:

To ensure that the Agilent 7500ce Inductively Coupled Plasma Mass Spectrometers (ICP-MS) #2 and #3 are operating within performance standards.

SCOPE:

This procedure applies to two (2) Agilent ICP-MS instruments described under equipment.

RESPONSIBILITY:

The instrument monitors for Agilent 7500ce #2 and #3 are John Urban and Nohora Shockey respectively. The alternates for Agilent 7500ce #2 are Barbara Barnes and Nohora Shockey. The alternates for Agilent 7500ce #3 are John Urban and Barbara Barnes. These individuals are responsible for seeing that the procedures described herein are performed. It is the responsibility of each analyst to verify that the instrument meets the performance standards prior to using it for sample analysis.

DEFINITIONS AND ACRONYMS:

ICP-MS – Inductively Coupled Plasma Mass Spectrometer

EM – Electron Multiplier. When the electron multiplier detector is “tuned” the applied voltages are optimized for maximum sensitivity without reducing the lifetime of the detector.

P/A – Pulse / Analog. The instrument operating software automatically switches the EM between the pulse and analog modes. P/A Factors are determined to assure calibration curve linearity across wide ranges of concentration, when both pulse and analog modes are used for detection.

SAFETY CONSIDERATIONS:

The instrument operates under high voltage and is a source of intense UV light. The instrument has integrated safety interlocks to protect the user from potential exposure which should not be defeated by the user. The instrument uses hydrogen gas (H₂). Caution must be used when hydrogen is in use. Refer to the instrument operating and hardware manuals for proper safety procedures.

The FCC Chemical Hygiene Plan and Material Safety Data Sheets should be consulted for pertinent information on the safe handling of reagents and standards used in conjunction with this equipment. The FCC Hazardous Waste Management Plan should be consulted for proper handling of wastes generated in conjunction with this equipment.

EQUIPMENT:

Agilent 7500ce ICP-MS	#2	#3
Serial Number	JP14101097	JP14101272
FDA Property Number	5105669	5105679

Logbooks: Daily Tuning Logbook, Instrument Maintenance/Repair Logbook

The equipment, logbooks, and manuals are located near instruments #2 and #3 in Rooms 117 and 121 respectively.

This document, when viewed electronically from the network file is the current issue. If you are using a printed copy, please confirm that it is the latest issue before use.

REAGENTS:

Reagent water and solutions of trace metals containing acids. Use reagents and acids which are trace metal grade or better.

Instrument Tuning Solution (1 ng/mL each of Li, Co, Y, Ce, Tl in 1% HNO₃ + 0.5% HCl): The instrument tuning solution is used when generating instrument tune reports and measuring stability. The tuning solution is prepared from stock standard solutions.

P/A Standard Solution: The P/A solutions are used to obtain the P/A factors. The P/A solutions are standard solutions that include many of the analytes in the sample to be analyzed, and cover the entire mass range. An optimal P/A tune will generate numerical factors for all masses in the acquisition method; however, it may be sufficient if factors are obtained for fewer masses, if factors are generated at several low, mid, and high masses. A standard containing 100 ng/mL is usually adequate for this purpose; however, additional solutions containing analyte at higher or lower concentration may need to be prepared. Refer to the ChemStation Operator's Manual for additional information.

QC ELEMENTS:

Analytical Limits

Tuning: Examine tuning report(s). The results must meet the limits contained in Attachment A. If not within limits, tune the ion optics and quadrupole analyzer as needed. Repeat the tune procedure.

Stability: The RSD obtained for each of the isotopes analyzed (⁷Li, ⁸⁹Y, and ²⁰⁵Tl) from the stability test must be less than 3%.

If any of the analytical limits cannot be met, consult the instrument operating manuals and/or the instrument monitor or alternate(s) for corrective action.

PROCEDURES (PROCESS DETAILS):

Check before each use:

Check the argon supply and replace if necessary. Check the liquid level in the waste container(s); dispose of waste properly as necessary. Check the condition of all glass components, sample tubing, and sampler and skimmer cones; clean or replace as necessary.

Daily with each use:

1. Start the instrument and allow to warm up approximately ½ hour. Load the method to be used for analysis. Load the appropriate tune file depending if the normal mode or the reaction cell modes (H₂ or He) will be used.

This document, when viewed electronically from the network file is the current issue. If you are using a printed copy, please confirm that it is the latest issue before use.

2. With the normal tune file loaded, tune EM prior to the first time use on each calendar week. This function is contained in the Autotune menu, which is part of the Tune screen. Save to the appropriate tune file.
3. While aspirating the 1ppb Tuning Solution, generate a tuning report for the normal mode. Compare the results of the test with the limits in Attachment A. If necessary, adjust the ion optics and quadrupole analyzer and generate a second tune report. Record the results on the Tuning Log Sheet.
4. In the normal (non-reaction cell) mode and using the method STABILIT, analyze the 1 ppb tuning solution. When the acquisition is complete, click on the button "Tabulate/Mass" and print the results. The RSD of the 10 replicates should be less than 3% for ^7Li , ^{89}Y , and ^{205}Tl . If the RSD is greater than 3% for any of the masses, repeat the acquisition one time. If the RSD's are still greater than 3%, consult the instrument monitor before proceeding with sample analyses. Record the results on the Tuning Log Sheet.
5. Load the tune file for the hydrogen mode (usually file "h2.U"). While aspirating the tuning solution, obtain values for average mean counts ($n = 200$) at m/z 56, 78 and 89. (Note that intensity at masses 56 and 78 are not analyte intensities, but from interfering polyatomic species.) Print the results from the screen and compare with the limits. Record the results on the Tuning Log Sheet. If any of the isotopes monitored do not meet the tuning specifications, adjust the tuning parameters or consult the instrument monitor or alternate for guidance.
6. Repeat, using the helium mode tune file (usually file "he.U") and isotopes 51, 75 and 89. (Note that intensity at masses 51 and 75 are not analyte intensities, but from interfering polyatomic species.) Record the results on the Tuning Log Sheet. If any of the isotopes monitored do not meet the tuning specifications, adjust the tuning parameters or consult the instrument monitor or alternate for guidance.
7. Aspirate the P/A Standard solution and generate a P/A report prior to the first time use on each calendar week. If enough factors are not generated using the P/A solution, repeat the procedure using a solution having the proper concentration of analyte. Be certain that the "Merge in the current data" box is checked when analyzing additional P/A solutions. Print the P/A factor report. Save the tune file. Copy the detector parameters to the hydrogen and helium tune files.
8. When these steps are completed, insert the reports in the appropriate Agilent 7500ce Daily Tuning Log binder.

Instrument Maintenance:

Periodic maintenance must be performed. Refer to Attachments C and D for the maintenance schedule and the maintenance checklist. The Instrument Maintenance

This document, when viewed electronically from the network file is the current issue. If you are using a printed copy, please confirm that it is the latest issue before use.

Checklist (Attachment D) must be filled out when doing the “Every Six Months” maintenance (See Attachment C) or whenever significant maintenance (such as changing pump oil or cleaning lenses) or repairs are performed.

RELATED PROCEDURES:

NA

RECORDS GENERATED:

Tuning reports, stability reports, P/A factor tuning report, and instrument maintenance checklist.

REFERENCES:

Agilent 7500 Series ICP-MS Installation Guide, Hardware Guide, and ChemStation Operator’s Manual, located near the instrument.

APPENDIXES (ATTACHMENTS):

Attachment A Daily Tuning Log

Attachment B Sample Tune Report

Attachment C Maintenance Procedures

Attachment D Instrument Maintenance Checklist

CHANGE HISTORY:

06/28/2005

1. Extend the scope of this SOP to Agilent 7500ce #3 ICP-MS
2. Include names of instrument monitor and alternates for Agilent 7500ce #3 ICP-MS.
3. Enter Agilent 7500ce #3 ICP-MS description under equipment.
4. Customize the order in which daily procedures are performed.
5. Change frequency for EM tune.
6. Remove references to the Daily Use Logbook in the “Procedures” section since this book is no longer in use. The Daily Tuning Log serves as the daily use log.
7. Generalize references to the names of the tune files to allow for use of alternate file names.
8. Under “Procedures” and in Attachment A, change the isotopes description to m/z only for the monitored interfering species in the hydrogen and helium modes.
9. Modify limits in Normal Mode (^7Li counts) and in Helium mode (counts at m/z 51 and ^{89}Y).
10. Modify limits for oxide and doubly charged in Attachment A.
11. Modify Attachment A to include results of the Stability test.
12. Under Equipment delete information on Recirculator and Autosampler since they are not unique parts of the instrument.

09/21/2005

1. Specify under Procedures the frequency in which P/A factors have to be performed.
2. Add P/A factor Tuning report to the Records Generated section.

This document, when viewed electronically from the network file is the current issue. If you are using a printed copy, please confirm that it is the latest issue before use.

Attachment A

Daily Tuning Log

Operator (Initials): _____

Date: _____

EM Tune performed? _____

P/A Factors adjusted? _____

Normal Mode

Sensitivity (1 ng/mL)	Specification		Achieved	
	Mean Counts *	% RSD	Mean Counts *	% RSD
⁷ Li	≥ 2500	< 15%		
⁸⁹ Y	≥ 8000			
²⁰⁵ Tl	≥ 4000			
Oxide (156/140)	≤ 5%			
Doubly Charged (70/140)	≤ 5%			
Axis ‡	± 0.1 of actual mass		⁷ Li	
			⁸⁹ Y	
			²⁰⁵ Tl	
W-10% ‡	0.65 to 0.80		⁷ Li	
			⁸⁹ Y	
			²⁰⁵ Tl	
Stability (n=10)	< 3%		⁷ Li	
			⁸⁹ Y	
			²⁰⁵ Tl	

Hydrogen Mode

m/z	Specification, Mean Counts *	Achieved, Mean Counts *
56	≤ 1000	
78	≤ 6	
⁸⁹ Y	≥ 3000	

Helium Mode

m/z	Specification, Mean Counts *	Achieved, Mean Counts *
51	≤ 50	
75	≤ 5	
⁸⁹ Y	≥ 1000	

* Mean Counts are the average of 200 - 0.1 second integrations of a 1 ng/mL (ppb) solution

‡ Indicate specification is achieved by checking appropriate box


Comments/Observations:

This document, when viewed electronically from the network file is the current issue. If you are using a printed copy, please confirm that it is the latest issue before use.

Attachment B SAMPLE TUNE REPORT


Tune Report

Tune File : Norm.U
 Comment : 12/03/04 no changes



Integration Time: 0.1000 sec
 Sampling Period: 0.6200 sec
 n: 200
 Oxide: 156/140 1.43%
 Doubly Charged: 70/140 2.10%

m/z	Range	Count	Mean	RSD%	Background
7	5,000	4514.0	4307.7 ✓	3.78	3.10
89	20,000	18912.0	18243.3 ✓	3.17	2.90
205	10,000	8014.0	7869.5 ✓	3.88	4.40
70/140	5	2.4%	2.1% ✓	7.76	
156/140	2	1.6%	1.4% ✓	11.73	



m/z:	7	89	205
Height:	4,299	18,419	7,942
Axis:	7.00 ✓	89.05 ✓	205.05 ✓
W-50%:	0.65 ✓	0.60	0.50
W-10%:	0.6500 ✓	0.700 ✓	0.6500 ✓

Integration Time: 0.1000 sec
 Acquisition Time: 22.7600 sec

Y axis : Linear

Attachment C

Maintenance Procedures

Routine maintenance procedures are summarized in Chapter 4 of the Agilent 7500 ICP-MS Hardware Manual, and should be consulted when there are questions about instrument maintenance.

Note that maintenance clocks have been preset in the ChemStation instrument operating software which will remind the user when major routine maintenance is due to be performed (such as checking or changing pump oil, checking lenses, etc.). When a particular procedure has been performed, the maintenance clock may be reset by the monitor or the alternate.

DAILY (OR BEFORE EACH USE): check the argon, hydrogen, and helium gas supplies and filters; check the waste vessels and empty when needed; check the condition of the peristaltic pump tubing; check the condition of the sampling and skimmer cones and their orifices. Clean or replace when necessary.

WEEKLY: Check the condition of the torch, spray chamber, end cap, and nebulizer; clean or replace if necessary. Check the water volume in the recirculator; replenish when needed. Check the oil level in the rotary pump; add oil if the level falls below the minimum.

MONTHLY: Check the level and condition of the oil in the rotary pump. Check the oil mist filter. Change oil and filter if necessary.

EVERY SIX MONTHS: Change the rotary pump oil. Check the condition of the plasma gas, auxiliary gas, and carrier gas tubing. Replace when necessary. Clean the pump strainer and vacuum the condenser fins on the recirculating chiller. Replace the water in the recirculator if algal growth is noticed in the reservoir. The Instrument Maintenance Checklist must be filled out (Attachment D).

YEARLY: Check or replace the oil mist filter of the rotary pump.

WHEN NEEDED: Evaluate and replace the electron multiplier. Clean or replace the ion lenses and the reaction lens assembly. Clean or replace the penning gauge.

Refer to the Hardware Manual for additional guidance.

This document, when viewed electronically from the network file is the current issue. If you are using a printed copy, please confirm that it is the latest issue before use.

Attachment D

INSTRUMENT MAINTENANCE CHECKLIST

AGILENT 7500ce # _____

Date _____

Operator's initials _____

Vacuum System

	Check	Clean	Adjust	Replace
Rotary Pump Oil				
Rotary Pump Mist Filters				
O Rings				
Penning Gauge				
Extraction-Omega Lens Assembly				
Octopole Assembly				

RF Generator and Sample Introduction

	Check	Clean	Adjust	Replace
Sampler and Skimmer Cones				
Gas Lines				
RF Contact Strip				
Torch Assembly				
Spray Chamber Cap & O-Ring				
Spray Chamber Drain & O-Ring				
Nebulizer Connector O-Rings				

Miscellaneous

	Check	Clean	Adjust	Replace
Water Lines				
Water Filter (at instrument inlet)				
Gas Filters				
Recirculator Fluid Level				
Recirculator Air Filter				
Recirculator Pump Strainer (Recirc. Inlet)				

Non-Routine Maintenance/Repairs:

This document, when viewed electronically from the network file is the current issue. If you are using a printed copy, please confirm that it is the latest issue before use.



A. Title: Speciation of arsenic compounds in liver resulting from roxarsone use in chickens

B. Validation History:

1. Original validation data collected: March through September 2010
2. Independent analyst validation completed: 15 Dec 2010
3. Method Trial completed:

C. Scope: This method is intended to speciate arsenic residues in chicken liver resulting from roxarsone use. In particular, it identifies and estimates roxarsone and its known metabolites, and quantifies the amount of inorganic arsenic (measured as arsenate) present above 0.6 ppb. Inorganic arsenic includes arsenite (AsIII) and arsenate (AsV); organic species such as monomethylarsinic acid (MMA) and dimethylarsinic acid (DMA) are also monitored as they may also present toxicity issues. Structures of these compounds are shown in Figure 1.




Figure 1. Structures of As compounds included in this speciation.



D. Principles: Tissue is extracted with an aqueous solution of tetramethylammonium hydroxide (TMAH), diluted with water, filtered to remove proteins and other macromolecules, and analyzed by ion chromatography-inductively coupled plasma-mass spectrometry (IC-ICP-MS). The MS is set to detect As ions (m/z 75). Peak identification is by retention time matching with external standards, with standard addition used when necessary. Quantification is by comparison to an external calibration curve.

E. Safety Precautions:

F. Reagents:

1. Water. MilliQ 18M Ω de-ionized water. Used for all water.
2. Methanol (MeOH). LC grade.
3. Nitric acid. Trace metals grade. Used to wash glassware and assist in 3-Amino dissolution.
4. Tetramethylammonium hydroxide (TMAH), 25% w/w aqueous solution. Reagent and Trace Metal grades. Prepare 0.625% solution by diluting 10 mL (reagent grade) to 400 mL with DI water. Store in polypropylene bottle at room temperature. Prepare weekly.
5. Mobile phases.
 - a. A-1% MeOH in water. Add 10.0 mL MeOH to a reservoir bottle. Bring to 1000g with water. Stir well. Use within 1 month.
 - b. B-100 mM TMAH, 1% MeOH in water. Add 36.46 g 25% Trace Metal grade TMAH, 10.0 mL MeOH, and water to 1000 g to a reservoir bottle. Stir well. Use within 1 month.
6. Standards.
 - a. *Roxarsone (4-hydroxy-3-nitrobenzenearsonic acid, CAS 121-19-7)*, >98% pure. (Rox). Acros Organics. Molecular weight = 263.03
 - b. *Arsentrioxide (CAS 85586-03-4)*. (AsIII). Spex Certiprep Speciation Standards Arsenic +3, 1000 μ g/g solution.
 - c. *Arsenic acid (CAS 1327-53-3)*. (AsV). Spex Certiprep Speciation Standards Arsenic +5, 1000 μ g/g solution.
 - d. *Dimethylarsinic acid, (CAS 75-60-5)*, 98.9% pure. (DMA). ChemService Inc. Molecular weight = 214.03



- e. *Monosodium acid methane arsenate sesquihydrate* (CAS 2163-80-6), 98.5% pure. (MMA). ChemService Inc. Molecular weight = 161.95
- f. *3-Amino-4-hydroxy-phenylarsonic acid* (CAS 2163-77-1), unknown purity. (3-Amino). Pfaltz & Bauer Rare and Fine Chemicals. Molecular weight = 233.06
- g. *N-acetyl-4-hydroxy-m-arsanilic acid ((3-acetamido-4-hydroxyphenyl)arsonic acid, CAS 97-44-9)*, unknown purity. (N-Acetyl). Pfaltz & Bauer Rare and Fine Chemicals. Molecular weight = 275.09

G. Equipment:

1. ICP-MS, Agilent 7500ce (upgraded from a 7500c), controlled by an HP Compact dc 7900 computer with Windows XP operating system and ICP-MS ChemStation, vers. B.04.00 instrumental control software. Chromatographic ICP-MS data was processed using PlasmaChrom ICP-MS Chromatographic Software, vers. C.01.00. Most chromatographic processing was done off-line on an IBM ThinkPad T series with operating system Windows XP vers. 5.1.
2. LC, Agilent 1200 Series with Instant Pilot control module (firmware B.02.07 [0001]), binary pump (firmware A.06.10 [005] with resident version A.06.10 [004]), vacuum degasser and refrigerated autosampler (firmware A.06.11 [001] with resident version A.06.10 [004]).
3. Omni-Prep Homogenizer, equipped with hard tissue disposable probes.
4. Single and multi-tube vortex mixers.
5. Centrifuges. Bench top centrifuge capable of 3000 rpm with buckets and carriers for 15mL and 50mL tubes (IEC Centra 8R with 218A rotor or equivalent).
6. Pipettors. Automatic pipettors capable of accurate delivery from 10 μ L up to 10.00 mL; repeat pipettor (Eppendorf) with assorted tips.
7. Centrifuge tubes. 15 mL polypropylene with plug-seal caps.
8. Ultrafilters. Centriprep with Ultracel-30 membrane (Millipore Corporation)
9. Autosampler vials and caps. Wash with 2% nitric acid (made from Trace Metal grade) and 4 rinses of water before using to remove trace AsIII and AsV.
10. Analytical column. IonPac AS18, 4.0 x 250 mm (Dionex cat. no. 0060549), equipped with IonPac AG18 guard column (Dionex cat. no. 060551)

H. Procedure:



1. Standard Preparation.

- a. All stock standard solutions are prepared in water except 3-Amino, which required acid to completely dissolve.
 - i. Accurately (4 significant digits) weigh 10-30 mg rox std. Dilute with 20.0 mL water, mix well and sonicate. Store refrigerated.
 - ii. Dilute 3-Amino (ca 20 mg, accurately weighed) with 20.0 mL water with 1-2 drops nitric acid added. Sonicate and mix well to dissolve. Store refrigerated. Bring to room temperature and mix well before using.
 - iii. Dilute N-Acetyl (ca 10 mg, accurately weighed) with 10.0 mL water. Sonicate and mix well to dissolve. Store refrigerated. Bring to room temperature and mix well before using.
 - iv. Dilute MMA (ca 10 mg, accurately weighed) with ca 5 g water (accurately weighed) and shake by hand to mix. (ca 1000 ppm)
 - v. Dilute MMA (ca 10 mg, accurately weighed) with ca 5 g water (accurately weighed) and shake by hand to mix. (ca 1000 ppm)
- b. For Rox, 3-Amino, N-Acetyl, MMA, and DMA, calculate stock solution concentration in mg of As per mL (or g \ddagger) as follows and verify concentration by total As determination. Use calculated concentration if experimental

$$\text{Conc} = (\text{purity} \times \text{weight in mg} \times \text{MW compound}) / (\text{volume in mL}\ddagger \times 74.92)$$

determination is within 5% of calculation. Otherwise, use experimentally determined concentration.

- c. For compounds with unknown purity (3-Amino and N-Acetyl), determine As concentration experimentally via nitric acid digestion and total As analysis. (For N-Acetyl, purity was determined to be 112%; for 3-Amino it was 78.7%.)
- d. *Working standard solutions.* Prepare in water. Calculate volumes of stock solutions of Rox, 3-Amino, and N-Acetyl which contain 2.00 mg As and dilute to 20.0 mL (g) with water. This is the 100 $\mu\text{g}/\text{mL}$ working solution. Dilute this solution 1:100 to make the 1 $\mu\text{g}/\text{mL}$ working solution. Calculate weight of MMA and DMA stock solutions which contain 0.1 mg As and dilute this amount to 100g with water to make 1 $\mu\text{g}/\text{g}$ working solutions. Prepare AsIII and AsV working solutions by dilution of the Spex Certiprep Arsenic +3 and +5, 1000 $\mu\text{g}/\text{g}$, solutions.
- e. *Calibration standards, mixed solutions at 100, 30, 10, 3, 1, 0.3, 0.1, 0.03 ng each compound as As/mL.* Not stable. Prepare daily in water. 100 ng/mL mixed



standard: Add 200 μL each 1 $\mu\text{g}/\text{mL}$ standard to a tube, bringing final volume to 2.00 mL with water. Mix 600 μL of this solution with 1400 μL water to make the 30 ng/mL solution. Prepare remaining standards by serial dilutions (200 μL + 1800 μL) from these two standards.

2. Controls and Fortified Samples (Quality Control Samples)

- a. *Method blank*: 0.5 mL water, 3 mL 0.625% TMAH, then 6.5 mL water.
- b. *Control tissue*: Either tissue from a chicken that has not been treated with Rox, or commercially purchased tissues that have been tested and shown to be free of residues > 0.6 ppb.
- c. *Fortified tissues*: 2000 ppb Rox fortified—add 10 μL of the 100 $\mu\text{g}/\text{mL}$ Rox standard to 0.5 \pm 0.05 g control tissue. 20, 4, or 2 ppb mix fortified—add 100, 20, or 10 μL 100 ng/mL mixed standard to 0.5 \pm 0.05 g control tissue. 1 ppb mix or single analyte fortified—add 50 μL 10 ng/mL mixed or single analyte standard to 0.5 \pm 0.05 g control tissue.
- d. *Fortification check samples*: Add the same volume of standard solution used to fortify tissues to a *method blank* sample. Subtract this volume from the amount of water used.

3. Sample Preparation

- a. Weigh 0.4 to 0.55 g tissue into 15 mL centrifuge tube. Record weight to nearest 0.001 g.
- b. Add 3.0 mL 0.625% TMAH.
- c. Homogenize ca 30s on Omni-Prep at 24,000 rpm.
- d. Cap tubes tightly; place on multi-tube vortexer or a shaker/rotator for at least 10 min. Ensure even mixing of all samples.
- e. Add 6.5 mL (minus standard fortification volume for QC samples) water to each sample. Cap tightly and mix well.
- f. Decant into Centriprep devices and centrifuge at 3000 rpm (2000 $\times g$) as long as needed to get >3 mL filtrate (about an hour).
- g. Transfer a portion of the filtrate solutions to autosampler vials.

4. Instrument Operating Parameters

- a. Set LC parameters as follows:



- i. Flow rate: 1 mL/min
 - ii. Injection volume: 50 μ L
 - iii. Autosampler temperature: 10°C
 - iv. Needle wash program: 3 sec using Mobile Phase B as the solvent
 - v. Runtime: 50 min
 - vi. Gradient: 45% B 0-17 min, linear ramp to 70% B 17- 17.1 min, 70% B 17.1-42 min, linear ramp to 45% B at 42-42.1 min.
- b. ICP-MS parameters:
- i. The following settings may be used as a starting point. ICP-MS tuning should be checked daily to ensure satisfactory performance. Save tune file as AsChrom.U.
 - (a) RF power: 1500 W
 - (b) Carrier gas: 1.1 L min⁻¹
 - (c) Makeup gas: 0.1 L min⁻¹
 - (d) Spray chamber temp: 2°C
 - (e) Nebulizer type: glass concentric
 - (f) Sampling depth: 8.5 mm
 - (g) Ions monitored: 75 (As), 77 (⁴⁰Ar³⁷Cl)
 - (h) Dwell time: 0.8 s (*m/z* 75), 0.2 s (*m/z* 77)
 - (i) Collision cell: ON, He mode
 - (j) Collision gas flow: 5.7 mL min⁻¹
 - ii. Set up a time-resolved method, monitoring ions 75 and 77 for 0.8 and 0.2s, respectively, and 1 replicate per ion. Acquire data for 3000 s (50 min).
5. Procedure for Instrumental Analysis of Samples, Controls and Standards.



- a. Turn on chiller and tighten peristaltic pump clamp onto spray chamber drain tubing. Open ICP ChemStation and ignite plasma. Tune according to lab or manufacturer procedure.
- b. Start LC. Connect cable to allow communication between ICP and HPLC. Ensure that backpressure is acceptable (2000 psi is normal, >3000 is problematic) and that ICP source is being drained by peristaltic pump (bubbles should be visible in drain tubing).
- c. Create/edit the sequence file on the ICP-MS data system. Make sure that the injection list on the LC controller matches the ICP-MS sequence.
- d. Inject a standard (100 ng/mL most commonly used) and check peak shape and response.
- e. Inject calibration standards set, method blank and fortification check samples, other QC samples, and sample extracts. Inject water blanks as needed. Follow with a repeat injection from one of the standards. Additional standards may be injected throughout the sequence to monitor retention time drift and instrument response.

I. Calculations:

1. Calculate the dilution factor (to be entered in ICP-MS sequence file) for each sample: $(9.5 + \text{sample wt}) / \text{sample wt}$. Use a dilution factor of 20 for fortified samples, method blanks, and fortification check samples.
2. Process data using ICP-MS Chromatographic software. Create a processing method with DMA, AsIII, MMA, AsV, 3-Amino, N-Acetyl, Rox, and unknown compounds. Use m/z 75 as the target ion, signal extraction time abt +/- 2 min, measure response by area, and identify peaks by best RT match. Construct standard curves for each known analyte using Linear Regression with $1/x^2$ weighting. Use Avg of Response Factors for the unknown compounds. Select as the response factor a number close to the slopes of responses for the nearest known analytes. Usually, this number is close to 10000.
3. Open a processing method and save it using the current date or date of analysis as part of the filename.
4. Starting with the calibration standards, open each data file, quantify it ("Quantitate," then "Calculate"), review the peaks, and save the review. For standards, select "Calibrate," "Update," "Update One Level," and replace the response and retention time with the current file's data. Save the method once all levels have been updated. Quantify, review, and save remaining data files from the sequence.



5. Export the quant results from these data files to a single csv file for further processing and summarization. The file will be named QUANTTAB.csv as the software default. Open the csv file in Excel. Add a new worksheet called Conc Summary, copy the concentrations (3rd column of numbers) from each with analysis, and paste it with transposition to a row on the new worksheet. Ensure each row has the correct sample name. Add all species columns to calculate total As species. Save the file, usually as *setname*QUANTTAB.xls, where *setname* follows the ChemStation convention of yyadd, with *a* being a single character code for month (A=January, B=February, etc.).

J. System Suitability and Quality Control:

1. Instrument Performance Specifications

- a. All peaks in the initial test standard (step H.5.d.) should be well-resolved except DMA, AsIII, and MMA, and appear similar to Fig. 1. Compounds shown are DMA (2.8 min), AsIII (3.0 min), MMA (3.2 min), AsV (7.0 min), 3-Amino (9 min), N-acetyl (10 min), and Rox (18 min).




Figure 2. Representative arsenic speciation ion chromatogram. This is a mix of 7 standards prepared in water at 10 ng/mL.

- b. The peak height of AsV should be approximately 10^5 in the 100 ng/mL standard.
- c. All peaks in the next to lowest standard (0.1 ng/mL) should be readily detectable ($S/N > 3$) by visual inspection. AsV should be readily detectable in the lowest standard (0.03 ng/mL).

2. Critical Points and Stopping Points

- a. Peak identity in sample extracts should be verified by individual standard addition if there is suspected retention time shift.

3. Stability of Analyte in Samples and Extracts



- a. Stability is to be determined in samples. Literature suggests overall stability is not likely to be a problem in tissues at -80° . Preliminary results suggest that samples should not be repeatedly thawed and refrozen.
 - b. Extracts in TMAH are not stable. We have observed some conversion of species in standard solutions and extracts after only 3 days at 10° . Extracts analysis should be completed within 48 hours of preparation.
4. Acceptance Tests for Critical Reagents
- a. The method blank should be free of peaks >0.6 ppb (0.03 ng mL^{-1}).
5. Acceptance Criteria for Results
- a. Coefficients of determination (r-square) for DMA, AsV, N-Acetyl, and Rox should exceed 0.99.

K. References:

1. Grant, T. 2004. Assessing The Environmental And Biological Implications Of Various Elements Through Elemental Speciation Using Inductively Coupled Plasma Mass Spectrometry: Chapter 4: Characterization of arsenic species in poultry tissue: Identification of 3-nitro-4-hydroxyphenylarsonic acid. Ph. Dissertation, University of Cincinnati. Inspiration for extraction.
2. Jackson, B. P. and Bertsch, P. M. (2001) "Determination of Arsenic Speciation in Poultry Wastes by IC-ICP-MS." *Environmental Science & Technology* **35**:4868-4873. Inspiration for chromatography.
3. Slingsby, Rosanne W., Al-Horr, Rida, Pohl, Christopher A, and Lee, Joung Hae (May 2007) "Use of Dual-Selectivity IC-ESI-MS for the Separation and Detection of Anionic and Cationic Arsenic Species." *American Laboratory*. Inspiration for choice of IonPac AS18 column.



L. Validation Data:

Accuracy and Precision from Fortified Controls

Fortification Level		DMA	AsIII	MMA	AsV	3-Amino	N-Acetyl	Rox	iAs (III+V)
2 ppm Rox n = 16	Recovery							73%	
	RSD							28%	
20 ppb mix n = 12	Recovery	106%	10%	82%	194%	83%	75%	72%	103%
	RSD	10%	275%	23%	15%	33%	33%	33%	12%
4 ppb mix n = 7	Recovery	116%	0%	94%	201%	103%	95%	86%	105%
	RSD	9%		17%	15%	32%	30%	28%	16%
2 ppb mix n = 10	Recovery	117%	2%	86%	185%	86%	81%	70%	103%
	RSD	7%	316%	20%	18%	32%	29%	33%	15%
1 ppb mix n = 10	Recovery	119%	0%	85%	157%	113%	78%	62%	103%
	RSD	23%		32%	23%	25%	42%	86%	19%
1 ppb AsV n = 5	Recovery				102%				93%
	RSD				13%				17%