

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01		Page 1 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Contents

A.	INTRODUCTION	2
B.	EQUIPMENT	2
C.	REAGENTS AND SOLUTIONS	3
D.	STANDARD(S)	4
E.	SAMPLE PREPARATION.....	20
F.	ANALYTICAL PROCEDURE	20
G.	CALCULATIONS / IDENTIFICATION	29
H.	SAFETY INFORMATION AND PRECAUTIONS	31
I.	QUALITY ASSURANCE PLAN	31
J.	APPENDIX.....	32
K.	APPROVALS AND AUTHORITIES.....	36

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 2 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

A. INTRODUCTION

1. Summary of Procedure

Animal drug residues are extracted from tissue using dispersive SPE for both extraction and sample clean up. The extracted residues are examined using UHPLC-MS-MS using a triple quadrupole mass spectrometer under electrospray ionization (ESI) conditions. Analytes are identified by comparison against matrix matched standards.

2. Applicability

This method is suitable for the screening and confirmation of animal drug residues in beef and porcine kidney and muscle tissue at the levels listed in Tables 18 and 19 in Appendix J.3.

Note: Refer to 21CFR for tolerance values set by FDA and 40CFR for tolerance values set by EPA.

B. EQUIPMENT

Note: Equivalent equipment may be substituted.

1. Apparatus

- a. Platform shaker - Cat. No. 6010, Eberbach
- b. Centrifuge - Thermo IEC, Sorvall RC-6 capable of 3720 rcf
- c. Balance - Mettler Top Loading Model PB300 Balance capable of weighing 2 ± 0.01 g
- d. Balance Analytical - Mettler Model X-205 Dualrange
- e. Turbovap LV Concentration Workstation - Biotage Corp
- f. C18 BakerBond Octadecyl (C18) 40 μ m Prep LC Packing
- g. Centrifuge tubes - Polypropylene (PP), 50 mL, Falcon Part number 352070
- h. Centrifuge tubes - Polypropylene (PP), 15 mL, Falcon Part Number 352096
- i. Whatman Mini-UniPrep Syringless filter vials - VWR 0.2 micron, PVDF, Cat. No. 12000-524.

Note: Avoid glass if the Mini-UniPrep filter vials are substituted with syringe filters and autosampler vials, and substitutes must be checked for possible retention of analytes.

- j. Magnetic stirrer and stirbars, freezer, volumetric flasks, graduated cylinders, Pasteur pipettes, repeating pipettes and tips, beakers, bottles, weigh boats, spatulas, funnels, bottle top volumetric dispensers, and other items.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 3 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

- k. LC vials with screw cap lids- Amber glass, 4 mL, Cat. No. 1232 P 49, Thomas Scientific
 - l. Plastic screw cap vials - Polypropylene, 4 mL, Cat. No. 1708 H 01, Thomas Scientific
2. Instrumentation
 - a. Waters UHPLC-MS/MS TQD system with MassLynx operating software.
 - b. UHPLC column - Waters UHPLC HSS T3, 2.1 x 100 mm, 1.8 μ m column with VanGuard Precolumn UHPLC HSS T3 2.1 x 5.0 mm, 1.8 μ m.

C. REAGENTS AND SOLUTIONS

Note: Equivalent reagents / solutions may be substituted. The stability time frame of the solution is dependant on the expiration date of the components used or the listed expiration date, whichever is soonest.

1. Reagents
 - a. Hexane - HPLC Grade, Fisher Optima Cat. No. H303-4
 - b. Acetonitrile (ACN) - HPLC Grade, Spectrum Chemical Co. Cat. No. HP412
 - c. Formic acid - Sigma Chemical Co., Cat. No. F0507-500ML
 - d. Water - LC Grade, House deionized water passed through an ELGA Pure Lab Ultra Filtration System.
 - e. Sodium hydroxide (NaOH) - pellets, Fisher, Item No. S318-3.
2. Solutions
 - a. 80:20 Acetonitrile/Water:

Measure 800 mL of acetonitrile using a graduated cylinder and transfer to a 1 L volumetric flask. Measure 200 mL of deionized water using a graduated cylinder and add to the volumetric flask containing the acetonitrile. Mix this solution and transfer to a dispenser bottle.
 - b. Hexane (saturated with acetonitrile):

Add 40 mL of acetonitrile to 1 L of hexane in a separatory funnel. Mix this solution vigorously. Allow the layers to separate and discard the acetonitrile layer. Transfer the saturated hexane to a dispenser bottle for storage. Mix this solution prior to use each day.
 - c. 0.1% Formic Acid in water:

Add 1.0 mL of formic acid to a 1 L volumetric flask. Dilute to volume with deionized water.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 4 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS	
Revision: .01	Replaces: CLG-MRM1.00
Effective: 8/06/2012	

- d. 0.03 M Sodium Hydroxide:
Add 0.12 g of NaOH to a 100 ml volumetric flask containing 80 mL of deionized water. Mix and allow solution to cool. Adjust to final volume using deionized water. Store in a plastic container.
- e. UHPLC Aqueous Mobile Phase (5% ACN, 95% Water, 0.1% Formic Acid):
Measure 50 mL of acetonitrile and add to a 1 L volumetric flask. Add 1.0 mL of formic acid to the flask. Bring to volume with deionized water. Mix and transfer to the aqueous reservoir of the LC.
- f. UHPLC Organic Mobile Phase (Acetonitrile, 0.1% Formic Acid):
Add 1.0 mL of formic acid into a 1 L volumetric flask. Bring to volume using acetonitrile. Mix and transfer to the organic reservoir of the LC.

D. STANDARD(S)

Note: Equivalent standards / solutions may be substituted. Purity and counterions are to be taken into account when calculating standard concentrations. The stability time frame of the solution is dependant on the expiration date of the components used or the listed expiration date, whichever ends sooner.

1. Standard Information

Table 1 – Standard information

Analyte Name	Manufacturer	Catalog number	Analyte Name	Manufacturer	Catalog number
Amoxicillin	US Pharmacopeia	1031503	Salbutamol	Sigma Aldrich	S8260
Ampicillin	Sigma Aldrich	A1593	Sarafloxacin	Abbott Labs ⁵	Not applicable
Beta/ Dexamethasone	MP Biochemicals	154853	Sulfachloro pyridazine	Sigma Aldrich	46778
Cefazolin	Sigma Aldrich	C5020	Sulfadiazine	Sigma Aldrich	S8626
Chloramphenicol	Sigma Aldrich	31667	Sulfadimethoxine	Sigma Aldrich	46794
Chlortetracycline	US Pharmacopeia	1129007	Sulfadoxine	US Pharmacopeia	1626500
Cimaterol	Tocris Bioscience	0435	Sulfaethoxy pyridazine	Fluka	02743
Ciprofloxacin	US Pharmacopeia	1134313	Sulfamerazine	Sigma Aldrich	S8876
Clindamycin	Sigma Aldrich	C5269	Sulfamethazine	Sigma Aldrich	S6256
Cloxacillin	US Pharmacopeia	1142005	Sulfamethizole	Sigma Aldrich	S5632
Danofloxacin	Pfizer ¹	Not applicable	Sulfamethoxazole	Sigma Aldrich	S7507
DCCD	Pfizer ¹	Not applicable	Sulfamethoxy pyridazine	Sigma Aldrich	S7257

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 5 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS	
Revision: .01	Replaces: CLG-MRM1.00
Effective: 8/06/2012	

Analyte Name	Manufacturer	Catalog number	Analyte Name	Manufacturer	Catalog number
Desethylene Ciprofloxacin	Bayer Healthcare ²	Not applicable	Sulfanilamide	US Pharmacopeia	1632004
Dicloxacillin	US Pharmacopeia	1189009	Sulfantran	Sigma Aldrich	46882
Difloxacin	Abbott Labs ⁵	Not applicable	Sulfapyridine	Sigma Aldrich	S6252
Enrofloxacin	Bayer Healthcare ²	Not applicable	Sulfaquinoxaline	Sigma Aldrich	45662
Erythromycin A	Sigma Aldrich	E0774	Sulfathiazole	Sigma Aldrich	S9876
Florfenicol	Sigma Aldrich	F1427	Tetracycline	US Pharmacopeia	1651009
Flunixin	US Pharmacopeia	1274607	Tilmicosin	Lilly ⁴	Not applicable
Gamithromycin	Hovione Farma ³	Not applicable	Tulathromycin A	Pfizer ¹	CP-472,295
Lincomycin	Sigma Aldrich	L6004	Tylosin	Sigma Aldrich	T6134
Melengestrol Acetate	MP Biochemicals	158952	Zeranol (B-Zearalanol)	Sigma Aldrich	Z0417
Nafcillin	US Pharmacopeia	1450007	Flunixin-d3	Sigma Aldrich ⁶	34083 (Vetranal)
Norfloxacin	Sigma Aldrich	N9890	¹³ C6 Sulfamethazine Phenyl	Sigma Aldrich ⁶	32519 (Vetranal)
Oxacillin	US Pharmacopeia	1481000	d7 Penicillin G	Toronto Research ⁶	B288600
Oxyphenylbutazone	Toronto Research	0876950			
Oxytetracycline	Sigma Aldrich	O5875			
Penicillin G	US Pharmacopeia	1502508			
Phenylbutazone	MP Biochemicals	153567			
Pirlimycin	Pfizer ¹	Not applicable			
Prednisone	Sigma Aldrich	P6254			
Ractopamine	Sigma Aldrich	34198			

¹ - Pfizer, Groton, CT

² - Bayer Healthcare, AG Business Group Pharma, PH-GDD-PT, Clinical Supplies Ops , Wuppertal, Germany.

³ - Hovione FarmaCiencia SA, Sete Casa, Loures, Portugal

⁴ - Lilly Corporate Center, Indianapolis, Indiana

⁵ - Abbott Labs, Chicago, IL

⁶ - Internal Standard (IS) – optional

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01		Page 6 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

2. Preparation of Standard Solution(s)

a. Animal Drug Stock Solutions and internal standard stock solutions

Prepare animal drug stock solutions and internal standard stock solutions at approximately 1.0 mg/mL when adequate material is available. Other concentrations are used based on two criteria:

- i. Solubility of the drug in the solvent
- ii. Cost and availability of the drug

For each stock solution, calculate the amount of base material needed (ex. accounting for purity and/or water and sulfate content) to prepare at the concentration listed below using the appropriate solvent listed.

Table 2 – Stock standard concentrations

Standard Analyte	Category	Solvent used	Stock Standard Concentration (ng/μL)
Cimaterol	Acetonitrile Mix	Acetonitrile	1000
Ractopamine	Acetonitrile Mix	Water	1000
Salbutamol	Acetonitrile Mix	Acetonitrile or methanol	1000
Chloramphenicol	Acetonitrile Mix	Acetonitrile	1000
Florfenicol	Acetonitrile Mix	Acetonitrile	1000
Ciprofloxacin	Acetonitrile Mix	0.03 M NaOH	1000
Danofloxacin	Acetonitrile Mix	0.03 M NaOH	1000
Difloxacin	Acetonitrile Mix	50%ACN/MeOH	500
Enrofloxacin	Acetonitrile Mix	Acetonitrile	500
Norfloxacin	Acetonitrile Mix	Acetonitrile	1000
Sarafloxacin	Acetonitrile Mix	Methanol	1000
Desethylene Ciprofloxacin	Acetonitrile Mix	0.03 M NaOH	300
Sulfachloropyridazine	Acetonitrile Mix	Acetonitrile	1000
Sulfadiazine	Acetonitrile Mix	Acetonitrile	1000
Sulfadimethoxine	Acetonitrile Mix	Acetonitrile	1000
Sulfadoxine	Acetonitrile Mix	Acetonitrile	1000
Sulfaethoxypyridazine	Acetonitrile Mix	Acetonitrile	1000
Sulfamerazine	Acetonitrile Mix	Acetonitrile	1000
Sulfamethazine	Acetonitrile Mix	Acetonitrile	1000
Sulfamethizole	Acetonitrile Mix	Acetonitrile	1000

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 7 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Standard Analyte	Category	Solvent used	Stock Standard Concentration (ng/μL)
Sulfamethoxazole	Acetonitrile Mix	Acetonitrile	1000
Sulfamethoxypyridazine	Acetonitrile Mix	Acetonitrile	1000
Sulfanilamide	Acetonitrile Mix	Acetonitrile	2000
Sulfantran	Acetonitrile Mix	Acetonitrile	500
Sulfapyridine	Acetonitrile Mix	Acetonitrile	1000
Sulfaquinoxaline	Acetonitrile Mix	Acetonitrile	500
Sulfathiazole	Acetonitrile Mix	Acetonitrile	1000
Clindamycin	Acetonitrile Mix	Acetonitrile	1000
Erythromycin A	Acetonitrile Mix	Acetonitrile	1000
Gamithromycin	Acetonitrile Mix	Acetonitrile	500
Lincomycin	Acetonitrile Mix	50%ACN/MeOH	500
Pirlimycin	Acetonitrile Mix	50%ACN/MeOH	1000
Tilmicosin	Acetonitrile Mix	Acetonitrile	1000
Tulathromycin A	Acetonitrile Mix	Acetonitrile	1000
Tylosin	Acetonitrile Mix	Acetonitrile	1000
Chlortetracycline	Acetonitrile Mix	Methanol	500
Oxytetracycline	Acetonitrile Mix	Methanol	1000
Tetracycline	Acetonitrile Mix	Methanol	500
Flunixin	Acetonitrile Mix	Methanol	1000
Phenylbutazone	Acetonitrile Mix	Acetonitrile	1000
Prednisone	Acetonitrile Mix	Methanol	1000
Oxyphenylbutazone	Acetonitrile Mix	Acetonitrile	1000
Beta/Dexamethasone	Acetonitrile Mix	Acetonitrile	1000
Melengestrol Acetate	Acetonitrile Mix	Acetonitrile	1000
Zeranol (B-Zearalanol)	Acetonitrile Mix	Methanol	1000
Amoxicillin	Beta Lactam Mix	Water	350
Ampicillin	Beta Lactam Mix	Water	1000
DCCD	Beta Lactam Mix	Water	1000
Cefazolin	Beta Lactam Mix	Water	1000
Cloxacillin	Beta Lactam Mix	Water	1000
Dicloxacillin	Beta Lactam Mix	Water	1000
Naficillin	Beta Lactam Mix	Water	1000
Oxacillin	Beta Lactam Mix	Water	1000

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 8 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Standard Analyte	Category	Solvent used	Stock Standard Concentration (ng/μL)
Penicillin G	Beta Lactam Mix	Water	1000
Flunixin-d3	Internal Standard Mix	Methanol	1000
13C6 Sulfamethazine Phenyl	Internal Standard Mix	Acetonitrile	1000
d7 Penicillin G	Internal Standard Mix	Water	500

“Acetonitrile Mix” stock standards will expire in 6 months when stored at ≤ -10 °C or at the time of the earliest expiring component. “Beta Lactam Mix” stock standards will expire in 2 months when stored at ≤ -10 °C or at the time of the earliest expiring component. “Internal Standard Mix” stock standards will expire in 2 months for d7 Penicillin G, or 6 months for Flunixin-d3 and 13C6 Sulfamethazine Phenyl when stored at ≤ -10 °C.

Note: Internal standards are optional in this method and can be used to monitor injection sequence performance within a set. For issues observed, such as inconsistent internal standard area counts, samples may be reinjected or reanalyzed as needed. If internal standards are not used, appropriate volumes and chemicals must be adjusted as outlined in appropriate method steps below.

b. Intermediate standard solutions

Prepare individual intermediate standard solutions as described for the analytes below in 10 mL volumetric flasks. Intermediate stock standards will expire in 6 months when stored at ≤ -10 °C or at the time of the earliest expiring component.

Table 3 – Intermediate standard solutions

Analyte	Stock Standard Concentration (ng/μL)	Solvent used	Volume Stock Standard (μL)	Intermediate Standard Concentration (ng/μL)
Cimaterol	1000	Acetonitrile	500	50
Ractopamine	1000	Acetonitrile	500	50
Salbutamol	1000	Acetonitrile	500	50
Chloramphenicol	1000	Acetonitrile	500	50
Flunixin	1000	Acetonitrile	1000	100
Zeranol (B-Zearalanol)	1000	Acetonitrile	500	50
Melengesterol Acetate	1000	Acetonitrile	500	50

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 9 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

- c. Antibiotic drug composite working (spiking) and internal standard (spiking) solutions
- i. Prepare the composite “Acetonitrile Mix” working solution(s) for the veterinary drugs contained in the acetonitrile spiking solutions using the stock and intermediate standard solutions above and the volumes listed in the tables below.
- Calculate or use the volume of stock or intermediate stock solution required to give the concentration listed for each standard type in the table below.
 - Pipet the calculated volume of stock into a 10 mL volumetric flask.
 - Dilute to 10 mL volume with acetonitrile.
 - Cap flask and mix.
 - Transfer solution into amber glass LC vials with screw cap lids.
 - Composite working solutions will expire in 6 months when stored at ≤ -10 °C or at the time of the earliest expiring component.

Table 4 – “Acetonitrile Mix” kidney working solution preparation

Standard Analyte	Stock or Intermediate Standard Concentration (ng/ μ L)	Bovine Kidney Screening Volume (μ L)	Bovine Kidney Confirmation Volume (μ L)	Porcine Kidney Screening Volume (μ L)	Porcine Kidney Confirmation Volume (μ L)
Beta/Dexamethasone	1000	25	25	25	25
Chloramphenicol	50	60		60	
Chlortetracycline	500	1000	1000	1000	1000
Cimaterol	50	120	120	30	30
Ciprofloxacin	1000	12.5	25	12.5	12.5
Clindamycin	1000	25	25	25	25
Danofloxacin	1000	12.5	12.5	12.5	12.5
Desethylene Ciprofloxacin	300	41.7	41.7	41.7	41.7
Difloxacin	500	25	25	25	25
Enrofloxacin	500	25	50	25	25
Erythromycin A	1000	25	25	25	25
Florfenicol	1000	50	50	50	50
Flunixin	100	62.5	62.5	62.5	62.5
Gamithromycin	500	50	100	50	100

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 10 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS	
Revision: .01	Replaces: CLG-MRM1.00
Effective: 8/06/2012	

Standard Analyte	Stock or Intermediate Standard Concentration (ng/μL)	Bovine Kidney Screening Volume (μL)	Bovine Kidney Confirmation Volume (μL)	Porcine Kidney Screening Volume (μL)	Porcine Kidney Confirmation Volume (μL)
Lincomycin	500	50	50	50	50
Norfloracin	1000	12.5	12.5	12.5	12.5
Oxyphenylbutazone	1000			25	
Oxytetracycline	1000	250	250	250	250
Phenylbutazone	1000			25	25
Pirlimycin	1000	125	125	125	125
Prednisone	1000	25	25	25	25
Ractopamine	50	30	120	30	
Salbutamol	50	60		30	
Sarafloxacin	1000	12.5	12.5	12.5	12.5
Sulfachloropyridazine	1000	25	25	25	100
Sulfadiazine	1000	25	25	25	25
Sulfadimethoxine	1000	25	25	25	25
Sulfadoxine	1000	25	25	25	25
Sulfaethoxypyridazine	1000	25	25	25	25
Sulfamerazine	1000	25	25	25	25
Sulfamethazine	1000	25	25	25	25
Sulfamethizole	1000	25	25	25	25
Sulfamethoxazole	1000	25	50	25	25
Sulfamethoxypyridazine	1000	25	25	25	25
Sulfanilamide	2000			25	50
Sulfanitran	500	50	100	50	50
Sulfapyridine	1000	25	25	25	25
Sulfaquinoxaline	500	50	50	50	50
Sulfathiazole	1000	25	25	25	25
Tetracycline	500	500	2000	500	500
Tilmicosin	1000	60		120	

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01		Page 11 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Standard Analyte	Stock or Intermediate Standard Concentration (ng/μL)	Bovine Kidney Screening Volume (μL)	Bovine Kidney Confirmation Volume (μL)	Porcine Kidney Screening Volume (μL)	Porcine Kidney Confirmation Volume (μL)
Tylosin	1000	50	50	100	100
Zeranol (B-Zearalanol)	50			120	120

Table 5 – “Acetonitrile Mix” muscle working solution preparation

Standard Analyte	Stock or Intermediate Standard Concentration (ng/μL)	Bovine Muscle Screening Volume (μL)	Bovine Muscle Confirmation Volume (μL)	Porcine Muscle Screening Volume (μL)	Porcine Muscle Confirmation Volume (μL)
Beta/Dexamethasone	1000	25	25	25	25
Chloramphenicol	50	60		60	
Chlortetracycline	500	1000	1000	1000	1000
Cimaterol	50	30		30	
Ciprofloxacin	1000	12.5	25	12.5	12.5
Clindamycin	1000	25	25	25	25
Danofloxacin	1000	12.5	12.5	12.5	12.5
Desethylene Ciprofloxacin	300	41.7	41.7	41.7	41.7
Difloxacin	500	25	25	25	25
Enrofloxacin	500	25	25	25	50
Erythromycin A	1000	25	25	25	25
Florfenicol	1000	50	50	50	50
Flunixin	100	62.5	62.5	62.5	62.5
Gamithromycin	500	50	50	50	200
Lincomycin	500	50	50	50	50
Melengestrol Acetate	50	100		100	
Norfloxacin	1000	12.5	12.5	12.5	12.5
Oxyphenylbutazone	1000	25	25	25	50
Oxytetracycline	1000	250	250	250	250
Phenylbutazone	1000	25	25	100	100

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01		Page 12 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Standard Analyte	Stock or Intermediate Standard Concentration (ng/μL)	Bovine Muscle Screening Volume (μL)	Bovine Muscle Confirmation Volume (μL)	Porcine Muscle Screening Volume (μL)	Porcine Muscle Confirmation Volume (μL)
Pirlimycin	1000	125	125	125	125
Prednisone	1000	25	50	25	25
Ractopamine	50	30	60	30	60
Salbutamol	50	30		30	
Sarafloxacin	1000	12.5	12.5	12.5	12.5
Sulfachloropyridazine	1000	25	25	25	25
Sulfadiazine	1000	25	25	50	50
Sulfadimethoxine	1000	25	25	25	25
Sulfadoxine	1000	25	25	25	25
Sulfaethoxypyridazine	1000	25	25	25	25
Sulfamerazine	1000	25	25	25	25
Sulfamethazine	1000	25	100	25	25
Sulfamethizole	1000	25	25	25	25
Sulfamethoxazole	1000	25	25	25	25
Sulfamethoxypyridazine	1000	25	25	25	25
Sulfanitran	500	50	200	50	100
Sulfapyridine	1000	25	25	25	25
Sulfaquinoxaline	500	50	50	50	50
Sulfathiazole	1000	25	25	25	25
Tetracycline	500	500	2000	500	1000
Tilmicosin	1000	30	30	30	30
Tulathromycin A	1000	500	500	500	500
Tylosin	1000	50	50	50	50

Using these volumes yield the following concentrations as found in the tables below:

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01		Page 13 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Table 6 – Acetonitrile Mix kidney working standard levels

Standard Analyte	Bovine Kidney Screening Level (ng/μL)	Bovine Kidney Confirmation Level (ng/μL)	Porcine Kidney Screening Level (ng/μL)	Porcine Kidney Confirmation Level (ng/μL)
Beta/Dexamethasone	2.5	2.5	2.5	2.5
Chloramphenicol	0.3		0.3	
Chlortetracycline	50	50	50	50
Cimaterol	0.6	0.6	0.15	0.15
Ciprofloxacin	1.25	2.5	1.25	1.25
Clindamycin	2.5	2.5	2.5	2.5
Danofloxacin	1.25	1.25	1.25	1.25
Desethylene Ciprofloxacin	1.25	1.25	1.25	1.25
Difloxacin	1.25	1.25	1.25	1.25
Enrofloxacin	1.25	2.5	1.25	1.25
Erythromycin A	2.5	2.5	2.5	2.5
Florfenicol	5	5	5	5
Flunixin	0.625	0.625	0.625	0.625
Gamithromycin	2.5	5	2.5	5
Lincomycin	2.5	2.5	2.5	2.5
Norfloxacin	1.25	1.25	1.25	1.25
Oxyphenylbutazone			2.5	
Oxytetracycline	25	25	25	25
Phenylbutazone			2.5	2.5
Pirlimycin	12.5	12.5	12.5	12.5
Prednisone	2.5	2.5	2.5	2.5
Ractopamine	0.15	0.6	0.15	
Salbutamol	0.3		0.15	
Sarafloxacin	1.25	1.25	1.25	1.25
Sulfachloropyridazine	2.5	2.5	2.5	10
Sulfadiazine	2.5	2.5	2.5	2.5
Sulfadimethoxine	2.5	2.5	2.5	2.5
Sulfadoxine	2.5	2.5	2.5	2.5

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 14 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Standard Analyte	Bovine Kidney Screening Level (ng/μL)	Bovine Kidney Confirmation Level (ng/μL)	Porcine Kidney Screening Level (ng/μL)	Porcine Kidney Confirmation Level (ng/μL)
Sulfaethoxypyridazine	2.5	2.5	2.5	2.5
Sulfamerazine	2.5	2.5	2.5	2.5
Sulfamethazine	2.5	2.5	2.5	2.5
Sulfamethizole	2.5	2.5	2.5	2.5
Sulfamethoxazole	2.5	5	2.5	2.5
Sulfamethoxypyridazine	2.5	2.5	2.5	2.5
Sulfanilamide			5	10
Sulfanitran	2.5	5	2.5	2.5
Sulfapyridine	2.5	2.5	2.5	2.5
Sulfaquinoxaline	2.5	2.5	2.5	2.5
Sulfathiazole	2.5	2.5	2.5	2.5
Tetracycline	25	100	25	25
Tilmicosin	6		12	
Tylosin	5	5	10	10
Zeranol (B-Zearalanol)			0.6	0.6

Table 7 – Acetonitrile Mix muscle working standard levels

Standard Analyte	Bovine Muscle Screening Level (ng/μL)	Bovine Muscle Confirmation Level (ng/μL)	Porcine Muscle Screening Level (ng/μL)	Porcine Muscle Confirmation Level (ng/μL)
Beta/Dexamethasone	2.5	2.5	2.5	2.5
Chloramphenicol	0.3		0.3	
Chlortetracycline	50	50	50	50
Cimaterol	0.15		0.15	
Ciprofloxacin	1.25	2.5	1.25	1.25
Clindamycin	2.5	2.5	2.5	2.5
Danofloxacin	1.25	1.25	1.25	1.25
Desethylene	1.25	1.25	1.25	1.25

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 15 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS	
Revision: .01	Replaces: CLG-MRM1.00
Effective: 8/06/2012	

Standard Analyte	Bovine Muscle Screening Level (ng/μL)	Bovine Muscle Confirmation Level (ng/μL)	Porcine Muscle Screening Level (ng/μL)	Porcine Muscle Confirmation Level (ng/μL)
Ciprofloxacin				
Difloxacin	1.25	1.25	1.25	1.25
Enrofloxacin	1.25	1.25	1.25	2.5
Erythromycin A	2.5	2.5	2.5	2.5
Florfenicol	5	5	5	5
Flunixin	0.625	0.625	0.625	0.625
Gamithromycin	2.5	2.5	2.5	10
Lincomycin	2.5	2.5	2.5	2.5
Melengestrol Acetate	0.5		0.5	
Norfloxacin	1.25	1.25	1.25	1.25
Oxyphenylbutazone	2.5	2.5	2.5	5
Oxytetracine	25	25	25	25
Phenylbutazone	2.5	2.5	10	10
Pirlimycin	12.5	12.5	12.5	12.5
Prednisone	2.5	5	2.5	2.5
Ractopamine	0.15	0.3	0.15	0.3
Salbutamol	0.15		0.15	
Sarafloxacin	1.25	1.25	1.25	1.25
Sulfachloropyridazine	2.5	2.5	2.5	2.5
Sulfadiazine	2.5	2.5	5	5
Sulfadimethoxine	2.5	2.5	2.5	2.5
Sulfadoxine	2.5	2.5	2.5	2.5
Sulfaethoxypyridazine	2.5	2.5	2.5	2.5
Sulfamerazine	2.5	2.5	2.5	2.5
Sulfamethazine	2.5	10	2.5	2.5
Sulfamethizole	2.5	2.5	2.5	2.5
Sulfamethoxazole	2.5	2.5	2.5	2.5
Sulfamethoxypyridazine	2.5	2.5	2.5	2.5
Sulfanitran	2.5	10	2.5	5
Sulfapyridine	2.5	2.5	2.5	2.5
Sulfaquinoxaline	2.5	2.5	2.5	2.5

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01		Page 16 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Standard Analyte	Bovine Muscle Screening Level (ng/μL)	Bovine Muscle Confirmation Level (ng/μL)	Porcine Muscle Screening Level (ng/μL)	Porcine Muscle Confirmation Level (ng/μL)
Sulfathiazole	2.5	2.5	2.5	2.5
Tetracycline	25	100	25	50
Tilmicosin	3	3	3	3
Tulathromycin A	50	50	50	50
Tylosin	5	5	5	5

- ii. Prepare the composite working solution(s) for the veterinary drugs contained in the "Beta Lactam mix" spiking solutions using the stock standard solutions above and the volumes listed in the table below.
 - (a) Calculate or use the volume of stock solution required to give the concentrations listed for each standard type in the table below.
 - (b) Pipet the calculated volume of stock into a 10 mL volumetric flask.
 - (c) Dilute to 10 mL volume with water
 - (d) Cap flask and mix
 - (e) Transfer solution into plastic screw cap vials.
 - (f) Composite working solutions will expire in 2 months when stored at ≤ -10 °C or at the time of the earliest expiring component.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 17 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Table 8 – Beta-lactam kidney working solutions preparation volumes

Standard Analyte	Stock or Intermediate Standard Concentration (ng/μL)	Bovine Kidney Screening Volume (μL)	Bovine Kidney Confirmation Volume (μL)	Porcine Kidney Screening Volume (μL)	Porcine Kidney Confirmation Volume (μL)
Ampicillin	1000	10	10	10	10
Cefazolin	1000	100	100	100	100
Cloxacillin	1000	10	10	10	10
DCCD	1000	100	100	100	200
Dicloxacillin	1000	100	100	100	100
Nafcillin	1000	100	100	100	100
Oxacillin	1000	100	100	100	100
Penicillin G	1000	50	50	50	50

Table 9 – Beta-lactam muscle working solutions preparation volumes

Standard Analyte	Stock or Intermediate Standard Concentration (ng/μL)	Bovine Muscle Screening Volume (μL)	Bovine Muscle Confirmation Volume (μL)	Porcine Muscle Screening Volume (μL)	Porcine Muscle Confirmation Volume (μL)
Amoxicillin	350	114		114	
Ampicillin	1000	10	10	10	10
Cefazolin	1000	100	100	100	100
Cloxacillin	1000	10	10	10	10
DCCD	1000	100	400	200	400
Dicloxacillin	1000	100	100	100	100
Nafcillin	1000	100	100	100	100
Oxacillin	1000	100	100	100	100
Penicillin G	1000	50	50	50	50

Using these volumes yield the following concentrations as found in the tables below:

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01		Page 18 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Table 10 – Beta-lactam kidney working solution levels

Standard Analyte	Bovine Kidney Screening Level (ng/μL)	Bovine Kidney Confirmation Level (ng/μL)	Porcine Kidney Screening Level (ng/μL)	Porcine Kidney Confirmation Level (ng/μL)
Ampicillin	1	1	1	1
Cefazolin	10	10	10	10
Cloxacillin	1	1	1	1
DCCD	10	10	10	20
Dicloxacillin	10	10	10	10
Nafcillin	10	10	10	10
Oxacillin	10	10	10	10
Penicillin G	5	5	5	5

Table 11 – Beta-lactam muscle working solution levels

Standard Analyte	Bovine Muscle Screening Level (ng/μL)	Bovine Muscle Confirmation Level (ng/μL)	Porcine Muscle Screening Level (ng/μL)	Porcine Muscle Confirmation Level (ng/μL)
Amoxicillin	4		4	
Ampicillin	1	1	1	1
Cefazolin	10	10	10	10
Cloxacillin	1	1	1	1
DCCD	10	40	20	40
Dicloxacillin	10	10	10	10
Nafcillin	10	10	10	10
Oxacillin	10	10	10	10
Penicillin G	5	5	5	5

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 19 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

- iii. Prepare the composite (spiking) solution for the isotopically-labeled veterinary drugs used for internal standards at 20 ppm if internal standards will be used. This solution can contain any number of the internal standards listed in the method.
 - (a) Calculate the volume of stock solution required to give the 20 ppm level (see the following table).
 - (b) Pipet the calculated volume of stock into a 5 mL volumetric flask.
 - (c) Dilute to 5 mL with acetonitrile.
 - (d) Cap flask and mix.
 - (e) Transfer 1.25 mL of solution into 4 amber glass LC vials.
 - (f) Store at < -10 °C. The stability of the solution is dependent on which internal standards are present. If the solution contains D7 Penicillin G, then the solution is stable for 2 months. Otherwise, the solution is stable for 6 months.

Table 12 – IS spiking solution

Standard Analyte	Stock or Intermediate Standard Concentration (ng/μL)	Volume (μL)	Final Volume (mL)	Working Standard Concentration (ng/μL)
13C6 Sulfamethazine Phenyl	1000	100	5	20
d7 Penicillin G	500	200	5	20
Flunixin-d3	1000	100	5	20

3. Preparation of External Calibration Curve (Optional)

Use the following table to prepare external standards.

Table 13 – Preparation of external standards

Target Conc.	Volume Acetonitrile Standard Mix (μL)	Volume Beta Lactam Standard Mix (μL)	Volume Internal Standard Mix (μL)	Volume Acetonitrile (μL)	Volume of 0.1% Formic Acid in water
0 X	0	0	20	120	860
1/4 X	10	10	20	110	850
1/2 X	20	20	20	100	840
1 X	40	40	20	80	820
2 X	80	80	20	40	780
3 X	120	120	20	0	740

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 20 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Note: The volume of 0.1% Formic Acid in water will be increased by 20 μ L to maintain a total volume of 1 mL if no internal standard solution is added.

E. SAMPLE PREPARATION

Samples collected fresh must be kept cold before and during shipping to the laboratory. Once received at the laboratory, samples must be frozen (< -10 °C) prior to grinding if they cannot be prepared on the day of receipt. Once frozen, the sample should be allowed to thaw, while keeping it as cold as possible. Dissect away fat and connective tissue. Grind tissue in blender or vertical cutter-mixer until homogeneous. Store samples frozen (< -10 °C) prior to analysis.

F. ANALYTICAL PROCEDURE

1. Preparation of Controls and Samples

- a. Weigh 2 ± 0.1 g of homogenized samples into labeled 50 mL polypropylene centrifuge tubes.

Note: Use corresponding blank tissue for controls for each specific species and tissue sample being analyzed.

- i. Screening - Weigh one portion each for a blank (negative control), a 1/2 X recovery (positive control), a matrix matched standard, and a check sample, if necessary.
- ii. Confirmation - Weigh six 2 g portions of blank tissue into 50 ml polypropylene centrifuge tubes. One for 1/2 X recovery (positive control), one for the blank (negative control) and four for the matrix matched standards (1/2X, 1X, 2X, and 3X). Weigh one additional portion for a check sample, if necessary.
- iii. Prepare recoveries, check samples, blank, matrix matched standard(s), and samples using the solutions and volumes in the table below:

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01		Page 21 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Table 14 – Preparation of controls and samples

Sample Type	Acetonitrile Standard Mix (μL)	Beta Lactam Standard Mix (μL)	Internal Standard Mix (μL)	Acetonitrile (μL)	Water (μL)
Samples and Negative Controls			40	160	160
Matrix Matched Standards				200	160
½ X Recovery	40	40	40	120	120
1X Recovery	80	80	40	80	80
2X Recovery	160	160	40		

Note: If no internal standard mix is used, the volume of acetonitrile added to each tube must be increased by 40 μL.

Note: The “1/2 X Recovery” sample type represents the minimum level of applicability (refer to Section J.3).

2. Extraction Procedure

- a. Vortex all uncapped tubes 10 seconds each to mix chemicals with matrix and allow to stand 5 minutes.
Note: Press down to start swirling action slowly and then increase speed to mix without splashing.
- b. Add 9.7 mL of 4/1 (v/v) acetonitrile/water to all tubes using a calibrated solvent dispenser. Cap tubes well.
- c. Place the racks of tubes in platform shaker on high for 5 minutes.
- d. Centrifuge the tubes at >3000 rcf for 5 minutes.
- e. Decant each extract into a pre-labeled 50 mL polypropylene centrifuge tubes containing 0.50 g of C18.
- f. Add 10 mL of hexane (saturated with acetonitrile) using a calibrated dispenser to all tubes. Cap all tubes well.
- g. Shake all tubes in the platform shaker on high for 1 minute.
- h. Centrifuge all tubes at approximately 3000 rcf for 5 minutes.
- i. Aspirate hexane to waste using a Pasteur pipette on a hose connected to a side-arm Erlenmeyer flask.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 22 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

- j. Pipet a 5 mL aliquot of the extract into a pre-labeled 15 mL polypropylene centrifuge tubes.
- k. Evaporate extract to ≤ 0.5 mL in TurboVap at $45 \pm 2^\circ\text{C}$ at 15 psi setting, changing to 20 psi after 10 minutes. Total time to evaporate to ≤ 0.5 mL is 45 to 60 minutes.
- l. Per the table below, add the volumes of the standard mixes and solutions to the matrix matched standards, add 140 μL acetonitrile to all other tubes, and dilute to 1.0 mL with 0.1% formic acid in water. Vortex all tubes for 5 seconds.

Note: After dilution to 1.0 mL with 0.1% formic acid in water, the extracts will contain $\approx 14\%$ acetonitrile by volume.

Table 15 – Preparation of matrix matched standards

Sample Type	Acetonitrile Standard Mix (μL)	Beta Lactam Standard Mix (μL)	Internal Standard Mix (μL)	Acetonitrile (μL)
Samples, Recoveries, and Negative Controls				140
1/2X Matrix Matched Standard	20	20	20	100
1X Matrix Matched Standard	40	40	20	80
2X Matrix Matched Standard	80	80	20	40
3X Matrix Matched Standard	120	120	20	0

Note: If no internal standard mix is used, the volume of acetonitrile added to each tube must be increased by 20 μL .

- m. Pipet 500 μL of all final extracts and calibration standards into bottom portions of pre-labeled Whatman Mini UniPrep Syringless Filter Vials, PVDF, 0.2 micron.
 - n. Pressed filter/caps onto the vials and placed in autosampler tray for UHPLC-MS/MS analysis.
3. Instrumental Settings
- Note: The instrument parameters may be optimized to ensure system suitability.*
- a. Instrument Operating Parameters - UHPLC system
 - i. Mobile phase for Residue analysis:
 - Mobile Phase A - 95% water / 5% ACN / 0.1% Formic Acid
 - Mobile Phase B - 100% ACN / 0.1% Formic Acid

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 23 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

Flush column with 1:1 A/B at a flow rate of 0.5 mL/min for three minutes. Change the mobile phase initial conditions to 100% A. Allow column to equilibrate until the "delta" value on the pressure reading is < 20.

ii. UHPLC gradient program:

Flow rate: 0.5 mL/min
Pressure Limits: 200 psi minimum; 15,000 psi maximum
Run time: 12.9 minutes

Table 16 – UHPLC gradient program

<i>Time (min)</i>	<i>% Mobile Phase A</i>	<i>% Mobile Phase B</i>	<i>Gradient</i>
0.00	100	0	none
0.10	80	20	linear
8.00	80	20	linear
9.50	60	40	linear
9.60	10	90	linear
12.90	10	90	linear

b. Autosampler program:

- i. Run time: 12.90 minutes
- ii. Injection loop: 20 µL
- iii. Loop option: Full Loop
- iv. Injection volume: 20 µL
- v. Weak wash solvent: 10/90 Acetonitrile/Water
- vi. Weak wash volume: 1200 µL
- vii. Strong wash solvent: 50/50 Acetonitrile/Water
- viii. Strong wash volume: 400 µL
- ix. Sample temperature: 10 °C

c. Column manager

- i. Column valve position: To match column location.
- ii. Column manager temperature: 40 °C
- iii. Use divert valve to divert eluant to waste 0.25 minutes prior to first peak and 0.25 minutes after last analyte peak.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 24 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

d. Instrument Operating Parameters - Mass Spectrometer

Mass Spectrometer calibration and resolution are to be done according to the manufacturer's specification using the manufacturer's supplied calibration solution.

- i. Type: MS/MS
- ii. Electrospray Source Parameters:
 - (a) Capillary (kV): 3.0
 - (b) Cone (V): Variable - analyte dependent
 - (c) Extractor (V): 3.0
 - (d) RF (V): 0.10
 - (e) Source Temperature (°C): 150
 - (f) Desolvation Temperature (°C): 450
 - (g) Cone Gas Flow (L/hr): 20
 - (h) Desolvation Gas Flow (L/hr): 900
 - (i) Collision Gas Flow (mL/min): 0.10
- iii. Analyzer Parameters:
 - (a) LM 1 Resolution: 10.6
 - (b) HM 1 Resolution: 14.8
 - (c) MSMS Mode Entrance: 1
 - (d) MSMS Mode Collision Energy: Variable – analyte dependent
 - (e) MSMS Mode Exit: 0.5
 - (f) LM 2 Resolution: 9.5
 - (g) HM 2 Resolution: 15.8
- iv. MS Method Parameters:
 - (a) Type: MRM
 - (b) Ion Mode: ES+
 - (c) Dwell (s): 0.005
 - (d) Start time (min): 0.4
 - (e) End time (min): 6.1

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 25 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

v. Instrumental Settings - Scan Parameters

Table 17 – Instrument settings

Analyte	RT (min)	window (min)	Dwell Times (ms)	Precursor Ion(m/z)	Product Ions (m/z)	Collision Energy (V)	Cone (V)
Sulfanilamide	1.06	0.42 - 2.00	5	172.87	92.88 75.90 65.90	36 22 20	44
Amoxicillin	1.43	1.25 – 1.90	5	366.14	114.00 349.34 208.09	22 10 14	20
Salbutamol	1.43	1.27 - 1.67	5	240.23	148.18 222.30 166.15	20 10 15	20
Cimaterol	1.47	1.32 - 1.72	5	220.04	142.97 115.88 88.91	24 16 10	20
DCCD	1.66	1.55 - 2.10	5	549.10	182.96 241.12 125.91	30 20 76	40
Lincomycin	1.82	1.69 - 2.09	5	407.31	126.10 359.22 389.21	35 20 15	40
Sulfadiazine	1.91	1.79 - 2.19	5	251.10	156.08 108.02 158.08	15 20 15	30
Ampicillin	1.95	1.82 - 2.35	5	350.14	106.07 114.00 160.07	24 30 24	26
Desethylene Ciprofloxacin	1.98	1.88 - 2.30	5	306.20	288.19 245.18 289.28	20 20 10	35
Sulfathiazole	2.00	1.90 - 2.30	5	256.07	156.08 108.04 101.03	15 25 25	25
Sulfapyridine	2.09	1.92 - 2.37	5	250.09	156.08 92.00 108.11	18 26 28	32
Norfloxacin	2.07	1.95 - 2.43	5	320.17	276.19 233.10 219.03	20 30 20	35
Tulathromycin A	2.08	1.93 - 2.38	5	806.75	72.01 577.47 116.04	35 25 30	40
Oxytetracycline	2.13	1.94 - 2.40	5	461.23	426.35 443.39 201.24	20 15 25	25

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 26 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS	
Revision: .01	Replaces: CLG-MRM1.00
Effective: 8/06/2012	

Analyte	RT (min)	window (min)	Dwell Times (ms)	Precursor Ion(m/z)	Product Ions (m/z)	Collision Energy (V)	Cone (V)
Ciprofloxacin	2.13	2.00 - 2.45	5	332.23	245.18 288.39 203.21	25 20 40	35
Ractopamine	2.19	2.08 - 2.48	5	302.20	164.03 106.98 120.98	15 30 30	20
Sulfamerazine	2.21	2.10 - 2.49	5	265.04	91.88 155.92 107.91	30 15 15	35
Danofloxacin	2.22	2.12 - 2.52	5	358.12	96.03 314.18 283.14	30 15 25	30
Tetracycline	2.27	2.16 - 2.56	5	445.23	154.13 410.22 427.32	30 20 15	30
Enrofloxacin	2.30	2.20 - 2.60	5	360.24	316.40 245.25 203.16	20 25 40	35
Sulfamethizole	2.44	2.30 - 2.70	5	271.06	156.08 91.97 107.98	16 30 26	28
Sulfamethazine	2.46	2.26 - 2.65	5	279.14	186.14 156.08 108.12	20 20 25	35
Sulfamethazine-6C13	2.46	2.24 - 2.73	5	285.17	186.06 98.02 124.05	18 32 26	32
Cefazolin	2.44	2.32 - 2.80	5	455.13	156.03 323.15 112.08	16 12 35	20
Sulfamethoxy pyridazine	2.48	2.34 - 2.72	5	281.12	156.09 126.17 108.06	20 20 25	30
Difloxacin	2.53	2.27 - 2.75	5	400.25	356.43 299.23 285.28	20 30 40	35
Sarafloxacin	2.48	2.35 - 2.75	5	386.09	342.18 299.17 270.20	20 25 45	40
Pirlimycin	2.67	2.44 - 2.84	5	411.29	112.20 363.33 110.34	40 20 45	30
Chlortetracycline	2.76	2.63 - 3.03	5	479.19	154.08 444.26 462.20	30 20 20	30

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 27 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS	
Revision: .01	Replaces: CLG-MRM1.00
Effective: 8/06/2012	

Analyte	RT (min)	window (min)	Dwell Times (ms)	Precursor Ion(m/z)	Product Ions (m/z)	Collision Energy (V)	Cone (V)
Clindamycin	2.81	2.68 - 3.08	5	425.30	126.22 377.40 124.32	40 20 45	45
Gamithromycin	2.85	2.70 - 3.10	5	777.79	83.02 115.97 158.10	48 36 46	66
Sulfachloropyridazine	2.84	2.29 - 3.00	5	285.00	156.08 107.98 91.97	16 26 30	28
Tilmicosin	3.00	2.83 - 3.25	5	869.79	174.16 132.23 696.63	35 35 35	45
Sulfadoxine	3.00	2.85 - 3.25	5	311.15	156.08 108.14 140.09	20 30 30	35
Sulfamethoxazole	3.00	2.89 - 3.29	5	254.02	92.06 155.97 107.89	25 20 20	30
Sulfaethoxypyridazine	3.03	2.91 - 3.31	5	295.13	156.09 140.18 108.09	20 20 25	30
Florfenicol	3.05	2.93 - 3.33	5	358.10	241.01 206.00 130.36	18 22 60	20
Chloramphenicol	3.25	3.00 - 3.60	5	323.11	274.97 164.99 118.78	15 30 40	15
Erythromycin A	3.45	3.28 - 3.68	5	734.75	158.15 115.89 576.48	30 40 20	30
Sulfadimethoxine	3.46	3.33 - 3.75	5	311.14	156.10 108.04 245.22	20 30 20	35
Sulfaquinoxaline	3.46	3.35 - 3.75	5	301.10	156.13 107.98 91.97	18 28 36	34
Prednisone	3.56	3.37 - 3.85	5	359.15	341.11 146.94 267.28	10 26 15	22
Tylosin	3.63	3.46 - 3.86	5	916.76	174.20 101.08 145.15	35 35 35	45
Pen G-d7	3.74	3.63 - 4.20	5	342.13	183.10 160.10 98.10	26 24 54	46

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 28 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS	
Revision: .01	Replaces: CLG-MRM1.00
Effective: 8/06/2012	

Analyte	RT (min)	window (min)	Dwell Times (ms)	Precursor Ion(m/z)	Product Ions (m/z)	Collision Energy (V)	Cone (V)
Penicillin G	3.77	3.61 - 4.03	5	335.13	176.02 160.10 114.01	16 18 30	18
Beta/Dexamethasone	4.01	3.89 - 4.29	5	393.22	373.24 147.07 355.29	10 28 14	20
Sulfanitran	4.05	3.93 - 4.33	5	336.15	156.08 134.13 92.71	14 28 38	30
Zeranol (B-Zearalanol)	4.26	4.13 - 4.53	5	323.21	305.20 189.09 149.02	10 24 30	16
Oxacillin	4.30	4.17 - 4.70	5	402.14	160.00 243.08 144.06	20 18 34	22
Cloxacillin	4.56	4.43 - 5.05	5	436.16	160.06 277.14 114.07	12 16 44	22
Nafcillin	4.69	4.55 - 5.10	5	415.22	199.17 171.07 115.11	16 42 70	20
Oxyphenylbutazone	4.74	4.59 - 4.99	5	325.18	120.06 148.16 204.12	24 30 16	26
Flunixin	4.82	4.61 - 5.03	5	297.05	279.05 109.00 264.04	22 50 32	42
Flunixin-d3	4.81	4.63 - 5.03	5	300.05	282.12 112.03 264.04	24 54 36	40
Dicloxacillin	4.94	4.80 - 5.35	5	470.19	160.06 311.11 114.07	14 16 48	22
Phenylbutazone	5.86	5.69 - 6.09	5	309.12	119.95 76.90 91.80	22 24 18	32
Melengesterol acetate	6.20	6.04 - 6.44	5	397.35	279.33 337.46 221.28	20 15 40	30

Note: Product ions are listed with the expected screening ion in bold (top) followed by diagnostic ions 1 and 2 (middle and bottom).

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 29 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

4. Sample Set
 - a. Screening Set
 - i. External Standard(s) (optional)
 - ii. Matrix matched standard
 - iii. Recovery(ies) (positive controls)
 - iv. Check sample (if necessary)
 - v. Blank (negative control)
 - vi. Up to 27 Samples
 - vii. External standard, matrix matched standard, or recovery
 - b. Confirmation Set
 - i. External Standard(s) (optional)
 - ii. Matrix matched standards
 - iii. Recovery(ies) (positive control)
 - iv. Check sample (if necessary)
 - v. Matrix Matched Blank (negative control)
 - vi. Up to 24 Samples
 - vii. External standard, matrix matched standard, or recovery

Note: Placing solvent blanks in the sample injection sequence is prudent in case a high finding leads to carry-over. Additionally, one may want to include an additional external standard, matrix matched standard, or recovery within the sample injection sequence to verify retention time and instrument response stability.

G. CALCULATIONS / IDENTIFICATION

1. Screening
 - a. The screening ion for a given analyte must be present. The required ion for each compound is listed in Table 13.
 - b. The retention times for the screening ion in the fortified recoveries must match the retention time of the screening ion in the matrix-matched standard within 5%. Retention time for the screening ions in the samples must match the retention time of the screening ions in a fortified recovery or the matrix matched standard within 5%.
 - c. The screening ion must have a signal-to-noise ratio ≥ 3 . This may be verified by visual inspection.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 30 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

- d. A sample is screened positive for an analyte if the following criteria are met:
 - i. The fortified recovery of the analyte must exceed 10% of the 1/2 X matrix matched standard level.
 - ii. The sample response equals or exceeds the 1/2 X (or level of interest) fortified recovery level.
- e. The level of the screening ion in the blank (negative control) must be less than 10% of the 1/2 X level for the matrix-matched standard.

Note: If a sample shows a positive response for a compound which did not meet screening criteria in the associated QC samples, then further testing of that sample is warranted.

2. Confirmation

- a. Monitored ions for each analyte will be assessed as follows:
 - i. Recovery retention times must match the retention time of the matrix matched standard within 5%. Retention time for the samples must match the retention time of the positive control or the matrix matched standard within 5%.
 - ii. All product ions specified for ratio matching are present with a signal-to-noise ratio ≥ 3 . This may be verified by visual inspection.
 - iii. One of the following ion ratio matching conditions is met:

Note: Ratios are calculated by dividing the area count of each diagnostic ion by the area count of the base ion. Ion ratios should be less than 1. If the ratio is not less than 1 for a sample set, the inverse of this ratio may be used.

 - (a) If two product ions are assessed, one sample ion ratio should match the calculated average ratio of the matrix-matched standards within a $\pm 10\%$ absolute difference.
 - (b) If three product ions are monitored, the presence of two sample ratios should match the calculated average ratio of the matrix-matched standards within a $\pm 20\%$ absolute difference.
 - iv. The fortified recovery of the analyte must exceed 10% of the 1/2 X matrix matched standard level.
 - v. The blank (negative control) must be less than 10% of the 1/2 X level for the matrix matched standard.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 31 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

H. SAFETY INFORMATION AND PRECAUTIONS

1. Required Protective Equipment -Safety eyewear, protective gloves, and lab coat.
2. Hazards

<i>Procedure Step</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
Antibiotic standards	Some individuals may have allergic reactions to certain β -lactams, sulfa, or other drugs.	Wear appropriate personal protective equipment to avoid dermal contact.
Acetonitrile, Methanol	Flammable	Keep in well-closed containers away from ignition sources. Avoid contact or prolonged exposure to vapors. Work in fume hood. Keep away from flame or heat.
Formic acid	Corrosive, Caustic	Wear personal protective equipment, avoid skin contact.

3. Disposal Procedures
Follow local, state and federal guidelines for disposal.

I. QUALITY ASSURANCE PLAN

1. Performance Standard
 - a. Screening Criteria
 - i. For set acceptance, 90% of the monitored analytes in the fortified recovery (positive control) must meet screening criteria. For sample reporting purposes, the analytes of interest in the fortified recovery (positive control) must meet screening criteria.
 - ii. The blank (negative control) must be negative using the criteria in Section G.
 - b. Confirmation Criteria
 - i. For set acceptance, nine of the following ten analytes must meet confirmation criteria (Beta-Dexamethasone, Desethylene Ciprofloxacin, Erythromycin A, Florfenicol, Flunixin, Oxytetracycline, Penicillin G, Prednisone, Sulfadimethoxine, and Sulfamethazine). For sample

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 32 of 36	
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

reporting purposes, the analytes of interest in the fortified recovery (positive control) must meet confirmation criteria.

- ii. The blank (negative control) must be negative using the criteria in Section G for the analytes of interest.

2. Critical Control Points and Specifications

Record

Acceptable Control

none known at this time

3. Intralaboratory Check Samples

a. System, minimum contents.

- i. Frequency: One per week per analyst when samples analyzed.
- ii. Records are to be maintained.

b. Acceptability criteria.

Refer to I. 1.

If unacceptable values are obtained, then:

- i. Investigate following established procedures.
- ii. Take corrective action as warranted.

4. Sample Condition upon Receipt

Cool or frozen

J. APPENDIX

1. References

[Reserved]

2. Chromatograms/spectra

[Reserved]

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01		Page 33 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS		
Revision: .01	Replaces: CLG-MRM1.00	Effective: 8/06/2012

3. Minimum Level of Applicability

Table 18 - Minimum Level of Applicability for kidney

Analyte (for kidney) <i>Level (ng/g)</i>	Screen		Confirmation	
	Bovine	Porcine	Bovine	Porcine
Ampicillin	20	20	20	20
Beta/Dexamethasone	50	50	50	50
Cefazolin	200	200	200	200
Chloramphenicol	6	6		
Chlortetracycline	1000	1000	1000	1000
Cimaterol	12	3	12	3
Ciprofloxacin	25	25	50	25
Clindamycin	50	50	50	50
Cloxacillin	20	20	20	20
Danofloxacin	25	25	25	25
DCCD (marker for Ceftiofur)	200	200	200	400
Desethylene Ciprofloxacin	25	25	25	25
Dicloxacillin	200	200	200	200
Difloxacin	25	25	25	25
Enrofloxacin	25	25	50	25
Erythromycin A	50	50	50	50
Florfenicol	100	100	100	100
Flunixin	12.5	12.5	12.5	12.5
Gamithromycin	50	50	100	100
Lincomycin	50	50	50	50
Nafcillin	200	200	200	200
Norfloxacin	25	25	25	25
Oxacillin	200	200	200	200
Oxyphenylbutazone		50		
Oxytetracycline	500	500	500	500
Penicillin G	100	100	100	100
Phenylbutazone		50		50
Pirlimycin	250	250	250	250
Prednisone	50	50	50	50
Ractopamine	3	3	12	
Salbutamol	6	3		

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 34 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS	
Revision: .01	Replaces: CLG-MRM1.00
Effective: 8/06/2012	

Analyte (for kidney)	Screen		Confirmation		
	Level (ng/g)	Bovine	Porcine	Bovine	Porcine
Sarafloxacin		25	25	25	25
Sulfachloropyridazine		50	50	50	200
Sulfadiazine		50	50	50	50
Sulfadimethoxine		50	50	50	50
Sulfadoxine		50	50	50	50
Sulfaethoxypyridazine		50	50	50	50
Sulfamerazine		50	50	50	50
Sulfamethazine		50	50	50	50
Sulfamethizole		50	50	50	50
Sulfamethoxazole		50	50	100	50
Sulfamethoxypyridazine		50	50	50	50
Sulfanilamide			100		200
Sulfanitran		50	50	100	50
Sulfapyridine		50	50	50	50
Sulfaquinoxaline		50	50	50	50
Sulfathiazole		50	50	50	50
Tetracycline		500	500	2000	500
Tilmicosin		120	240		
Tylosin		100	200	100	200
Zearalanol			12		12

Table 19 - Minimum Level of Applicability for muscle

Analyte (for muscle)	Screen		Confirmation		
	Level (ng/g)	Bovine	Porcine	Bovine	Porcine
Amoxicillin		80	80		
Ampicillin		20	20	20	20
Betamethasone		50	50	50	50
Cefazolin		200	200	200	200
Chloramphenicol		6	6		
Chlortetracycline		1000	1000	1000	1000
Cimaterol		3	3		
Ciprofloxacin		25	25	50	25
Clindamycin		50	50	50	50

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 35 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS	
Revision: .01	Replaces: CLG-MRM1.00
Effective: 8/06/2012	

Analyte (for muscle)	Screen		Confirmation	
	Level (ng/g)	Bovine	Porcine	Bovine
Cloxacillin	20	20	20	20
Danofloxacin	25	25	25	25
DCCD	200	400	800	800
Desethylene Ciprofloxacin	25	25	25	25
Dicloxacillin	200	200	200	200
Difloxacin	25	25	25	25
Enrofloxacin	25	25	25	50
Erythromycin A	50	50	50	50
Florfenicol	100	100	100	100
Flunixin	12.5	12.5	12.5	12.5
Gamithromycin	50	50	50	200
Lincomycin	50	50	50	50
Melengestrol Acetate	10	10		
Nafcillin	200	200	200	200
Norfloxacin	25	25	25	25
Oxacillin	200	200	200	200
Oxyphenylbutazone	50	50	50	100
Oxytetracycline	500	500	500	500
Penicillin G	100	100	100	100
Phenylbutazone	50	200	50	200
Pirlimycin	250	250	250	250
Prednisone	50	50	100	50
Ractopamine	3	3	6	6
Salbutamol	3	3		
Sarafloxacin	25	25	25	25
Sulfachloropyridazine	50	50	50	50
Sulfadiazine	50	100	50	100
Sulfadimethoxine	50	50	50	50
Sulfadoxine	50	50	50	50
Sulfaethoxypyridazine	50	50	50	50
Sulfamerazine	50	50	50	50
Sulfamethazine	50	50	200	50
Sulfamethizole	50	50	50	50
Sulfamethoxazole	50	50	50	50

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

CLG-MRM1.01	Page 36 of 36
Title: Screening and Confirmation of Animal Drug Residues by UHPLC-MS-MS	
Revision: .01	Replaces: CLG-MRM1.00
Effective: 8/06/2012	

Analyte (for muscle)	Screen		Confirmation	
	Bovine	Porcine	Bovine	Porcine
<i>Level (ng/g)</i>				
Sulfamethoxypyridazine	50	50	50	50
Sulfanitran	50	50	200	100
Sulfapyridine	50	50	50	50
Sulfaquinoxaline	50	50	50	50
Sulfathiazole	50	50	50	50
Tetracycline	500	500	2000	1000
Tilmicosin	60	60	60	60
Tulathromycin A	1000	1000	1000	1000
Tylosin	100	100	100	100

K. APPROVALS AND AUTHORITIES

1. Approvals on file.
2. Issuing Authority: Director, Laboratory Quality Assurance Division.