

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 1 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

Contents

A.	INTRODUCTION	2
B.	EQUIPMENT.....	3
C.	REAGENTS AND SOLUTIONS.....	4
D.	STANDARDS.....	5
E.	SAMPLE PREPARATION.....	8
F.	ANALYTICAL PROCEDURE	8
G.	CALCULATIONS	12
H.	HAZARD ANALYSIS.....	13
I.	QUALITY ASSURANCE PLAN.....	15
J.	WORKSHEET	17
K.	APPENDIX.....	19

United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science

SOP No: CLG- RAC1.01	Page 2 of 22	
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

A. INTRODUCTION

1. Theory

Ractopamine, a β -adrenergic agonist, is extracted from swine and bovine liver and muscle with methanol. An aliquot of the extract is evaporated, borate buffer is added, and ractopamine is extracted into ethyl acetate by liquid/liquid partition. The ethyl acetate extract is further purified by passing it through an acidic alumina solid phase extraction column. Ractopamine is eluted from the column with methanol. The methanol extract is evaporated to dryness, then dissolved in dilute acetic acid, filtered and analyzed for ractopamine using high performance liquid chromatography (HPLC) with fluorescence detection.


2. Applicability

This method has been validated for the determination of ractopamine in swine and bovine liver at 25 - 300 ppb and in swine and bovine muscle tissues at 1.0 - 300 ppb.


3. Structure

Ractopamine HCl is a mixture of four stereoisomers in approximately equal proportions (RS, SR, RR, and SS). This HPLC method does not distinguish between these stereoisomers, and thus results in a single peak for all four stereoisomers present in a particular sample.

Structure of ractopamine HCl.



Structure of Ritodrine



**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 3 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

B. EQUIPMENT

1. Apparatus

Note: Equivalent instrument and apparatus may be substituted.

- a. HPLC pump - Agilent 1100 series quaternary pump, G1311A.
- b. HPLC autosampler - Agilent 1100 series ALS, G1313A.
- c. HPLC variable wavelength fluorescence detector - Agilent FLD, G1321A.
- d. HPLC detector recording device - Agilent Chem Station, Agilent Technologies, Wilmington, Delaware.
- e. HPLC mobile phase filtering and degassing apparatus - Millipore Microfiltration Assembly, 47 mm, filtering with a type GV 0.2 micron filter.
- f. Analytical balance (± 0.0001 g) - Mettler Toledo AG204.
- g. Top loading balance (± 0.01 g) - Mettler Toledo PG5002.S Delta Range.
- h. Magnetic stirrer and stir bars Corning Stirring/Hot Plate - PC 420, and VWR, 5/16" Diameter x 2" stir bars, 58949-038.
- i. Meat grinder or food processor - Rival model 2100 M/2.
- j. Branson Sonifier 450 ultrasonic generator with a 1/4 inch micro tip or a mechanical blender such as an UltraTurrax No.T25 equipped with an S25N-10G dispersing rotor.
- k. Vortex-2 Genie test tube vortexer - Scientific Industries, Bohemia, NY.
- l. Polypropylene centrifuge tubes - 50 mL conical, with closures, Blue Max 352070, 30 x 115 mm, 25/pk, Becton Dickerson, Franklin Lakes, NJ.
- m. Centrifuge - IEC Centra-8R centrifuge.
- n. Volumetric flasks - 50 and 100 mL volumetric flasks.
- o. Micropipettes - 500 - 3000 μ L.
- p. Test tubes - 16 x 100 mm test tubes with closures - CMS, Inc., Cat. Nos. 339-309 and 270-671, respectively (Fisherbrand Disposable Culture Tubes, Borosilicate glass, 16 x 125 mm, catalog #14-961-30).
- q. Volumetric pipettes - 10 mL and 15 mL, Class A.
- r. Serological pipettes - 10 mL, disposable, polystyrene, Falcon 357530.
- s. Mixing cylinders - 100 mL graduated mixing cylinders, stoppered, CMS, Inc., Cat.

United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science

SOP No: CLG- RAC1.01		Page 4 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

No. 101-527.

- t. Glass bottles - 1000 mL glass bottles.
- u. Pasteur pipettes - Disposable glass Pasteur pipettes.
- v. Membrane filters - Disposable 13 mm PVDF Membrane 0.22 µm filters, HPLC Certified Minispike Outlet, Gelman Acrodisc PM4450T.
- w. Small disposable syringes - Becton-Dickerson, 1 mL syringe, 309602.
- x. Spatulas - Metal or Teflon coated.
- y. N-Evap air/nitrogen stream evaporator - Meyer N-Evap Analytical Evaporator, Organomation Associates, Inc. South Berlin, MA.
- z. Vacuum apparatus for solid phase cartridges, or syringes for sample application to cartridges - Supelco Visiprep, Supelco, 595 North Harrison Road, Bellefonte, PA.
- aa. Ultrasonic water bath - Branson model 2200, 125 watts.
- bb. pH meter - Orion 611, ThermoOrion, 500 Cummings Center, Beverly, MA.
- cc. Solid Phase Extraction (SPE) Cartridges - Acidic alumina (Alumina A) approximately 2 g packing, Activity Grade I, Waters Sep-Pak Classic, No. 51800, Waters Sep-Pak Plus. No. 20500.

Note: Not all brands of acidic alumina produce acceptable results. See Section K.1, Appendix, SPE CARTRIDGE TESTING.
- dd. HPLC Column - C18 reversed phase such as Supelcosil LC-18-DB, 5 µm, 4.6 x 250 mm, Part No. 58355.

C. REAGENTS AND SOLUTIONS

1. Reagents

Note: An equivalent solution or reagent may be substituted.

- a. Water (H₂O) - HPLC grade or distilled, deionized.
- b. Methanol (MeOH) - HPLC grade, Fisher A452-4.
- c. Acetonitrile (ACN) - HPLC grade, Fisher A998-4.
- d. Ethyl acetate - HPLC grade, Fisher E195-4.
- e. Glacial acetic acid - Reagent grade, Fisher A38-212.
- f. Sodium hydroxide (NaOH) - Reagent grade, Fisher S318-1.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 5 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

- g. 1-Pentanesulfonic acid sodium salt - HPLC grade, Baker 2841-05.
- h. Sodium borate decahydrate - Reagent grade, Baker 3570-01.

2. Solutions

- a. 1 N sodium hydroxide solution:

Dissolve 40 g sodium hydroxide in approximately 800 mL deionized water. Mix well. Cool to room temperature and dilute to 1 L.

Note: Exothermic reaction. Recommend the flask be placed in an ice bath when initially dissolving the sodium hydroxide.

- b. Borate Buffer (0.025M, pH 10.3):

Add 9.54 ± 0.05 g sodium borate to 900 mL of HPLC grade water in a graduated cylinder or glass bottle, and dissolve by mixing, add 1 N sodium hydroxide (approximately 40 mL) until the pH is 10.3 ± 0.1 . Dilute to 1 L with HPLC water. Store at room temperature.

Note: Check buffer monthly. Acceptable pH range is 9.5 -11.0.

- c. HPLC Mobile Phase:

Add 320 mL HPLC acetonitrile to 680 mL HPLC water. Then add 20 mL of glacial acetic acid and 0.87 ± 0.05 g 1-pentanesulfonic acid. Mix well, filter through a 0.45 μ m filter if necessary, and degas.

- d. Sample Diluent (acetic acid - 2 percent v/v):

Add 20 mL reagent grade glacial acetic acid to 980 mL HPLC water and mix well.

- e. Mobile phase, for alternative HPLC column:

Prepare as directed in section C.2.c., except use 250 mL of acetonitrile and 750 mL of water.

D. STANDARDS

1. Source

Note: The analyst may prepare different standard volumes and/or concentrations to cover the range of interest.

- a. Ractopamine HCl Reference Standard: Eli Lilly and Company, Indianapolis, Indiana. Store at room temperature.
- b. Ritodrine HCl Reference Standard: US Pharmacopeia Convention, 12601

United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science

SOP No: CLG- RAC1.01		Page 6 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

Twinbrook , Rockville, Maryland. Catalog Number, 1604701. Store at room temperature.

2. Reference Standard Preparation

a. Stock solution (1.00 mg/mL):

Note: The stock solution must be adjusted for purity during preparation.

Prepare a ractopamine hydrochloride standard stock solution by adding 100 ± 1 mg of ractopamine hydrochloride reference standard to a 100 mL volumetric flask and diluting to volume with methanol. This solution is stable for three months at 2 - 8°C.

CAUTION: Wear gloves when handling reference standard. Do not inhale the dust of the primary reference standard.

b. Intermediate standard (10 µg/mL):

Pipet 1.0 mL of standard stock solution into a 100 mL volumetric flask and dilute to volume with sample diluent. Mix well. This solution is stable for one month at 2 - 8°C.

c. Prepare the following standard solutions for the analysis of liver at 25 - 300 ppb:

i. Fortification standard (1.5 µg/mL):

Pipet 15 mL of 10 µg/mL intermediate standard solution into a 100 mL volumetric flask and dilute to volume with sample diluent. Mix well. This solution is stable one month at 2 - 8 °C.

ii. Ractopamine HCl external standard curve solutions (25, 50, 75, 150 and 300 ng/mL):

Prepare volumetric dilutions of the 10 µg/mL intermediate standard using sample diluent. For 25, 50, 75, 150 and 300 ng/mL solutions make 250 µL to 100 mL, 500 µL to 100 mL, 750 µL to 100 mL, 1.5 mL to 100 mL and 3.0 mL to 100 mL volumetric dilutions, respectively. These solutions are stable for one month at 2 - 8 °C.

d. Prepare the following standard solutions for the analysis of muscle at 0.5 - 2.0 ppb:

i. Intermediate Standard (100 ng/mL):

Pipet 1.0 mL of 10 µg/mL intermediate standard solution into a 100 mL volumetric flask and dilute to volume with sample diluent. Mix well. This solution is stable one month at 2 - 8 °C.

ii. Fortification standard (10 ng/mL):

Pipet 10 mL of 100 ng/mL intermediate standard solution into a 100 mL

United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science

SOP No: CLG- RAC1.01		Page 7 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

volumetric flask and dilute to volume with sample diluent. Mix well. This solution is stable one month at 2 - 8 °C.

- iii. Ractopamine HCl external standard curve solutions (0.5, 0.75, 1.0, 1.5, 2.0 ng/mL):

Prepare volumetric dilutions of the 100 ng/mL intermediate standard using sample diluent. For 0.5, 0.75, 1.0, 1.5, and 2.0 ng/mL solutions make 500 μ L to 100 mL, 750 μ L to 100 mL, 1000 μ L to 100 mL, 1500 μ L to 100 mL, and 2000 μ L to 100 mL volumetric dilutions, respectively. These solutions are stable for one month when maintained at 2 - 8 °C.

3. Resolution Standard Preparation:

Note: Prepare these solutions as necessary to check column resolution.

- a. Ritrodine HCl stock solution (1 mg/mL):

Weigh 50 ± 0.1 mg of ritodrine hydrochloride reference standard into a 50 mL volumetric flask and dilute to volume with methanol.

- b. Ritrodine HCl intermediate standard (10 μ g/mL):

Pipet 1.0 mL of standard stock solution into a 100 mL flask and dilute to volume with sample diluent. Mix well.

- c. For the analysis of liver at 25 - 300 ppb.

- i. Resolution solution, mixed external standard (25 ng/mL):

Pipet 250 μ L of the 10 μ g/mL ritodrine HCl intermediate solution and 250 μ L of the 10 μ g/mL ractopamine HCl intermediate solution into a 100 mL flask and dilute to volume with sample diluent. Mix well.

- d. For the analysis of muscle at 0.5 - 2.0 ppb:

- i. Ritrodine HCl intermediate standard (100 ng/mL):

Pipet 1 mL of 10 μ g/mL intermediate standard solution into a 100 mL volumetric flask and dilute to volume with sample diluent. Mix well.

- ii. Resolution solution, mixed external standard (1 ng/mL):

Pipet 1 mL of the 100 ng/mL ritodrine HCl intermediate solution and 1 mL of the 100 ng/mL ractopamine HCl intermediate solution into a 100 mL flask and dilute to volume with sample diluent. Mix well.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 8 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

E. SAMPLE PREPARATION

1. Preparation and Storage of Tissues

- a. Initial processing includes grinding or blending of the tissues using a food grinder (or cryogenic grinding) to produce homogenous samples. Grind a minimum 500 g sample of tissue when possible.
- b. Store all tissues at < -10 °C when not processing or sub-sampling. Ractopamine has been shown to be stable in frozen tissue for one year.

Note: Extreme care should be taken to make sure all tissue residue containing ractopamine is cleaned from glassware and other laboratory items in contact with samples and standards. It is recommended that disposable items be used whenever possible and that labware used with standards and other sources containing high levels of ractopamine be kept separate from that used to prepare samples in the low ppb range.

F. ANALYTICAL PROCEDURE

1. Tissue Extraction

- a. Weigh 10.0 ± 0.2 g of frozen or partially thawed ground sample tissue into a suitable container such as a 50 mL polypropylene centrifuge tube.

Note: Prepare blank and recovery samples at this time by weighing two 10 g blank tissues as part of the sample set.

For the analysis of liver at 25 - 300 ppb, prepare a 150 ppb recovery by adding 1 mL of 1.5 µg/mL fortification standard (D.2.c.i) to one of the tissue blanks.

For the analysis of muscle at 0.5 - 2.0 ppb, prepare a 1 ppb recovery by adding 1 mL of the 10 ng/mL fortification standard (D.2.d.ii) to one of the tissue blanks.

- b. Add 20 ± 1 mL of methanol to the sample.
- c. Homogenize the tissue slurry for approximately one minute using an ultrasonic cell disrupter equipped with a 1/4 inch micro tip. Alternatively, the tissue may be blended for approximately one minute using a suitable blender to produce a uniform slurry. The probe must be cleaned between samples with methanol and water rinses, and given a final methanol rinse. A detergent may also be used to assist in cleaning mechanical homogenizers (e.g. UltraTurrax/Polytron). Let the sample stand at room temperature for 10 - 15 minutes to enhance solvent contact with tissue.
- d. Centrifuge the tissue slurry at approximately 3000 rpm (approximately 1500 g) for

United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science

SOP No: CLG- RAC1.01		Page 9 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

10 minutes. Exact speed and centrifugal force is not critical provided a good sediment pack is obtained. Refrigeration may be used, but is not necessary.

- e. Decant the supernatant into a 100 mL graduated mixing cylinder or other appropriate graduated glassware.
- f. Add a second 20 mL of methanol to the tissue, vigorously suspend the centrifuged pack with a spatula, centrifuge as in step d above, and add the second supernatant to the first. Combined extracts will be cloudy.
- g. Repeat Step f, adding the third supernatant to the first and second.
- h. Dilute the combined supernatants to 60 ± 1 mL with methanol and mix well.

Note: This is a suitable stopping point. Extracts may be stored for 7 days at 2 - 8 °C.

- i. Pipette 8.0 ± 0.1 mL aliquot of the combined supernatant into a 16 x 100 mm test tube. Evaporate the sample to less than 0.5 mL under air or nitrogen at 49 ± 2 °C.

Note: Dry the sample aliquot until less than 0.5 mL remains. This is sufficient to remove most of the methanol. A thin film of oily residue will remain on the side of the test tube.

- j. Add 2 mL of borate buffer and mix briefly. A repipet dispenser or disposable pipette is sufficiently accurate to use for the buffer addition.

2. Liquid/Liquid Extraction

- a. Add 7 mL of ethyl acetate to the test tube containing the sample and vortex for at least 30 seconds. A repipet dispenser or disposable pipette is sufficiently accurate to use for the ethyl acetate addition.
- b. Centrifuge the tube for 5 minutes at 2000 rpm (approximately 560 g).
- c. Transfer the upper layer (ethyl acetate) into a clean 16 x 100 mm test tube using a disposable pipette or other suitable means, taking care to remove as much of the layer as possible *without removing any of the lower fraction* (borate buffer).
- d. Repeat steps 2.a. and 2.b. with fresh ethyl acetate. To save time, do not transfer upper layer. Proceed to solid phase extraction.

3. Solid Phase Extraction

Note: Steps a. - c. should be performed continuously without letting the cartridges run dry.

- a. Wet an acidic alumina SPE cartridge using approximately 5 mL of ethyl acetate. Let the solvent drain to the surface of the cartridge. Flow rate is not important.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 10 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

- b. Transfer the second ethyl acetate fraction from step 2.d. to the cartridge; follow with the first ethyl acetate fraction from step 2.c. Drain the combined ethyl acetate fractions to the surface of the SPE cartridge at a flow rate of approximately 2-4 mL/minute using vacuum if necessary.
- c. Wash the cartridge with approximately 5 mL of ethyl acetate at approximately the same flow rate as the sample application and stop the flow when the liquid reaches the surface of the cartridge packing. Discard the cartridge effluent.
- d. To elute ractopamine from the cartridge, add approximately 10 mL of methanol to the cartridge and collect the effluent in a 16 x 100 mm test tube or equivalent vessel. Force the liquid completely from the cartridge using either pressure or vacuum if necessary. Flow rate of methanol should be no greater than approximately 5 mL/minute.

Note: This is a suitable stopping point. The methanol effluent may be stored for 14 days at 2 - 8 °C if needed, before continuing with the method.

- e. Evaporate the sample to dryness using an air or nitrogen stream and a water bath or heater set at 49 ± 2°C.
- f. Dissolve the sample in 1.0 mL of sample diluent by swirling the tube or vortexing vigorously for 15 - 30 seconds, and sonicating for approximately 15 seconds in an ultrasonic water bath.
- g. Filter the sample through a PVDF 0.22 µm x 13 mm syringe filter using a small disposable syringe and collect the filtered sample in a HPLC autosampler vial.

Note: This is a suitable stopping point. Samples may be stored for 4 days at 2 - 8 °C if needed, until the sample analysis is completed.

4. Sample Analysis

a. Operating Conditions

Note: Instrument conditions may be modified to optimize chromatograph, if necessary.

- i. Column: 4.6 mm i.d. x 25 cm Supelcosil LC-18-DB,
(guard column may be used)

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 11 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

- ii. Excitation: 226 nm
 - iii. Emission: 305 nm
 - iv. Flow Rate: 1.0 mL/minute
 - v. Injection Volume: 100 µL
 - vi. Column Temperature: Ambient (20 - 25 °C)
 - vii. Run Time: 10 minutes
 - viii. Mobile Phase: Refer to section C.2.c.
- b. Initial test for HPLC system suitability (to be used when initiating analysis for the first time, or whenever instrumentation, HPLC column, or elution parameters are changed, or degradation of instrument performance is suspected).
- i. Inject an appropriate resolution solution and full standard curve curves for liver and muscle analysis, and verify system meets requirements specified in items ii – vi below.
 - ii. Ractopamine must be baseline resolved from ritodrine. See Figure 1, Section K.2, for example chromatogram.
 - iii. The signal/noise ratio of the 25 ng/mL standard for a liver set, or the 0.5 ng/mL standard for a muscle set, should be at least 5.
 - iv. The ractopamine peak should be well resolved from the solvent front and elute between 4 - 8 minutes.
 - v. The ractopamine retention time range for all injected standards (minimum to maximum) should not exceed 3% of the average value.
 - vi. The correlation coefficient that is calculated for the standard curve(s) must be ≥ 0.995 .
- Note: Decreasing acetonitrile concentration in the mobile phase increases the resolution of ractopamine from ritodrine and may separate impurities co-elution with ractopamine. Mobile phase composition should be optimized for each HPLC/column system.
- c. Daily HPLC suitability test (to be run before analyzing samples)
- i. Check system suitability by injecting standard curve(s) appropriate for the day's analysis. Verify that system meets requirements specified above in section b. iii – vi.
- d. Analyze Sample Set
- i. If above items are satisfactory, continue with injection scheme.

United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science

SOP No: CLG- RAC1.01		Page 12 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

A recommended injection sequence is as follows:

Resolution standard (optional)
Standards
Recovery
Blank
Samples
At least one standard

Note: It is recommended that the column be flushed with a strong solvent/water mixture (acetonitrile or MeOH/water, 80:20) after the end of the analytical batch to remove residual matrix from the column.

- e. Evaluate Sample Set Data
 - i. Measure the HPLC retention times and peak areas for ractopamine HCl in all sample set chromatograms.
 - ii. Identify a peak in any sample as ractopamine HCl if its retention time matches that of an appropriate reference standard within $\pm 3.0\%$ relative.
 - iii. If the peak area for ractopamine HCl in a sample exceeds the high end of the standard curve, the final extract should be diluted and re-injected along with one standard curve set. Another option is to take a smaller aliquot of the initial methanol extract and reprocess the sample. If the amount found exceeds known SPE cartridge capacity, the sample should be reprocessed using a smaller aliquot (step F.1.i) in order to assure a consistent recovery.

Note: The capacity of Waters Sep-Pak Plus Alumina A cartridges has been verified with 150 ppb recovery samples using 4x the recommended aliquot of the combined supernatant (step F.1.i).
 - iv. Verify that the standard injected at the end of the set meets the system suitability criterion b.iii. listed above.

5. Chromatograms

See Figures 1, 2 and 3 in section K.2

G. CALCULATIONS

- a. Using linear regression analysis, calculate the slope, intercept, and correlation coefficient of a standard curve constructed by plotting peak areas versus concentration (ng/mL) for all of the injected standards.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 13 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

Note: Standard curve must have a correlation coefficient (r-value) ≥ 0.995 over the concentration range used for quantitation.

- b. The concentration of ractopamine can be calculated using the following equation:

$$\text{ppb Ractopamine HCl} = \frac{[(A - B) \times E]}{(C \times F)}$$

A = HPLC peak area of sample injection

B = Intercept from the calibration curve

C = Slope of the calibration curve (area/ng/mL)

E = Total volume (mL) = (Initial volume/aliquot volume) x final volume

F = Weight of tissue sample (g)

Note: Area response should be linear with respect to the concentration of ractopamine. (Refer to section K.1. in the event acceptable recovery and blank control results are not obtained).

- c. Results may be reported only if the following conditions have been met:
- i. The correlation coefficient calculated for the standard curve is ≥ 0.995 .
 - ii. No quantifiable peak eluting within the elution time window for ractopamine is detected in the blank sample.
 - iii. The recovery of the positive control falls within the limits specified in section I.1.

H. HAZARD ANALYSIS

1. Required Protective Equipment —Safety eyewear, protective gloves, and lab coat.
2. Hazards

<i>Reagent</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
Glacial acetic acid	Strong acid	Wear protective equipment, avoid contact with skin.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 14 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

Ractopamine HCl	Eye irritant and exposure may increase heart rate.	Wear protective equipment, avoid breathing powder.
Ritodrine HCl	Irritant, fast or irregular heartbeat, nausea, shortness of breath.	See ractopamine HCl above.
Methanol, Ethyl Acetate	Flammable	Keep in well-closed containers in a cool place and away from fire. Use it in well-ventilated hood.

3. Disposal Procedures

<i>Reagent</i>	<i>Hazard</i>	<i>Recommended Safe Procedures</i>
Glacial acetic acid	Strong acid, burns	Collect waste in a tightly sealed container and store away from non-compatibles in a cool, well ventilated, acid liquid storage area/cabinet for disposal in accordance with local, State, and Federal regulations.
Ractopamine HCl	Eye irritant and exposure may increase heart rate.	Collect waste in a tightly sealed container and store in a cool, well-ventilated storage area/cabinet for disposal in accordance with local, State, and Federal regulations.
Ritodrine HCl	Irritant, fast or irregular heartbeat, nausea, shortness of breath.	See ractopamine HCl above.
Organic solvents	Flammable	Collect waste in a tightly sealed container and store away from non-compatibles in a cool, well ventilated, flammable liquid storage area/cabinet for disposal in accordance with local, State, and Federal regulations.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 15 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

I. QUALITY ASSURANCE PLAN

1. Performance Standard

<i>Analyte</i>	<i>Tissue</i>	<i>Analytical Range</i>	<i>Acceptable Recovery</i>	<i>Acceptable Repeatability (CV)</i>
Ractopamine HCl	Liver	25 - 300 ppb	60 - 115 %	≤ 20 %
Ractopamine HCl	Muscle	1.0 - 300 ppb	60 - 115 %	≤ 20 %

Acceptable correlation coefficient for standard curve: ≥0.995.

2. Critical Control Points and Specifications

<i>Record</i>	<i>Acceptable Control</i>
Sample weight	10.0 ± 0.2 g
Methanol	20 ± 1 mL
Combined methanol supernatant	60 ± 1 mL
Aliquot Volume	8 ± 0.1mL
Water bath temperature	49 ± 2 °C

3. Readiness To Perform (FSIS Training Plan)

a. Familiarization

i. Phase I: Standards - Duplicate standard curve on each of 3 consecutive days, which will include the following:

- (a) Level 1: 0, 25, 50, 75, 150, and 300 ng/mL
- (b) Level 2: 0, 0.5, 0.75, 1.0, 1.5, 2.0 ng/mL

ii. Phase II: Fortified Samples:

- (a) For liver: Three recovery curves in liver, on three separate days at 0, 25, 50, 150, and 300 ppb.
- (b) For muscle: Three recovery curves in muscle, on three separate days at 0, 1.0, and 2.0 ppb.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 16 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

Note: Phase I and II may be performed concurrently.

- iii. Phase III: Check samples given by the supervisor or designee.
 - (a) Seven unknown samples split between liver and muscle including one blank.
 - (b) Approval from the Supervisor of Record and the Laboratory Quality Assurance Manager (QAM) is required to commence official analysis.
 - b. Acceptability criteria.
Refer to section I.1. above.
4. Intralaboratory Check Samples
- a. Frequency:
 - i. One sample per week per analyst as samples analyzed.
 - ii. Records are maintained.
 - b. Acceptability criteria:
If unacceptable values are obtained, then:
 - i. Stop all official analyses by that analyst.
 - ii. Take corrective action.
5. Sample Acceptability and Stability
- a. Matrix: swine and bovine liver and muscle
 - b. Sample receipt size: approximately 500 g
 - c. Condition upon receipt: not spoiled or rancid
 - d. Sample storage:
 - i. Time: 2 months (This is the prepared sample.)
 - ii. Condition: Frozen <-10°C.
6. Sample Set
- Note: Each sample set must include:
- a. Blank tissue.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 17 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

- b. Fortified control
 - i. For liver
One blank tissue fortified with 1 mL of the 1.5 µg/mL fortification standard (addition of the fortification solution D.2.c.i. to tissue yields a 150 ppb concentration of ractopamine).
 - ii. For muscle
One blank tissue fortified with 1 mL of the 10 ng/mL fortification standard (addition of the fortification solution D.2.d.ii. to tissue yields a 1 ppb concentration of ractopamine).
 - c. Samples
7. Sensitivity
- a. Minimum proficiency level (MPL): 1.0 ppb for muscle; 25 ppb for liver

J. WORKSHEET

The worksheet, on the following page, is only an example and can be removed for photocopying.

United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science

SOP No: CLG- RAC1.01		Page 19 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

K. APPENDIX

1. SPE cartridge testing

In the event that acceptable standard recovery (> 60 %) and control results are not attainable using the method, the following steps may be taken to determine the suitability of the acidic alumina SPE cartridges.

Note: 2 g acidic alumina activity I from Waters, Alltech Isolute, and ICN proved to be acceptable using this test.

a. Preparation of Test Solutions

i. Ractopamine HCl fortification solution (0.2 µg/mL):

Pipet 2 mL of the 10 µg/mL intermediate solution into a 100 mL volumetric flask and dilute to volume with sample diluent. Mix well. (This solution is stable for one month when stored at 2 - 8 °C.)

ii. Ractopamine HCL external standard (2.5 ng/mL):

Pipet 10 mL of the 25 ng/mL standard, or 2.5 mL of the 100 ng/mL intermediate standard, into a 100 mL volumetric flask and dilute to volume with sample diluent.

b. Test Procedure 1

i. Fortify 2 mL of Borate Buffer with 0.1 mL of the 0.2 µg/mL fortification standard (20 ng).

ii. Perform steps F.2.a through F.3.g and analyze as described in the method.

iii. At least 85% (17 ng) of the fortification should be recovered.

iv. If 85% is not recovered, a new source of acidic alumina should be tested.

c. Test Procedure 2

The suitability may further be evaluated using the following procedure, which will determine the SPE performance in the presence of control and fortified liver tissue extracts.

i. Prepare control extracts of tissue following steps F.I.a through j.

ii. Fortify the desired number of samples with 0.1 mL of a 0.2 µg/mL fortification standard (20 ng).

Note: 20 ng fortified at the borate buffer stage is equivalent to 15 ppb in tissue.

iii. Perform steps F.2.a through F.3.g and analyze as described in the

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**


SOP No: CLG- RAC1.01		Page 20 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

method.

- iv. At least 80% (16 ng) of the fortified analyte should be recovered.
- v. The area of the tissue blank at the retention time of ractopamine should be less than 20% of the area obtained when injecting 100 µl of the 2.5 ng/mL external standard.
- vi. If 80% is not recovered, a new source of acidic alumina should be tested.

2. Chromatograms

Figure 1. Resolution of ractopamine and ritodrine (2.5 ng each injected on column) using a Supelcosil LC-18-DB column.



**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 21 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

Figure 2. Chromatograms showing blank and fortified liver tissues, respectively





Figure 3. Chromatograms showing blank and fortified, 1 ppb, muscle tissues, respectively



- Alternative HPLC column - C18 reversed phase HPLC column Beckman Ultrasphere IP, 5 µm, 4.6 x 250 mm, Part No. 235335.

**United States Department of Agriculture
Food Safety and Inspection Service, Office of Public Health Science**

SOP No: CLG- RAC1.01		Page 22 of 22
Title: Determination of Ractopamine Hydrochloride by High Performance Liquid Chromatography		
Revision: .01	Replaces: .00	Effective: 07/06/04

Approved By:	Date Approved:
Stephen Powell	06/23/04
Leon Ilnicki	06/23/04
David Martin	06/24/04
Jess Rajan	06/24/04
Phyllis Sparling	06/24/04
Charles Pixley	06/23/04