
Dispase Splitting Protocol	Berggren, Travis	
	tberggren@salk.edu	
Date Submitted	April 4, 2012	
Submitted by	Berggren, Travis	
	tberggren@salk.edu	
Adapted from	Salk STEM Cell Core in-house protocols	
Contributors	Lutz, Margaret. Modesto, Veronica.	
Affiliation	The Salk Institute	

Introduction:

This protocol is used for general maintenance and passaging of hES and iPS cells. It assumes that the cells are grown in 6 well format. The well to be split should be close to confluence. The colonies should be large and touching each other.

Flow Chart:



Protocol Steps:

Prepare the Matrigel Plate:

- 1. Thaw a 0.5mg vial of Matrigel on ice.
- 2. Bring up the thawed 0.5mg of Matrigel in 12 mls of ice cold DMEM/F12.
- 3. Immediately plate 2 mls per well of diluted Matrigel on a 6 well plate.
- 4. Plate must sit at room temperature for 1 hour or in a 37C incubator for at least 20 minutes before it can be used.
- 5. Matrigel plates can be made up to one week ahead of time and stored in the 4C refrigerator wrapped in Parafilm. Note: After removing from the 4C refrigerator, the plates must sit at room temperature for one hour or in a 37C incubator for 20 minutes before being used.
- 6. After Matrigel plate has sat for appropriate amount of time, aspirate off DMEM/F12
- 7. Place 1.5 mls of warm hES media into each well to be split into. Note: We recommend using mTeSR1 but any media formulated for use in feeder free conditions can be used.

Prepare the Dispase Solution:

- 8. Dissolve 2mg of Dispase into 1 ml of DMEM/F12.
- 9. Alternatively, make a 10x stock of Dispase (20mg/ml) and store in 1 ml aliquots in the -20C freezer. These aliquots can be thawed and diluted with 9mls of DMEM/F12.
- 10. Once in solution, Dispase can be stored in the 4C refrigerator for up to 2 weeks.

Passage the hES cells:

- 11. Mark areas of differentiation on the well to be split using the microscope objective marker.
- 12. Aspirate off spent media.
- 13. Place 1 ml of warm Dispase in each well to be split.

- 14. Incubate at 37C for 5-7 min. Note: After incubation, the edges of the colonies should be starting to curl up.
- 15. Aspirate off the Dispase.
- 16. Rinse the well gently 3 times with 1 ml warm DMEM/F12. Aspirate after each wash.
- 17. After the third rinse, aspirate off previously marked spots of differentiation with a Pasteur pipette.
- 18. Add desired amount of media for scraping off cells. For example: For a 1:4 split, add 2mls of media.
- 19. While holding pipette at a 90° angle, scrape a glass pipette back and forth across well while slowing expelling media to rinse off colonies.
- 20. Repeat this until at least 90% of colonies are detached from the well.
- 21. Pipette cells gently up and down twice to break up the colonies into smaller pieces.
- 22. Place appropriate amount of cell/media mixture into each well of the previously prepared Matrigel plate. For example: For a 1:4 split, scrape off in 2mls of media and plate 0.5mls into each well of the previously prepared plate.
- 23. Place in 37C incubator.
- 24. Give the plate a gentle shake back and forth and side to side to help distribute the colonies evenly in the wells. Tip: Remember to feed the wells not split in the old plate. These wells can be used to split the next day if this split did not go well.
- 25. Change the media the next day.

Materials:

Product	Supplier	Catalog #
Growth Factor Reduced	BD	354230
Matrigel		
Dispase	Life Technologies	17105041
DMEM/F12	Life Technologies	11330057
mTeSR1	Stemcell Technologies	05850
6 well plate	USA Scientific	
Microscope Objective	Nikon	MBW10020
Marker		

Serological Pipettes	Various Sources	
5 ml Glass Pipette	Fisher Scientific	13-678-27E
Pasteur Pipettes	Fisher Scientific	13-678-20D