


Title	Differentiating Neural Stem Cells into Oligodendrocytes
Date Submitted	May 5, 2012
Submitted by -	Efthymiou, Anastasia - anastasia.efthymiou@nih.gov
Adapted from -	Gibco Protocol
Contributors -	Efthymiou, Anastasia
Affiliation(s) -	NIH CRM - NIAMS – Laboratory of Stem Cell Biology

### ❖ Introduction:



Oligodendrocyte stained for PLP (green) and an intracellular marker (red) - Prof. Klaus-Armin Nave, Max-Planck-Institute for Experimental Medicine

### ❖ Protocol:

Neural stem cells (NSCs) will proliferate as progenitors a few times even after the complete growth medium is replaced with the appropriate differentiation medium. If the cells reach 90% confluency, it might be necessary to split the cells at a 1:2 ratio. However, do not split the cells once they reach day 9-10 of differentiation when they can get damaged during the passaging process.

1. Plate the NSCs on a polyornithine and laminin- coated culture dish in complete StemPro NSC SFM at  $2.5 \times 10^4$  or  $5 \times 10^4$  cells/cm<sup>2</sup>.
2. After 2 days, change the medium to oligodendrocyte differentiation medium. Change the spent medium every 3 to 4 days.

### ❖ Materials:

polyornithine and laminin-coated culture dish
StemPro NSC SFM

oligodendrocyte differentiation medium		
<b>StemPro NSC SFM Complete Media</b>		
Component	Final concentration	Amount
KnockOut™ D-MEM/F-12	1X	97 mL
GlutaMAX™-I Supplement	2 mM	1 mL
bFGF (prep as 100 µg/mL stock)	20 ng/mL	20 µL
EGF (prep as 100 µg/mL stock)	20 ng/mL	20 µL
StemPro® Neural Supplement	2%	2 mL
<b>Oligodendrocyte Differentiation Medium</b>		
Component	Final concentration	Amount
Neurobasal® Medium	1X	97 mL
B-27® Serum-Free Supplement	2%	2 mL
GlutaMAX™-I Supplement	2 mM	1 mL
T3	30 ng/mL	0.1 mL

❖ Troubleshooting:

❖ **References:**