



Title	Differentiating Neural Stem Cells into Neurons
Date Submitted	May 5, 2012
Submitted by -	Efthymiou, Anastasia - anastasia.efthymiou@nih.gov
Adapted from -	Gibco Protocol
Contributors -	Efthymiou, Anastasia
Affiliation(s) -	NIH CRM - NIAMS – Laboratory of Stem Cell Biology

❖ **Introduction:**



Human iPSC-derived neuronal culture: Neurons/Tuj1 (Green); Glia/GFAP (Red) Nuclei/DAPI (Blue) - Rakesh Karmacharya, MD, PhD (Broad/McLean) from Harvard Stem Cell Institute

❖ Protocol:

Neural stem cells (NSCs) will proliferate as progenitors a few times even after the complete growth medium is replaced with the appropriate differentiation medium. If the cells reach 90% confluency, it might be necessary to split the cells at a 1:2 ratio. However, do not split the cells once they reach day 9-10 of differentiation when they can get damaged during the passaging process.

1. Plate neural stem cells on a polyornithine and laminin- coated culture dish in complete StemPro NSC SFM at 2.5×10^4 or 5×10^4 cells/cm².
2. After 2 days, change the medium to neural differentiation medium. Change the spent medium every 3 to 4 days.
3. If expedited differentiation is desired, add 0.5 mM of dibutyryl cAMP (Sigma, Cat. no. D0627) to the differentiation medium daily starting at day 7 of differentiation for 3 days.

IMPORTANT! Do not expose cells to air at any time after they have differentiated into neurons.

❖ Materials:

polyornithine and laminin-coated culture dish		
complete StemPro NSC SFM medium		
neural differentiation medium		
dibutyryl cAMP		
StemPro NSC SFM Complete Media		
Component	Final concentration	Amount
KnockOut™ D-MEM/F-12	1X	97 mL
GlutaMAX™-I Supplement	2 mM	1 mL
bFGF (prep as 100 µg/mL stock)	20 ng/mL	20 µL
EGF (prep as 100 µg/mL stock)	20 ng/mL	20 µL
StemPro® Neural Supplement	2%	2 mL
Neural Differentiation Medium		
Component	Final concentration	Amount
Neurobasal® Medium	1X	97 mL
B-27® Serum-Free Supplement	2%	2 mL
GlutaMAX™-I Supplement	2 mM	1 mL

❖ Troubleshooting:

❖ References: